Zentral CENTRAL SERVICE STERILISATION

Security- Check

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Competency offensive – knowledge the key to success

Series of lectures by the KLS Martin Group, April 2013

Gudrun Westermann

The complexity of medical devices that are routinely cleaned and sterilized has grown dramatically in recent years. This relentless pace of change means that Central Sterile Supply Department (CSSD) personnel have to contend with endless new challenges. Medical device users must therefore bear this in mind and keep abreast of such changes. To that effect, the KLS Martin Group has launched a «competency offensive» to assist CSSD managers. This training event, which is held at various locations throughout Germany, was opened by Andreas Daum, CSSD manager at Odenwaldkreis Health Centre, with a lecture on how to deal with changes to instrument surfaces.

Instrument surface changes

He said that if instruments developed discolorations, this would not only negatively impact on their value and functionality, it would also make it more difficult to evaluate the cleaning results. Depending on the nature of the surface changes, there could also be negative implications for the patient; but it was not easy to determine this. Daum went on to say that there was no point in simply removing surface deposits without investigating their cause. Indeed the opposite was true – treatment measures that were too aggressive could damage the instruments' passive layer which was needed to protect their surfaces. Daum stated that to generate an adequate passive layer, new as well as repaired or refurbished instruments should be washed three times before being first used.

When investigating the cause of these problems, the water and steam quality should also be taken into account, in addition to the chemicals used (detergents and their residues, care agents) and the general method of handling instruments, in the operating room (OR) too.

Daum gave an overview of the various types of corrosion as well as of silicate deposits and discolorations arising from

rouging. He gave tips on preventing and removing these. While rouging was a topical issue, it occurred primarily at temperatures above 160 °C and at pressure levels higher than 4 bar. As such, stated Daum, this problem was not too commonly encountered in medical device reprocessing. But cleaning the sterilizers and checking the steam quality were of paramount importance.

Routine reprocessing activities

The next lecture, delivered once again by Andreas Daum, focused on routine reprocessing activities. He spoke about the various steps of which the reprocessing cycle consisted, while emphasizing important aspects of the different phases. It was advisable to compile standard operating procedures (SOPs) for management of medical devices after use at the site of use. OR staff needed commensurate training to ensure e. g. devices were properly pre-cleaned and dismantled.

If devices were manually reprocessed (as admissible for Critical A devices, e. g., if there was no washer-disinfector (WD) available, such as in medical practitioners' offices), detailed SOPs were needed and should be observed. Disinfection should be carried out with virucidal agents covering the entire spectrum of viruses.

But preference should always be given to automated reprocessing. It was important to check the validation report and precisely document the results of routine checks and report these in the event of any discrepancies. More in-depth tests were needed for equipment that did not conform to the valid regulations. That could be accomplished for instance by using cleaning indicators and conducting measurements based on thermologgers and pressure loggers.

If rinse aids were used, their compatibility with the medical device materials had to be ensured. Such incompatibilities could arise, in particular, in the case of synthetic materials, which were particularly suscep-

tible since these chemical rinse aids were applied to the medical devices in the last rinse cycle.

The packaging must be chosen in line with the sterilization process. Daum stressed that every type of packaging had to be taken into account when validating the sterilizers. No other types of packaging could be used for sterilization.

In the case of containers, the requirements governing the filter system used had to be taken into account. For reusable filters that could be used only for a certain number of reprocessing cycles, the respective number of cycles should always be documented. But no such time limit applied for valve systems, however, these required servicing.

Only qualified personnel were authorized to release reprocessed medical devices. Daum stated that this called for commensurate expertise.

Inspections – be prepared!

In the second part of his talk Daum spoke about the inspections carried out by the supervisory authorities. What was important here was to be well prepared, e. g. to check that all documentation was complete and up to date (e. g. manufacturer's instructions). The employees, too, needed to be prepared because they could be asked various questions by the supervisory authorities. Prior training in the most important procedural instructions was advisable.

It was also beneficial to have the CSSD management inspect their own premises in advance, paying special attention to cleanliness and safety aspects, as well as to the shelf-life data of process chemicals. Daum finished off by stressing that the authorities could only demand things that were demonstrably required in writing

*: Recommendation for hygienic reprocessing practices for medical devices, jointly compiled by the Commission for Hospital Hygiene and Infection Prevention at the Robert Koch Institute (RKI) and the Federal Institute for Drugs and Medical Devices (BfArM)

(e. g. in the KRINKO-BfArM Recommendation*). In the event of any discrepancies, he recommended contacting the hospital infection control consultant. Adopting an open and friendly approach was helpful in all matters.

I Process validation – what are the requirements

Ralf Wesch, from the firm Vanguard, dealt with the requirements for process validation. He explained its various components, i. e. installation, operational and performance qualification and when these were to be carried out. He used photos of routine validation activities to demonstrate various points.

To test the cleaning efficacy, instruments harbouring defined test soils as well as everyday soils should be used. The maximum period of time instruments were likely to be left unattended after use in routine practice until they were reprocessed was important. If any changes were

made to the cleaning detergents, then the process had to be revalidated immediately (performance requalification), since otherwise reprocessing would be deemed to have been conducted with a non-validated process.

For sterilization Wesch pointed out that routine operations could not be conducted concurrently with validation. He advocated carrying out servicing before validation (no more than six weeks before). The prescribed annual calibration of sterilizer sensors could also be performed at the same time.

All types of packaging used had to be taken account of at the time of validation, if necessary running different programmes. Wesch also elaborated on the microbiological validation needed additionally in certain cases, e. g. if the structure of the sterile supplies did not permit the use of sensors at critical sites. That meant that the instruments used for validation would not be available for reuse for around two weeks.



Wesch finished off by summarizing the various sections of a validation report and how these were best organized. He recommended having the entire reprocessing process assessed at the time of validation and taking advantage of the presence of the validation inspector by asking him to clarify any uncertainties.

The RKI-BfArM Recommendation was also intended as a means of enhancing reproducibility. Wesch cited by way of example some of the changes made to the amended version of the recommendation published in 2012.

I Choosing a sterile barrier system

The focus of a second lecture, given also by Ralf Wesch, was on comparing soft packaging with hard packaging. The user should in principle think carefully when choosing a sterile barrier system. How much space was available in the OR, in the sterile supply warehouse and in the CSSD? For sheet packaging materials, for example, larger

The evolution of the revolution – 2nd Generation

Record Keeping

Reprocessing: objective – physical – digital

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3M™ Electronic Test System (ETS)

- Objective Monitoring
- Electronic Archiving & Controlling
- Unambiguous Result (»Pass« or »Fail«)
- Reusable
- Professional Service

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3M – The Partner for your SSD
Call us for a free Trial:
Customer Services: 0 15 09 - 61 33 27 (UK)

packing tables were needed. Was there enough reprocessing capacity available to warrant the use of containers? Economic and ecological aspects also played a role. Wesch demonstrated the differences between the packaging types used for packing, sterilization and storage. Whereas containers could be easily stacked, were stable and rugged, soft packaging was less resistant to damage and needed more space – when stored, or even for the packing process itself. The «human factor» also played a pivotal role in packing and was a potential source of errors. And, not least, this type of packaging resulted in significantly more waste which, besides, had to be treated as contaminated, in particular when used devices were wrapped in non-woven materials after use. On the other hand, the investments needed for soft packaging were lower. There was no need for servicing or for documentation of cycles. The packaging sizes were variable and easier to handle, not least because of the lower weight. Conclusion: both systems had their advantages and disadvantages that had to be taken into consideration in the individual case.

I Trends and perspectives

In the last lecture, Johannes Gulde, product manager of the KLS Martin Group, reported on trends and perspectives in the CSSD. First, he spoke about sustainability. But what did that mean in the CSSD, an area where high energy and water consumption as well as large amounts of waste could, in principle, be expected? But Gulde continued by stating that despite all these reservations, ecological gains could be made if resources were used more efficiently, including the use of a more environmentally friendly generation and recycling of the resources used. This also involved economic aspects. Alone the optimization of communication could save time and resources. A more efficient process throughput could be achieved through astute organization of workplaces and appropriate space allocation. Gulde stressed that keeping operating costs under control was vital. These were often overlooked, but could be the source of potential savings. Top priority had to be accorded to patient protection. But responsibility had also to be borne for the employees, e. g. by assuring effective occupational safety, continuing professional development and an ergonomic workplace.

I Ecological assessment of a CSSD

Gulde presented a project called «Ecological assessment of a CSSD». Ecological assessment is aimed at identifying the environmental pollution caused by products throughout their entire life cycle and analysing the impact of such influences.

Based on data from two hospitals, Gulde demonstrated the workflow patterns within a CSSD. The data collected were based on the performance characteristics of large apparatuses, such as e.g. their water and power consumption, but also the proportion of non-woven packaging used in the process and ensuing amount of waste.

Armed with these data, effective optimization measures could be identified and implemented. The most important points mentioned by Gulde are as follows:

1. Large sterilizers. Appliances with greater capacities consume on average fewer resources per batch.
2. A more circumspect management of packaging materials. Consumption should be kept to a minimum. Container systems reliant on consumables, e. g. paper filters and seals, account for more than one-third of all waste. Containers that do not need consumables produce only one-tenth of waste. Hence such a system helps to cut down by around 2.1 tons the amount of waste that would otherwise be generated during its service life.
3. Optimization of programmes for washing containers so as to use fewer resources. This is greatly underpinned by the use of a special tunnel washer for containers, but this is often not used because of lack of space and it is not economical for smaller establishments that have a less fast container turnaround time.
4. Formulation of sustainable economic policies. These imply that, on the part of all employees, responsibility for resources is presupposed, underpinned and required.
5. Verification of the energy mix in place. Is sustainable production of this energy mix assured or not?

I RFID and automation

Other trends reviewed by Gulde in his talk included radio frequency identification (RFID) and automation. First of all he explained the characteristics of RFID, comparing this method also with optical identification systems such as DataMatrix or Barcode, and reported on experiences of clinical application at a university hospital

in Osaka, Japan. There all medical devices were marked with DataMatrix, while trays and containers were additionally equipped with RFID tags. At the same time the system was linked to a suitable software facility permitting exploitation of important advantages:

- Display of shelf-life data, with corresponding priority assignment of trays
- Display of tray availability
- Identification of the localization of trays in real time
- Recording of use frequency with statistical evaluation

The results showed that RFID gave rise to improved process sequences, efficient management of repairs and to reduction of «unplanned» costs. Search and handling charges were reduced and disintegration of instrument trays kept to a minimum. Against a background of some 6,000 operations, it was possible to save almost 140,000 euros in this way.

Gulde finished off by stating that RFID was highly cost intensive at present and would definitely have to be equipped with automation facilities in order to realize its full potential.

In the last part of his talk Gulde spoke about the prospects held out by automation, while basing his insights on a CSSD in Gentofte, Denmark. But not all processes could be automated, e. g. requests from the OR, provision of supplies and return and storage after sterilization. What kind of experiences had been gathered from such an automated CSSD? Process optimization measures, e. g. through shorter supply times and optimization of warehouse management permit savings of 0.5 million euros for the entire process per year. To begin with, these are associated with high installation costs of 3.8 million euros. Apart from the long-term financial benefits, the system also had implications for staff and patients: contact with contaminated materials was reduced and manual processes kept to a minimum. Complete digital management of the process documentation and tracking was possible – a plus point for patient safety, which after all was the top priority.

Conclusion: development of sustainability in the CSSD is a topical issue. Whether the prospects outlined here come to fruition remains to be seen. ■