
UHP

Herbert-Ulnakopfprothese®
Herbert Ulnar Head Prosthesis®

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Timothy -J. Herbert/Jörg van Schoonhoven

Disorders of the Distal Radio-Ulnar Joint and the Need for a New Ulnar Head Prosthesis

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Fig. 1: Rotation of the radius on the ulna.
 From: Bowers WH The distal radio-ulnar joint.
 In: Green DP, ed Operative Hand Surgery.
 2nd. ed. New York: Churchill Livingstone, 1988: 939-990.

Introduction

Disorders of the distal radio-ulnar joint (DRUJ) are common and are usually associated with pain, weakness, instability and loss of forearm rotation. The aim of surgical treatment is to restore a painfree range of forearm rotation with a stable ulna and wrist. Today, there is no single surgical procedure which can be guaranteed to restore painfree rotation and at the same time correct any underlying instability of length discrepancy. For this reason treatment of the DRUJ remains controversial.

The ulnar head prosthesis has been designed to restore pain-free forearm rotation at the same time as any length discrepancy is corrected, stability is achieved by means of a simple soft tissue repair (so called "ulnar-flap"), which is combined with a reconstruction of the triangular fibro-cartilage complex (TFCC).

Anatomy

The DRUJ consists of a shallow "ball and socket" articulation between the head of the ulna and the sigmoid fossa of the radius. This articulation allows for 180 degrees of forearm rotation (radius on ulna, Fig. 1), which is essential for normal function of the hand.

The TFCC acts as a cushion between the head of the ulna and the proximal row of the carpus. By means of its strong soft tissue attachments to the apex of the ulnar head and the capsule of the DRUJ, it is also the major stabiliser of the DRUJ. Similarly, its distal expansions merge with the wrist capsule to stabilise the ulnar side of the carpus on the wrist. It allows for a free range of flexion, extension, ulnar deviation and radial deviation at the ulno-carpal joint.

Bowers (1991) provides an excellent review of the anatomy and function of the DRUJ in the normal and unstable wrist.

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Disorders of the DRUJ

Any loss of congruity between the sigmoid fossa and the ulnar head will result in painful loss of forearm rotation. Causes include congenital abnormalities such as Madelung's deformity, radial fractures, inflammatory arthritis or degenerative arthritis as a consequence of tears of the TFCC.

Tears of the TFCC are extremely common, particularly in patients with ulnar plus variance (81%, Palmer and Werner 1981). Tears may be due to acute or chronic trauma, and lead to pain and loss of motion in the ulnar side of the wrist; loss of integrity of the TFCC results in instability of both the DRUJ and the ulno-carpal joint. In time, this is likely to become a fixed deformity leading to secondary osteoarthritis.

As the TFCC fails, the normal alignment between the distal radius and ulna is lost and x-rays will show evidence of increasing ulnar plus deformity to the extent that the ulna starts to impact on the lunate and/or triquetrum. This impaction may result in triquetro-lunate instability and ulno-carpal osteoarthritis.

The commonest cause of these problems is fracture of the distal radius. This may be associated with avulsion fracture of the ulnar styloid process, leading to functional instability of the TFCC. More commonly, radial malunion results in a secondary ulnar plus deformity. Occasionally direct damage to the sigmoid fossa may lead to the development of secondary osteoarthritis.

Furthermore trauma may result in acute dislocation of the DRUJ itself, often in association with radial head fracture and tear of the interosseous membrane (Essex-Lopresti lesion).

Other conditions affecting the DRUJ and TFCC include growth disorders (Madelung's deformity), primary osteoarthritis, metabolic disorders (e.g. gout), rheumatoid arthritis. Finally, many patients present with painful instability following previous surgical procedures on the DRUJ.

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Current operations for DRUJ disorders


Abb. 1: Ulnakopfresektion nach Darrach. R=Radius, U=Ulna, T=Triquetrum
Fig. 1: Darrach's procedure. R=radius, U=ulna, T=triquetrum

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1. Excision of the ulnar head (Darrach's procedure Tulipan et al. 1991, Fig. 1) destabilizes the TFCC leading to ulno-carpal and radio-ulnar instability. Absence of the ulnar head often results in painful radio-ulnar impingement.

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2. Bowers procedure (Bowers 1985, Fig. 2) attempts to overcome the problems associated with Darrach's procedure as outlined above. By preserving the attachment of the TFCC to the ulnar styloid process, some stability is retained. However, soft tissue interposition between the sigmoid fossa and the resected surface of the ulna is seldom strong enough to prevent some degree of ulno-carpal and radio-ulnar instability, due to loss of the normal alignment between the radius and the ulna. Furthermore, in cases of ulnar plus variance, ulno-carpal impaction will occur unless the ulna is shortened at the same time.

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3. The Sauve-Kapandji procedure (Sanders et. al. 1991, Fig. 3), if performed correctly, will relieve any ulno-carpal impaction whilst maintaining function of the TFCC. However, rotation is only restored at the expense of an unstable ulnar pseudoarthrosis which commonly results in painful radio-ulnar impingement. Furthermore the exact positioning of the radio-ulnar fusion is difficult, with the result that many patients still have ulno-carpal impaction and TFCC problems.


Abb. 4: Der „ulnargestielte Weichteilappen“. Aus: „Stanley D. und T.J. Herbert, The Swanson Ulnar Head Prosthesis for post-traumatic disorders of the distal radio-ulnar joint“, Journal of Hand Surgery, 1992, 17b, S. 682-688.

Fig. 4: The "ulnar-flap". From Stanley D. Herbert T.J. "The Swanson Ulnar Head Prosthesis for post-traumatic disorders of the distal radio-ulnar joint". J. Hand Surg. 1992; 17b: 682-688

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4. Silastic ulna head replacement (Swanson 1973) was designed to prevent radio-ulnar instability following ulnar head resection. With modifications previously described (Stanley and Herbert 1992, Fig. 4) it may also allow for reconstruction of the TFCC and restoration of ulno-carpal instability. However, the silastic prosthesis is subject to excessive wear which may be complicated by the development of silicone synovitis. Many of these prostheses have had to be removed resulting in painful ulno-carpal instability.

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At the present time, most surgeons prefer either the Bowers or the Sauve-Kapandji procedure. Darrach's procedure is normally reserved for low demand patients (e.g. rheumatoid arthritis), whilst the Swanson prosthesis is no longer available and has fallen out of favour. When correctly carried out, and in the appropriate patients, any of the before mentioned procedures can produce reasonably satisfactory results. However, there are many patients who are unhappy with the results of surgery and who are seeking a suitable revision procedure to restore stability and painfree rotation at the DRUJ.

It is for these reasons that we sought to develop a new ulnar head prosthesis which, combined with a careful reconstruction of the TFCC and the use of an ulnar based capsulo-retinacular flap, could overcome all of the shortcomings associated with previous surgical procedures.

Design of the prosthesis

Preliminary biomechanical studies were carried out by the Author in conjunction with the Department of Biomechanical Engineering at the University of New South Wales. Stem and head dimensions were determined by x-ray and CT-Scan measurements on 100 normal wrist joints. It was found that the full range could be covered by having 3 stem diameters and 3 head sizes. Further, it was recognized that there would be a need for a revision stem with a built up neck to compensate for excessive shortening associated with previous surgical procedures.

It was felt that the stem should be designed to have a "contact fit" within the medullary canal of the ulna, with a porous surface for bone ingrowth so that cement fixation would not be required. Interchangeable stem and head sizes demanded the use of a "morse lock" between the two components. The head is designed to conform with the diameter of the sigmoid fossa, restoring the normal alignment between the radius and the ulna, and providing a support underlying the soft tissue flap repair. The distal end of the head is concave in order to avoid excessive pressure on the underside of the TFCC. Ceramic was chosen as the most biocompatible material with the best wear characteristics for use as a hemiarthroplasty.

Finally, it was not felt necessary to have any soft tissue attachment points on the prosthesis, (i.e. to reproduce an ulnar styloid process), as the flap repair has been shown to provide adequate tension of the TFCC and give excellent stability both in the distal radio-ulnar and ulnocarpal joints (Fig. 5).

A simple and elegant set of instrumentation has been developed to facilitate insertion of the prosthesis, to ensure accurate length correction and sizing of the prosthesis.

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Fig. 5: The KLS Martin prosthesis in situ, stabilised by the "ulnar-flap"

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Reference List

- Herbert TJ, van Schoonhoven J:
Ulnar Head Replacement
Techniques in Hand and Upper Extremity Surgery 2007;
11(1):98-108
- Van Schoonhoven, Herbert TJ:
The Dorsal Approach to the Distal Radioulnar Joint
Techniques in Hand and Upper Extremity Surgery 2004;
8(1): 11-15
- Grechening W, Peicha G, Fellingner M:
Primary Ulnar Head Prosthesis for Treatment of
an Irreparable Ulnar Head Fracture Dislocation.
Journal of Hand Surgery 2001; 26B: 3: 269-271.
- Herbert TJ, van Schoonhoven J:
Ulnar Head Prosthesis: A new solution for problems
at the distal radioulnar joint.
Hand Anthroplasties 2000; 145-149.
- van Schoonhoven J, Fernandez DL, Bowers WH, Herbert TJ:
Salvage of Failed Resection Arthroplasties of the Distal Radio-
ulnar Joint using a new Ulnar Head Prosthesis.
Journal of Hand Surgery 2000; 25-A: 3: 438-446.
- Fernandez DL:
Acute and Chronic Derangement of the Distal
Radio Ulnar Joint after Fractures of the Radius.
EFORT 1999; 41-53
- Stanley D, Herbert TJ:
The Swanson ulnar head prosthesis for
post-traumatic disorders of the distal radio-ulnar joint.
Journal of Hand Surgery 1992; 17B: 682-688.
- Bowers WH:
Instability of the distal radio-ulnar articulation.
Hand Clinics 1991; 7: 311-327.
- Sanders RA, Frederick HA, Hontas RB:
The Sauve-Kapandji procedure: A salvage operation
for the distal radio-ulnar joint.
Journal of Hand Surgery 1991; 1125-1129.
- Tulipan DJ, Eaton RG, Eberhart RE.:
The Darrach procedure defended: Technique redefined
and long-term follow up.
Journal of Hand Surgery 1991; 16A: 438-444.
- Bowers WH:
Distal radio-ulnar joint arthroplasty:
The hemiresection-interposition technique.
Journal of Hand Surgery 1985; 10A: 169-178.
- Palmer AK, Werner FW:
The triangular fibrocartilage complex
of the wrist – Anatomy and function.
Journal of Hand Surgery 1981; 6: 153-162.
- Swanson AB:
Implant arthroplasty for disabilities of the distal radio-ulnar joint:
Use of a silicone rubber capping implant
following resection of the ulnar head.
Orthop. Clinics 1973; 4: 373-382.

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Ulnar Head Prosthesis for Disorders of the Distal Radio-Ulnar Joint

Indications

- Salvage of:
 - Darrach
 - Bowers
 - Sauve-Kapandji
- Primary osteoarthritis
- Post-traumatic osteoarthritis:
 - Radial fractures
 - TFCC tears
 - Ulnar impingement
- Rheumatoid Arthritis
- Tumors

Contraindications


- Excess radial deformity
- Inadequate bone stock
- Inadequate soft tissues

Pre-operative assessment

Careful clinical examination is required in order to assess acutely the extent of instability of the ulna and/or the carpus. Similarly it is important to recognize any fixed deformity or dislocation, as this will need to be corrected at the time of surgery. If dorsal or volar dislocation of the ulnar head is due to an underlying radial deformity, this should be corrected by means of an appropriate radial osteotomy before proceeding to ulnar head replacement.

Standard 90/90 degree x-rays of both forearms and wrists have to be obtained preoperatively. The appropriate x-ray template, (allow for magnification factor, normally 1:1.1), should then be used to assess optimum resection level and likely size of head and stem required. Comparison with the opposite normal wrist is the most helpful way to determine the appropriate resection level.

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Set-up

The operation is carried out with the arm extended on an operating side table in full pronation and under tourniquet control. X-ray control, using a fluoescan is highly desirable.


Incision (Fig. 1)

A 6–8 cm dorsal incision is made, centred over the distal radio-ulnar joint (DRUJ).

In the case of a revision procedure, the incision should take account of previous scars and also extend proximally to allow adequate exposure of the ulna shaft.

The dorsal sensory branch of the ulnar nerve is found and retracted together with the volar skin flap. Be aware of the possibility of a transverse branch of this nerve passing towards the radial side of the wrist at the level of the ulno-carpal joint. The 5th extensor compartment is opened longitudinally and the Extensor Digiti Minimi (EDM) tendon is mobilized and retracted radially.

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
*Raising the flap
(Figure 2)*

A capsulo-retinacular flap is marked out as shown in Fig. 2. The base of the flap is centred on the pisiform bone with its apex at the mid-point of the DRUJ. Start raising the flap by incising through the capsule where it attaches to the rim of the sigmoid fossa. The flap is then carefully raised by sharp dissection off the dorsal part of the TFCC distally and off the neck of the ulna proximally.

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
Stay sutures are used to retract the flap, allowing exposure of both the DRUJ and ulno-carpal joints. Note that the 6th extensor compartment (Extensor Carpi Ulnaris – ECU) should not be opened as this forms an integral part of the base of the flap (Fig. 3).

If the floor of the 6th compartment is deficient (degenerative process or previous surgery), this should be repaired using fine interrupted mattress sutures so that the ECU tendon is no longer visible but is contained within the base of the flap.

For revision procedures the soft tissues between the ulna proximally and the near edge of the flap distally, are incised longitudinally and reflected as one layer to expose the ulnar stump.

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Ulnar head resection (Figure 4)

The neck of the ulna is exposed using Hohmann retractors. The resection level is marked using the resection guide (26-265-05-07). This ensures that the head of the prosthesis will lay in the optimum position, 2 mm proximal to the TFCC. The tip of the resection guide is carefully inserted over the distal end of the radius at its extreme ulnar border (Fig. 4).

Level 1 resection is used for the standard prosthesis and is the preferred resection level.

Level 2 may be used when there is insufficient bone stock to allow for seating the prosthesis more distally (e.g. after Bowers' procedure). The standard+ stem should be used at this level in order to ensure optimum length of the prosthesis (at 2 mm ulnar minus).

Because of the length discrepancy between the three head sizes (+2mm/-2mm) and the length discrepancy of the standard and standard+ stem (2 mm) all combinations can be performed with resection level 1 and 2, as shown in table 1.

Level 3 is used only with the revision stem, in case of extreme ulnar shortening (e.g. after previous surgery such as Darrach's procedure). Resection should be carried out at the appropriate level, depending on the size of head selected.

- s: small
- m: medium
- l: large

The ulnar head is then grasped with the holding forceps and is rotated externally and distally so that the soft tissue attachments can be carefully peeled away using a combination of blunt and sharp dissection. In this way the attachment of the apex of the TFCC to the underside of the flap will be preserved. If there is an ununited ulnar styloid fragment, this may need to be removed separately and any resulting deficiency carefully repaired.


Kopf	klein	mittel	groß
Schaft	26-220-01-04	26-220-03-04	26-220-05-04
Standard		Niveau 1	Niveau 2
Standard+	Niveau 1	Niveau 2	

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Head	small	medium	large
Stem	26-220-01-04	26-220-03-04	26-220-05-04
Standard		Level 1	Level 2
Standard+	Level 1	Level 2	

Table 1: Possibilities of combinations at the different resection levels

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
Reconstruction of triangular fibrocartilage (Figure 5)

Following resection of the ulnar head, the triangular fibrocartilage (TFC) is carefully inspected and palpated on both its proximal and distal surfaces. Any tears or deficiencies are repaired using fine 4/0 non-absorbable sutures. Where necessary, local soft tissue flaps may be used to reinforce the repair of the triangular fibro cartilage.

Assessment of joints

The triquetro-lunate joint is checked and stabilised as necessary. The sigmoid fossa is then examined and cleared of any scar tissue or osteophytes which may interfere with seating of the head of the prosthesis. If necessary the sigmoid fossa may be deepened with a power reamer in order to improve seating at the head of the prosthesis and/or to enhance stability.

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Reaming of ulna (Figure 6)

Place the appropriate sided (right/left) soft tissue retractor (26-265-20-07, 26-265-22-07) under the distal end of the ulna and open the medullary cavity using the sharp awl (22-130-14-07). The standard small reamer is hammered into its maximum depth, at which stage it should be a tight fit; if not, progress up to the medium or large reamer as appropriate. The slot of the hammer (23-264-19-07) is used to remove the reamer, which should always have been inserted to its full depth.

In the case of a revision procedure, use only the special revision reamers in similar fashion.

At this stage at least two small drill holes should be made in the dorsal rim of the sigmoid notch, to facilitate insertion of perosseous sutures during later flap repair.

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Trial reduction and assessment

Insert the appropriate trial stem and head. Check that the stem fits snugly and that the head lies at the appropriate level (approx. 2 mm) below the TFCC. This should be checked fluoroscopically with the wrist in full pronation.

Advance the flap radially over the ulnar head to check that this can be re-attached to the radius and provide adequate stability.

If the flap is too tight, a smaller head should be tried. If this is selected, then the standard stem should be changed to a standard+ stem in order to ensure optimum length of the prosthesis.

If the flap is too loose, a larger head may be selected. In this case, the ulna must be shortened back to level 2 in order to ensure optimum length of the prosthesis (2 mm ulnar minus).

Range of forearm rotation is assessed at this stage, and if it remains limited, the cause should be sought and dealt with as appropriate (e.g. osteophyte block; contracture of interosseous membrane, etc.).

The trial prosthesis is removed, using the explantation chisels (26-265-35-07, 26-265-37-07) and the wound and medullary canal are irrigated.

Inserting the definitive prosthesis

The appropriate stem (standard, standard+ or revision) and head (small, medium or large) are selected. The stem is carefully hammered into position using the conical impactor (26-265-30-07). The conical end of the stem should be clean and dry before the ceramic head is impacted.

At this stage, intra-operative x-rays are used to confirm that the prosthesis is in optimum position.


Warnung:

Excessive force should never be used during insertion of the definitive prosthesis. If the stem appears to "lock" during insertion, it should be removed and further reaming carried out.

Removal of prosthesis

If for any reason it is necessary to remove the prosthesis, this can be done using the stem extractor (26-265-07-07) provided in the set.

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Flap Closure (Fig. 7)

The dorsal aspect of the TFCC is reattached to the underside of the flap using one or two fine non-absorbable mattress sutures. The flap is then carefully advanced until it is tight enough to provide complete stability whilst still allowing full range of rotation. Reduction is facilitated if the assistant pronates the carpus (pushes up on the pisiform) and reduces the ulna (pushes down on the ulnar shaft).

The flap should be reattached, under the appropriate tension, to the dorsal rim of the sigmoid notch, using a perosseous suture technique (size 0 non-absorbable suture material is recommended). The remainder of the flap is sutured using an overlap (2 layers) repair where appropriate. The EDM tendon should be left free and no attempt is made to reconstruct the 5th extensor compartment. The proximal extension of the flap along the ulnar neck is sutured as is the distal extension overlying the ulno-carpal joint. At this stage a final assessment is made to check that there is a full range of rotation with no abnormal crepitus or instability in either the radio-ulnar or ulno-carpal joint.

Closure

Suction drains are inserted. The wound is closed in an appropriate fashion and a firm wool and crepe bandage applied. A 6-inch plaster slab is then moulded around the ulnar side of the wrist in order to protect the soft tissue repair.

Post operative management

Dressings are normally left intact until the sutures are removed at 12–14 days. However, where necessary, the dressings may be reduced and the appropriate splints and hand therapy commenced after 48 hours. Where stiffness is likely, early active and passive exercises should be commenced, alternatively, if instability is more likely to be a problem, a “sugar tong” type splint is fitted as soon as possible in order to control forearm rotation and protect the repair. Full clinical and radiological assessment is carried out after 5–6 weeks. At this stage splintage should be discarded and active mobilising and strengthening exercises encouraged. Where appropriate an elastic ulnar carpal support may be used as protection particularly for heavy workers.

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A case of radial malunion treated by hemiresection arthroplasty of distal radio-ulnar joint. Note persistent ulno-carpal impaction. Patient complained of chronic pain and instability at DRUJ and wrist.

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X-Ray following revision with standard ulnar head prosthesis.

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Normal radio-ulnar alignment in lateral view indicates stability. Pain and instability relieved and the patient regained normal function at the DRUJ.

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Unsatisfactory Darrach's procedure with painful instability.

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
Revision prosthesis in place. Note correction of ulnar length which, combined with the use of the ulna sling, resulted in restoration of normal stability and range of motion.

Implantate
Implants

STERILE !

Standard-Kragen
Standard collar

26-210-01-09

klein
small

Ti

26-210-03-09

mittel
medium


Ti

26-210-05-09

groß
large

Ti

Standard-Kragen +
Standard collar +

26-210-11-09

klein
small

Ti

26-210-13-09

mittel
medium


Ti

26-210-15-09

groß
large

Ti

Revisions-Kragen
Revision collar

26-210-21-09

klein
small

Ti

26-210-23-09

mittel
medium

Ti

26-210-25-09

groß
large

Ti

Keramik-Köpfe
Ceramic heads

26-220-01-04

klein
small

Ce

26-220-03-04

mittel
medium

Ce

26-220-05-04

groß
large

Ce

alternativ: Kobalt-Chrom-Köpfe (CoCr)
alternative: Cobalt-chromium heads (CoCr)

26-230-01-04

klein
small

CCr

26-230-03-04

mittel
medium

CCr

26-230-05-04

groß
large

CCr

Probeimplantate
Trial implants

Ti St Ce CCr STERILE ↓
 Titan Stahl Keramik Kobalt-Chrom dampfsterilisiert
 Titanium Steel Ceramic Cobalt-chromium steam sterilized

Standard-Kragen
Standard collar

26-211-01-05

klein
small

St

26-211-03-05

mittel
medium

St

26-211-05-05

groß
large

St

Standard-Kragen +
Standard collar +

26-211-11-05

klein
small

St

26-211-13-05

mittel
medium


St

26-211-15-05

groß
large

St

Revisions-Kragen
Revision collar

26-211-21-05

klein
small

St

26-211-23-05

mittel
medium

St

26-211-25-05

groß
large

St

Probe-Köpfe
Trial heads

26-221-01-05

klein
small

St

26-221-03-05

mittel
medium

St


26-221-05-05

groß
large

St

26-265-13-07
 15 cm/ 6"
 Knochenhaltezange
 Bone holding forceps


23-730-17-07
 17 cm/ 6¾"
 Knochenhaltezange
 Bone holding forceps

26-265-15-07
 21 cm/ 8¼"
 Scharfer Löffel
 Sharp spoon

22-130-14-07
 14 cm/ 5½"
 Knochen-Pfriem
 Sharp awl

26-265-05-07
 16 cm/ 6¼"
 Resektions-Messinstrument
 Resection guide

23-264-19-07
 19 cm/ 7¾"
 Hammer
 Hammer

26-265-20-07 rechts/right
26-265-22-07 links/left
 11,5 cm/ 4¾"
 Weichteilretractor
 Soft tissue retractor

Standard-Raspel
Standard reamer

1/2

26-250-01-07
klein
small

26-250-03-07
mittel
medium

26-250-05-07
groß
large

Revisions-Raspel
Revision reamer

1/2

26-260-01-07
klein
small

26-260-03-07
mittel
medium

26-260-05-07
groß
large

1/2

26-265-30-07
9,5 cm/ 3 3/4"
Einschlaginstrument
Conical impactor

1/2

26-265-35-07 klein/small
26-265-37-07 groß/large
24,5 cm/ 9 1/2"
Explantationsmeißel
Explantation chisel

1/2

26-265-07-07
22 cm/ 8 3/4"
Extraktor
Extractor

26-200-00-04 Set

Implantate

1 x 26-210-01-09	UHP-Standard-Schaft, klein
1 x 26-210-03-09	UHP-Standard-Schaft, mittel
1 x 26-210-05-09	UHP-Standard-Schaft, groß
1 x 26-210-11-09	UHP-Standard+-Schaft, klein
1 x 26-210-13-09	UHP-Standard+-Schaft, mittel
1 x 26-210-15-09	UHP-Standard+-Schaft, groß
1 x 26-210-21-09	UHP-Revisions-Schaft, klein
1 x 26-210-23-09	UHP-Revisions-Schaft, mittel
1 x 26-210-25-09	UHP-Revisions-Schaft, groß
1 x 26-220-01-04	UHP-Prothesenkopf Keramik, klein
1 x 26-220-03-04	UHP-Prothesenkopf Keramik, mittel
1 x 26-220-05-04	UHP-Prothesenkopf Keramik, groß

Probe-Implantate

1 x 26-211-01-05	Probe-Standard-Schaft, klein
1 x 26-211-03-05	Probe-Standard-Schaft, mittel
1 x 26-211-05-05	Probe-Standard-Schaft, groß
1 x 26-211-11-05	Probe-Standard+-Schaft, klein
1 x 26-211-13-05	Probe-Standard+-Schaft, mittel
1 x 26-211-15-05	Probe-Standard+-Schaft, groß
1 x 26-211-21-05	Probe-Revisions-Schaft, klein
1 x 26-211-23-05	Probe-Revisions-Schaft, mittel
1 x 26-211-25-05	Probe-Revisions-Schaft, groß
1 x 26-221-01-05	Probe-Prothesenkopf Metall, klein
1 x 26-221-03-05	Probe-Prothesenkopf Metall, mittel
1 x 26-221-05-05	Probe-Prothesenkopf Metall, groß

Instrumente

1 x 22-130-14-07	Knochen-Pfriegel
1 x 23-264-19-07	Herbert-Hammer
1 x 23-730-17-07	Knochen-Haltezange
1 x 26-265-13-07	Knochen-Haltezange
1 x 26-250-01-07	Standard-Raspel, klein
1 x 26-250-03-07	Standard-Raspel, mittel
1 x 26-250-05-07	Standard-Raspel, groß
1 x 26-260-01-07	Revisions-Raspel, klein
1 x 26-260-03-07	Revisions-Raspel, mittel
1 x 26-260-05-07	Revisions-Raspel, groß
1 x 26-265-05-07	Resektionsmessinstrument
1 x 26-265-07-07	Extraktor
1 x 26-265-15-07	Scharfer Löffel
1 x 26-265-20-07	Weichteilretraktor, rechts
1 x 26-265-22-07	Weichteilretraktor, links
1 x 26-265-30-07	Einschlaginstrument
1 x 26-265-35-07	Explantationmeißel, klein
1 x 26-265-37-07	Explantationmeißel, groß

Lagerung

1 x 55-443-22-04	Codierschild mit Beschriftung, ohne Loch
1 x 55-443-13-04	Farbreiter, blau
1 x 55-442-13-04	MicroStop®-Container 140 x 300 x 600 mm
1 x 55-910-16-01	Revisions-Implantate-Rack
1 x 55-910-17-01	Standard-Implantate-Rack
1 x 55-910-18-01	Instrumenten-Rack
1 x 90-666-52-21	Röntgen-Schablone Standard
1 x 90-667-52-21	Röntgen-Schablone Standard+
1 x 90-668-52-21	Röntgen-Schablone Revision

26-200-00-04 Set

Implants

1 x 26-210-01-09	UHP standard stem, small
1 x 26-210-03-09	UHP standard stem, medium
1 x 26-210-05-09	UHP standard stem, large
1 x 26-210-11-09	UHP standard+ stem, small
1 x 26-210-13-09	UHP standard+ stem, medium
1 x 26-210-15-09	UHP standard+ stem, large
1 x 26-210-21-09	UHP revision stem, small
1 x 26-210-23-09	UHP revision stem, medium
1 x 26-210-25-09	UHP revision stem, large
1 x 26-220-01-04	UHP ceramic head, small
1 x 26-220-03-04	UHP ceramic head, medium
1 x 26-220-05-04	UHP ceramic head, large

Trial Implants


1 x 26-211-01-05	Trial standard stem, small
1 x 26-211-03-05	Trial standard stem, medium
1 x 26-211-05-05	Trial standard stem, large
1 x 26-211-11-05	Trial standard+ stem, small
1 x 26-211-13-05	Trial standard+ stem, medium
1 x 26-211-15-05	Trial standard+ stem, large
1 x 26-211-21-05	Trial revision stem, small
1 x 26-211-23-05	Trial revision stem, medium
1 x 26-211-25-05	Trial revision stem, large
1 x 26-221-01-05	Trial metallic head, small
1 x 26-221-03-05	Trial metallic head, medium
1 x 26-221-05-05	Trial metallic head, large

Instruments

1 x 22-130-14-07	Sharp awl
1 x 23-264-19-07	Herbert hammer
1 x 23-730-17-07	Bone holding forceps
1 x 26-265-13-07	Bone holding forceps
1 x 26-250-01-07	Standard reamer, small
1 x 26-250-03-07	Standard reamer, medium
1 x 26-250-05-07	Standard reamer, large
1 x 26-260-01-07	Revision reamer, small
1 x 26-260-03-07	Revision reamer, medium
1 x 26-260-05-07	Revision reamer, large
1 x 26-265-05-07	Resection guide
1 x 26-265-07-07	Extractor
1 x 26-265-15-07	Sharp spoon
1 x 26-265-20-07	Soft tissue retractor, right
1 x 26-265-22-07	Soft tissue retractor, left
1 x 26-265-30-07	Cone impactor
1 x 26-265-35-07	Explantation chisel, small
1 x 26-265-37-07	Explantation chisel, large

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 UHP-Endoprothesenpass (90-226-59-10)

Sterilisation

The implants are delivered sterile and singly packed. Resterilisation of the implants after the expiry date has run out is possible through the manufacturer, as long as the packaging is still sealed.

Used materials:

Protheses	Titanium
Coating	Titanium
Ceramic heads	Zirkonia ceramic
Cobalt-Chrome heads	CoCr ISO 5832/4

 UHP CD-ROM (90-738-31)

 UHP Endoprosthesis passport (90-226-59-10)

KLS Martin Group

KLS Martin France SARL
68200 Mulhouse · France
Tel. +33 3 89 51 31 50
france@klsmartin.com

KLS Martin UK Ltd.
Reading RG1 3EU · United Kingdom
Tel. +44 1189 000 570
uk@klsmartin.com

KLS Martin do Brasil Ltda.
CEP 04.531-011 São Paulo · Brazil
Tel.: +55 11 3554 2299
brazil@klsmartin.com

Gebrüder Martin GmbH & Co. KG
Representative Office
121471 Moscow · Russia
Tel. +7 499 792-76-19
russia@klsmartin.com

Martin Italia S.r.l.
20864 Agrate Brianza (MB) · Italy
Tel. +39 039 605 67 31
italia@klsmartin.com

Nippon Martin K.K.
Tokyo 113-0024 · Japan
Tel. +81 3 3814 1431
nippon@klsmartin.com

KLS Martin Australia Pty Limited
Artarmon NSW 2064 · Australia
Tel.: +61 2 9439 5316
australia@klsmartin.com

Gebrüder Martin GmbH & Co. KG
Representative Office
201203 Shanghai · China
Tel. +86 21 5820 6251
china@klsmartin.com

Martin Nederland/Marned B.V.
1271 AG Huizen · The Netherlands
Tel. +31 35 523 45 38
nederland@klsmartin.com

KLS Martin L.P.
Jacksonville, FL 32246 · USA
Tel. +1 904 641 77 46
usa@klsmartin.com

KLS Martin SE Asia Sdn. Bhd.
14100 Simpang Ampat, Penang · Malaysia
Tel.: +604 505 7838
malaysia@klsmartin.com

Gebrüder Martin GmbH & Co. KG
Representative Office
Dubai · United Arab Emirates
Tel. +971 4 454 16 55
middleeast@klsmartin.com

Gebrüder Martin GmbH & Co. KG
A company of the KLS Martin Group
KLS Martin Platz 1 · 78532 Tuttlingen · Germany
Postfach 60 · 78501 Tuttlingen · Germany
Tel. +49 7461 706-0 · Fax +49 7461 706-193
info@klsmartin.com · www.klsmartin.com