Nervstimulator NervScout

Nerve stimulator NervScout

Estimulador de nervios NervScout

Neurostimulateur NervScout

Neurostimolatore NervScout

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Instructions for Use

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Mode d'emploi

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CE 0297

DE-EN-ES-FR-IT

V. 1.2- (2013-04)

REF 90-160-58-12

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Valid for hardware status HW 01

KLS martin
GROUP

Deutsch.....	4
English	34
Español	63
Français	93
Italiano	124


Symbolerklärung




Symbol explanation

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Légende

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	Die Pfeile zeigen nach oben, Sendung aufrecht lagern und transportieren Arrows pointing up, transport and store upright Las flechas indican hacia arriba, almacenar y transportar en posición vertical Les flèches indiquent vers le haut, stocker et transporter le colis à la verticale Trasportare e immagazzinare il collo in posizione verticale, con le frecce rivolte verso l'alto
	Zerbrechliches Packgut, Druck und Erschütterung vermeiden Fragile, avoid pressure and shock Mercancía frágil, rogamos evitar compresiones y sacudidas Fragile, éviter les chocs et les secousses Prodotto imballato fragile, evitare pressioni e scosse
	Trocken aufbewahren Store in a dry place Almacenar en un lugar seco Conserver à l'abri de l'humidité Conservare in un luogo asciutto

	Umgebungstemperatur für Lagerung und Transport Ambient temperature for storage and transportation Temperatura ambiente para el almacenamiento y el transporte Température ambiante pour stockage et transport Temperatura ambiente per magazzino e trasporto
	Angaben zur minimalen und maximalen Luftfeuchtigkeit für Lagerung und Transport Information on minimum and maximum atmospheric humidity for storage and transportation Indicaciones sobre la humedad atmosférica mínima y máxima para el almacenaje y el transporte Chiffres indiquant l'humidité relative de l'air minimale et maximale pour les stockage et transport Indicazioni relative all'umidità dell'aria minima e massima per il magazzino e il trasporto
	Angaben zum minimalen und maximalen Luftdruck für Lagerung und Transport Information on minimum and maximum atmospheric pressure for storage and transportation Indicaciones sobre la presión atmosférica mínima y máxima para el almacenaje y el transporte Chiffres indiquant la pression atmosphérique minimale et maximale pour les stockage et transport Indicazioni relative alla pressione atmosferica minima e massima per il magazzino e il trasporto
	CE-Konformitätskennzeichnung CE mark of conformity Marca CE de conformidad Marquage CE Contrassegno CE

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Contents

1	Product Liability and Warranty	36
1.1	General Information	36
1.2	Scope of Delivery	36
1.3	Warranty	37
1.4	User's Inspection	37
1.5	Hotline.....	38
2	Notices concerning this Document	38
2.1	Validity of this Document	38
2.2	Symbols Used in this Document	39
3	Functional Principle and Intended Use	40
4	Putting the Unit into Service	41
4.1	Functions of Controls, Connections and Indicators.....	41
4.1.1	Front View	41
4.1.2	Rear View.....	42
4.2	Setup with maxium® Electrosurgical Unit.....	43
4.3	Mains Connection	44
4.4	Turning the Unit ON/OFF.....	44
4.5	Description of Controls & Optical Indicators	45
4.5.1	Stimulation Current Amplitude	45
4.5.2	Pulse Mode Selection	45
4.5.3	Activation Indicator	45
4.5.4	Error Message.....	45
4.6	Stimulation Current Mode	46
4.6.1	Stand-Alone Operation.....	46
4.6.2	Using the NervScout in Conjunction with an Electrosurgical Unit.....	46
4.6.3	Operation with NeuroPen Handle	49
5	Safety Measures.....	50
6	Safety Checks.....	52
7	Consumables and Accessories	53
8	Requirements for Third-Party Accessories.....	54
9	Care Instructions.....	55
9.1	Cleaning and Disinfecting	55
9.2	Sterilization of Reusable Accessories	56
9.3	Non-Sterilizable Accessories	56

10	Technical Data	57
11	Guidelines and Manufacturer's Declaration Regarding Electromagnetic Compatibility (EMC)	58
12	Ecological Information	62
12.1	Disposal of Packing	62
12.2	Ecological Aspects of Operation.....	62
12.3	Disposal of the Unit	62

1 Product Liability and Warranty

1.1 General Information

We thank you for having decided to buy a KLS Martin product. This product carries the CE mark, which means that it satisfies the essential requirements laid down in the EC Directive concerning medical devices.

We are the manufacturer of this product:

Gebrüder Martin GmbH & Co. KG

A company of the KLS Martin Group

KLS Martin Platz 1 · D-78532 Tuttlingen · Germany

Postfach 60 · D-78501 Tuttlingen · Germany

Tel. +49 7461 706-0 · Fax +49 7461 706-193

info@klsmartin.com · www.klsmartin.com

1.2 Scope of Delivery

- NervScout, item no. 80-065-00-04, 100–240 V is delivered with the following components:
 - Mains cable, country-specific
 - Instructions for Use

Optional components (to be ordered separately):

- Electrode handle with 1 pushbutton for electrodes with shaft Ø 4 mm
- NeuroPen handle with 3 pushbuttons for electrodes with shaft Ø 4 mm
- Electrode handle without pushbutton for electrodes with shaft Ø 4 mm
- Stimulation electrode (silver-coated), Ø 1 mm
- Stimulation electrode (silver-coated), Ø 2 mm
- Stimulation electrode (silver-coated), Ø 3 mm
- Connecting cable for disposable neutral electrodes
- Disposable neutral electrodes
- Single-pedal foot switch for activation via pedal
- Mounting set for use in conjunction with a maxium[®] electrosurgical unit

Connects the nerve stimulator to a KLS Martin electrosurgical unit

- MABS connecting cable, foot switch signal line, 1.2 m
- NervScout neutral electrode connecting cable, "I", (1.2 m)
- NervScout neutral electrode connecting cable, "M", (1.2 m)

1.3 Warranty

Our Standard Terms and Conditions of Sale effective at the time shall apply. Agreements diverging from these Standard Terms and Conditions do not restrict the legal rights of the buyer.

Any warranty exceeding the above provisions shall require a contractual form and shall exclude component-related vandalism, software updates and consumables.

Important Notices

The product may only be repaired by Gebrüder Martin or a qualified person or firm expressly authorized by Gebrüder Martin to perform such work.

If the repair is carried out by a person or firm specially authorized by Gebrüder Martin, the operator of the product is required to obtain from the repairer a certificate with details about the nature and scope of the repair work done. This certificate must show the date of the repair and the details of the person or firm carrying out the work and must be signed.

In all cases where a party other than the product manufacturer performed the work, repaired products must be additionally marked with the repairer's ID label.

Improper interventions or alterations performed by third parties during the period of limitation shall void any and all warranty claims. Unauthorized actions performed on the product shall invalidate any liability claims against Gebrüder Martin.

1.4 User's Inspection

Immediately upon receipt, the goods must be checked for completeness and potential damage in transit. Notice of any such damage must be given immediately.

1.5 Hotline

- Should you have any questions on how to handle the product or use it for clinical applications, please do not hesitate to contact the Product Management:
Tel: +49 7461 706-243
Fax: +49 7461 706-190
 - Should you have any technical questions, please do not hesitate to turn to our Martin Service Center:
Tel: +49 7461 706-343
Fax: +49 7461 706-408
E-mail: service@klsmartin.com
 - Should you have any questions concerning maintenance contracts or training courses, please contact our Technical Service Manager:
Tel: +49 7461 706-332
E-mail: service@klsmartin.com
-

NOTICE

To answer your technical questions as efficiently as possible, our service technicians require the serial number of the product. Therefore, please have this number at hand when contacting our hotline. It is part of the information provided on the rating plate; see section 4.1.2 "Rear View", page 42.

2 Notices concerning this Document

Non-observance of this document can lead to serious or even lethal patient injury!

Improper handling and care as well as non-intended use can lead to premature wear and / or pose a risk to patients and users!

Be sure to read, understand and follow the instructions given below!

- Every user is required to read this document completely and follow them carefully.
- In particular, be sure to heed all cautions, warnings and danger notices.
- Keep this document accessible to users at all times.
- If the NervScout is operated in conjunction with an electrosurgical unit, the Instructions for Use supplied with the electrosurgical unit must be observed as well.

2.1 Validity of this Document

This document is valid for hardware status HW 01. If older hardware versions are used, it may be that functions described in this document are not available.

2.2 Symbols Used in this Document

Throughout this document, important information (such as general or safety-related notices) is marked with the following symbols and signal words:

WARNING

Danger of death or serious injury!

Indicates a situation which, if not avoided, **could** result in death or serious injury!

CAUTION

Danger of minor injury!

Indicates a situation which, if not avoided, could result in minor or moderate injury!

NOTICE

Risk of material damage!

Indicates a situation which, if not avoided, could lead to material damage (loss of time, data loss, device failure, etc.)!

3 Functional Principle and Intended Use

The NervScout nerve stimulator is used for intraoperative differentiation of nerves and muscles in neuromuscular tissue structures that can be stimulated electrically.

Via an active (positive) electrode, negatively polarized electric current pulses (cathodic stimulation) are periodically delivered into the tissue. The current then takes a divergent path from the point of contact to the neutral (negative) electrode, with current density decreasing in the tissue as the distance from the active electrode increases. In a first approximation, this can be defined as an “inverse square law relationship”: as electric current passes through tissue, it should be expected that the current density decreases by a factor of 4 – down to just 25% – when the distance increases by a factor of 2.

Nerve pathways can be identified by assessing their transmission characteristics. This is done by comparing the observed muscle responses with expected maximum contractions. Using a process of “iterative spatial approximation” (i.e., approaching a nerve repeatedly with the stimulation electrode using a lower stimulation current intensity each time while checking for muscle response), it is possible to localize nerves and minimize damage to sensitive structures during operations. The NervScout can thus be used for intraoperative nerve function diagnosis and assessment and for nerve localization to prevent unintended nerve damage to the greatest possible degree.

Cathodic stimulation means using negative rectangular pulses that are delivered into the tissue by electrodes with a diameter of 1, 2 or 3 mm in monopolar application. Activation is possible by finger switch or foot switch. The NervScout can be used as a stand-alone device or in conjunction with an electrosurgical unit (combination mode).


NOTICE

Please note that the KLS Martin NervScout is a pure nerve stimulator. It is not possible to use it for neuromonitoring.

4 Putting the Unit into Service

4.1 Functions of Controls, Connections and Indicators

4.1.1 Front View

1 Connecting socket for electrode handles and NeuroPen handle

2 Connecting socket for neutral electrode

3 Optical activation indicator (white)

4 Power switch (ON/OFF)

5 Dosing error and malfunction indicator (red)

6 Pulse mode selector button (toggle function)

7 Pulse mode indicator (0-3)

8 Change buttons for stimulation current amplitude preselection
(1-35 mA, selectable via UP/DOWN buttons)

9 Stimulation current amplitude indicator

Floating output (isolated from ground) type CF; the device is defibrillation-proof

Observe Instructions for Use!

4.1.2 Rear View


- 1 Neutral electrode output socket for NE connecting cable, for connection to a KLS Martin electrosurgical unit used in combination mode
- 2 Foot switch output socket for foot-switch signal cable item no. 80-181-51-04, for connection to a KLS Martin electrosurgical unit used in combination mode
- 3 Foot switch input (e.g. for item no. 80-811-50-04, 80-811-04-04, 80-824-01-04)
- 4 Foot switch input (e.g. for item no. 80-821-02-04, 80-821-03-04, 80-811-30-04, 80-811-32-04)
- 5 Rating plate (ID label)
- 6 Equipotential bonding contact acc. to IEC 60601-1
- 7 Power input socket (IEC C14 inlet connector) for line voltages ranging between 100 and 240 VAC


Rating plate of KLS Martin nerve stimulator NervScout

4.2 Setup with maxium® Electrosurgical Unit

Usually, the two system components (maxium® electrostimulation unit and NervScout nerve stimulator) are stacked, placing one on top of the other and connecting them with ball pins (1) and two interlocking rails (2). This connection system is available as an optional accessory for the nerve stimulator and therefore not part of the scope of delivery of the mounting set (item no. 08-474-00-10). Before starting the mounting process, the two ball pins must be screwed into the top of the NervScout (3). The rail must be fixed to the underside of the maxium® (4), making sure that the wedge-type (tapered) side (5) is pointing backwards. The counterpart can be mounted in either direction on the housing top of the NervScout.

If you want to install the two units on a pendant system or on an inclined surface, we recommend securing them not only with ball pins, but using connection rails (item no. 08-022-00-18) in addition. In this case, you may need to tap holes for mounting these fixation elements.

If you prefer to install the total system (NervScout, maxium® electrostimulation unit, accessories) as a closed unit, we recommend using the maxium® Cart II.

In this setup, the maxium® electrostimulation unit is fixed in place on top of the NervScout nerve tester and subsequently connected to the maxium® Cart / maxium® Cart II using the connection elements specified above. The KLS Martin maxium® Beamer can also be integrated into this combination system. In this case, the order of the setup is as follows:

- maxium® Cart II
- NervScout
- Beamer
- maxium®

A combination system including the KLS Martin marVac® smoke evacuator is set up as follows:

- maxium® Cart II
- NervScout
- marVac®
- maxium®

4.3 Mains Connection

Prior to switching on the unit for the first time, verify that your supply system matches the voltage specified on the rating plate (located next to the power input socket). If this is not the case, please contact our Martin Service Center.

The unit incorporates a so-called universal PSU allowing operation on voltages ranging from 100–240 V, 50/60 Hz. It may only be operated with power supply systems providing protective grounding.

4.4 Turning the Unit ON / OFF


Once the nerve stimulator has been properly prepared and installed, it can be turned on with the ON/OFF (I/O) switch.

For safety reasons, the NervScout always starts with the “zero setting”: no stimulation current pulse preset (digit “0” displayed), and stimulation current intensity set to the minimum of 1 mA. To make the unit ready for operation, you first need to select a pulse value and set the desired stimulation current intensity.

4.5 Description of Controls & Optical Indicators

4.5.1 Stimulation Current Amplitude

The change buttons (up/down) allow you to preset the stimulation current intensity across a range of 1–35 mA. The selected setting is then shown on the two-digit 7-segment display. Each time the unit is turned on, the last-used value is automatically reset to 1 mA.

4.5.2 Pulse Mode Selection


The desired pulse length is selected via the pulse mode selector button. Note that once the initial self-test has been completed after turning on the unit, the NervScout is still not in operating condition, as indicated by the digit "0" shown on the display. To set the desired pulse mode, operate the toggle button repeatedly to move from "0" to "1", "2" or "3" as required. Thereafter, the cycle starts again with "0".

4.5.3 Activation Indicator

Immediately to the right of the connecting socket for the NervScout handle, the unit has an activation indicator that lights up white during stimulation current activation.

4.5.4 Error Message

If an error or malfunction occurs, e.g. in the event of misdosage due to an internal fault, the error indicator lights up red continuously and stimulation current output is stopped at once. If activation is continued in this situation, the activation alarm tone is sounded continuously.

4.6 Stimulation Current Mode

4.6.1 Stand-Alone Operation

Making the unit ready for operation

- Turn on the unit with the power switch.
- The unit comes up with the following settings: stimulation current pulse form “0”, stimulation current amplitude “1 mA”. It is still not ready for use at this point.
- Connect the neutral electrode to be used on the patient (the NE connecting cable of KLS Martin electrosurgical units (“I” version) is compatible).
- Attach the NE to the patient’s skin.
- Use the pulse mode selector button to select the stimulation current pulse form of your choice:
 - “1” 50 ms / 550 ms pulse / pause (e.g. for TUR)
 - “2” 330 μ s / 550 ms pulse / pause (e.g. for OMF applications, open surgery)
 - “3” 250 μ s / 9 ms pulse / pause (for tetanic stimulation / permanent contraction)
- Connect sterile NeuroPen to NervScout.
- Insert sterile stimulation electrode (ball electrode) into NeuroPen.

Activation

Upon activation, the white activation indicator lights up and an acoustic signal is heard continuously. As soon as stimulation current is flowing through the tissue (closed patient circuit), the continuous alarm tone is switched to pulsing signal tones emitted in sync with the current pulses actually delivered into the tissue. The pitch of the acoustic signal correlates with the preselected stimulation current amplitude. Potential neuromuscular responses will take place in sync with the acoustic signal.


4.6.2 Using the NervScout in Conjunction with an Electrosurgical Unit

Making the unit ready for operation

Apart from the points specified in section 4.6.1, this mode of operation requires additional connections between the electrosurgical unit and the NervScout for the neutral electrode and, if applicable, for foot switch operation. As regards the neutral electrode, the NE connecting cable is used to connect the NE output socket on the rear panel of the NervScout to the NE input socket of the electrosurgical unit. If a foot switch is required, the foot switch signal is passed on in the same way from the foot-switch output socket on the rear panel of the NervScout to the foot-switch input socket of the electrosurgical unit using the MABS connecting cable.


When using a double pedal foot switch in stimulation current mode, activation is possible with either pedal.

Use in stand-alone operation

Use in combination mode

Active instrument: Foot switch activates HF current

⚠ CAUTION

In combination mode, monopolar HF generator operation is possible only if the stimulation-current patient cable is disconnected from the NervScout.

The KLS Martin NervScout connecting socket features a “Connected” identification function. As long as no stimulation patient cable has been plugged into the NervScout, the neutral electrode is properly connected from the patient to the electrosurgical unit via the nerve stimulator and the NE connecting cable linking the two units. There are no restrictions to electrosurgery under this condition. However, as soon as a patient cable is connected to the NervScout, the neutral electrode connection to the electrosurgical unit is automatically interrupted and the neutral electrode is exclusively assigned to the nerve stimulator in this state. As a result, only the bipolar mode is available for electrosurgical operation.

So when using the NervScout in combination mode, the patient cable must be replugged from the electrosurgical unit to the NervScout and vice versa if a single neutral electrode is used jointly for stimulation and electrosurgery.

In this way, the use of a joint handle or patient cable for stimulation and electrosurgery is therefore possible.

NOTICE

When using the NervScout jointly with the KLS Martin maxium® electrosurgical unit, NervScout parameter setting must nonetheless be done separately. Parametrization via the maxium® unit is not possible.

4.6.3 Operation with NeuroPen Handle

The NeuroPen handle allows you to adjust the stimulation current amplitude directly on the handle. Using the “+” (2) and “-” (3) change buttons, you can adjust across a range of 1-35 mA.

Button (1) is used for activating the device.

- 1 Activation button
- 2 Change button “+”
- 3 Change button “-”

NOTICE

Please note that, unlike the electrode handles with one pushbutton and no buttons, the NeuroPen handle cannot be used for electrosurgery.

Please also observe the processing instructions for accessories, available at the KLS Martin website: <http://kls-martin.com/phnet/>

5 Safety Measures

CAUTION

Danger of injury!

The use of non-functional, damaged or defective accessories can pose a danger to patients or users and compromise the specified functions of the nerve stimulator.

- The use of non-approved accessories from other manufacturers can also pose a hazard. If in doubt, please contact Gebrüder Martin.
 - Accessories marked for “single use only” (such as disposable neutral electrodes) must neither be reprocessed nor reused.
 - Accessories unfit for use must be withdrawn from service.
-

CAUTION

General safety instructions

- The NervScout may only be used by persons that have been properly trained in operating the device!
 - Be sure to read the Instructions for Use prior to using the product for the first time.
 - The NervScout may be used only within the scope of its intended purpose.
 - Stimulation currents $> 2 \text{ mAeff/cm}^2$ must be used with special care.
 - Do not apply the neutral electrode and the active electrode contralaterally in the head and thoracic areas.
 - Stimulation can cause unexpected muscle responses. To minimize the risk of mechanical injury, keep sharp-edged objects away from the patient.
 - Do not use stimulation currents with an intensity higher than necessary.
 - When using the NervScout jointly with an electrosurgical unit, make sure that the application rules governing HF surgery are duly observed and a minimum distance of approx. 15 cm is maintained between the neutral electrode of the nerve stimulator and the active electrode of the electrosurgical unit.
 - Care must be taken to prevent any accidental contact between the stimulation electrode and other conductive parts (probes, electrodes, other metal parts).
-

CAUTION
--

Contraindications

- The nerve stimulator must not be used in the following cases:
 - Patients with cardiac pacemakers *
 - Patients with implanted cardioverter-defibrillator (ICD) *
 - Patients with arrhythmia
 - Pregnant patients
 - Epileptic patients
 - NervScout must not be used for cardiac stimulation!

* only under the supervision of a cardiologist

NOTICE

Please note that the KLS Martin NervScout is a pure nerve stimulator. Neuromonitoring is not possible with this device.

6 Safety Checks

At least every 24 months, the following checks and tests must be performed on the unit by qualified persons who have the necessary training, knowledge and practical experience to carry out such safety-related work properly, responsibly and independently.

- Visual inspection: check the unit and its accessories for mechanical and functional damage/defects.
- Check all safety-related labels for legibility.
- Check the fuses for compliance with rated current and prearcing time/current characteristic.
- Perform functional test in accordance with the Instructions for Use.
- Check that the stimulation current amplitude can be adjusted continuously across the range of 1–35 mA.
- Check optical indicator for correct functioning during activation.
- Check that continuous tone is sounded during activation in no-load condition. Pitch continuously increasing between 5 and 25 mA depending on preset current intensity.
- Check that continuous tone is switched to pulsing tone from a current intensity of 1 mA, up to a load resistance $\leq 2.7 \text{ k}\Omega$.
 - Check pulse lengths and pulse pauses using an oscilloscope
- Carry out electrical checks in accordance with the test report sheet for periodic safety checks.
 - Leakage currents must not exceed 1.5 times the value measured initially and at the same time not exceed the specified limit value, either.

The values measured first are available from the test report attached, which was established when the unit was first installed.

We recommend the user to record all safety checks and related results in a medical device logbook.

If the unit is not fully reliable and/or safe to operate, it must be withdrawn from service and marked as “not ready for operation”. A technical check is always required in such cases.

NOTICE

The above checks may be carried out only by qualified Gebrüder Martin service personnel or by a service representative that has been specially authorized to carry out such work.

7 Consumables and Accessories

To ensure reliable stimulation, KLS Martin offers the following accessories/consumables:

Description	Item number
Electrode handle with 1 pushbutton for electrodes with shaft Ø 4 mm	80-217-01-04
NeuroPen handle with 3 pushbuttons for electrodes with shaft Ø 4 mm	80-217-03-04
Electrode handle without pushbutton for electrodes with shaft Ø 4 mm	80-220-00-04
Stimulation electrode silver-coated, Ø 1 mm	80-565-01-04
Stimulation electrode silver-coated, Ø 2 mm	80-565-02-04
Stimulation electrode silver-coated, Ø 3 mm	80-565-03-04
MABS connecting cable, foot switch signal line, 1.2 m *	80-181-51-04
NervScout neutral electrode connecting cable, "I", (1.2 m) *	80-065-01-04
NervScout neutral electrode connecting cable, "M", (1. m) *	80-065-02-04
Connecting cable for disposable neutral electrodes	80-294-44-04
Disposable neutral electrode	80-344-09-04
Single-pedal foot switch for activation via pedal	80-811-50-04
Mounting set for use in conjunction with a maxium® electrosurgical unit	08-474-00-10

* Connects the nerve stimulator to a KLS Martin electrosurgical unit

8 Requirements for Third-Party Accessories

NOTICE

The unit may only be used in combination with accessories, parts subject to wear and disposables that have been evidenced as safe for use by way of a Declaration of Conformity.

For accessories certified for use with the KLS Martin nerve stimulator, please refer to the KLS Martin Accessories Catalog, which is available from Gebrüder Martin or can be directly downloaded at www.klsmartin.com.

⚠ CAUTION**Risk of malfunction and danger of injury!**

Improper handling by the user can lead to failure of the device and cause injury! Therefore, please heed the following:

- All accessories used must be suitable for the voltages specified in section 10 "Technical Data", page 57.
 - Use only accessories that are in perfect working order!
-

9 Care Instructions

NOTICE

Be sure to observe all local hygiene and processing regulations.

9.1 Cleaning and Disinfecting

The responsibility for cleaning and disinfection lies with the operator/user of the product. Be sure to observe your local cleaning and hygiene regulations!

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⚠ WARNING

Danger of serious injury when using defective accessories!

Using a unit that is not in perfect working order or using defective accessories can pose a hazard to patients or users and compromise the proper functioning of the device.

- Be sure to maintain the unit and its accessories in perfect operating condition. Perform a visual inspection and check the system for insulation damage, cleanliness and integrity.
 - A unit or accessories found to be unfit for use must be withdrawn from service.
-
- Never sterilize the unit!
 - Always disconnect the unit from the supply before cleaning and disinfecting it.
 - No liquid must be allowed to enter the unit. Therefore clean and disinfect it with care.
Do not use sprays!
 - Never clean the unit with abrasive cleaners, disinfectants or solvents that could scratch the housing or damage the device in other ways.
 - For surface cleaning and surface disinfection, use the procedure recommended by your hospital or use another nationally recognized and approved method.
 - When disinfecting accessories (by surface disinfection or immersion), be sure to observe the instructions of the manufacturer of the disinfectant regarding material compatibility, dosage and exposure time.
 - Accessories not suitable for sterilization, such as the foot switch, must be cleaned with a disinfecting detergent.
 - All disinfectant residues must be completely removed before putting the unit in operation.
 - If accessories must be disinfected or sterilized, be sure to observe the respective Instructions for Use.

9.2 Sterilization of Reusable Accessories

The handles and the stimulation electrodes must be processed in accordance with the instructions provided in the respective Instructions for Use.

(Steam sterilization at 134 °C [273.2 °F] / 20 min).

Please also observe the processing instructions for accessories, available at the KLS Martin website: <http://kls-martin.com/phnet/>

<i>NOTICE</i>

The responsibility for proper cleaning, disinfection and sterilization lies with the operator or product user. Be sure to observe all local regulations (including potential restrictions)!

9.3 Non-Sterilizable Accessories

- Accessories such as the foot switch must be cleaned and disinfected.

10 Technical Data

Description	Nerve stimulator NervScout	
Power supply:	100–240 V $\pm$ 10%; 50/60 Hz	
Power input:	15 VA	
Protection class:	I	
Degree of protection:	IP 21	
Line fuses:	2 x "T 1.25 AH, 250 V"	
Equipotential bonding:	Acc. to IEC 60601-1	
Dimensions:	(W x H x D) 390 x 108 x 425 mm	
Weight:	3.5 kg	
Type of stimulation:	Constant current	
Pulse form:	Negative rectangle	
Pulse length/pause:	Mode "1": 50 ms / 550 ms Mode "2": 330 μ s / 550 ms Mode "3": 250 μ s / 9 ms	
Pulse amplitude:	1–35 mA (continuously adjustable in 1-mA increments)	
Pulse voltage:	max. -60 V	
Classification acc. to 93/42 EEC:	IIb	
Type acc. to EN 60601-1:	Applied part type CF, defibrillation-proof	
Electromagnetic compatibility:	Acc. to EN 55011 and IEC 60601-1-2 Safety requirements acc. to IEC 801	
Environmental conditions for transportation and storage:	Ambient temperature Relative humidity Atmospheric pressure	-25°C to +70°C (-13°F to 158°F) 10% to 100% 500 hPa to 1,060 hPa
Environmental conditions for operation:	Ambient temperature Relative humidity Atmospheric pressure	+10°C to +40°C (+50°F to 104°F) 30% to 75% 700 hPa to 1,060 hPa
 0297	Conforms to Directive 93/42/EEC on medical devices	

11 Guidelines and Manufacturer's Declaration Regarding Electromagnetic Compatibility (EMC)

Guidance and manufacturer's declaration according to IEC 60601-1-2, clause 6.8.3.201 a) 3) Table 201: Electromagnetic emission

This NervScout device is intended for use in an electromagnetic environment as specified below. The user should assure that it is used in such an environment.

Emission test	Compliance level	Electromagnetic environment – Guidance
Harmonic emissions according to IEC 61000-3-2	Class A	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Voltage fluctuations / flicker emissions according to IEC 61000-3-3	Complies	

**Guidance and manufacturer's declaration according to IEC 60601-1-2,
clause 6.8.3.201 a) 3) Table 201: Electromagnetic emission**

The device is intended for use in an electromagnetic environment as specified below. The user should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – Guidance
Electrostatic discharge (ESD) according to IEC 61000-4-2	±6 kV contact discharge ±8 kV air discharge	±6 kV contact discharge ±8 kV air discharge	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Fast electrical transients / bursts according to IEC 61000-4-4	±2 kV for mains power lines ±1 kV for input and output lines	±2 kV for mains power lines ±1 kV for input and output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge according to IEC 61000-4-5	±1 kV differential-mode voltage ±2 kV common-mode voltage	±1 kV differential-mode voltage ±2 kV common-mode voltage	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short power interruptions and supply voltage fluctuations according to IEC 61000-4-11	<5% U_T (>95% dip in U_T) for ½ period 40% U_T (60% dip in U_T) for 5 periods 70% U_T (30% dip in U_T) for 25 periods <5% U_T (>95% dip in U_T) for 5 seconds	0% U_T (100% dip in U_T) for ½ period 40% U_T (60% dip in U_T) for 5 periods 70% U_T (30% dip in U_T) for 25 periods 0% U_T (100% dip in U_T) for 5 seconds	Mains power quality should be that of a typical commercial or hospital environment. If the user requires continued operation during mains power interruptions, it is recommended to power the device from an uninterruptible power supply.
Power frequency (50 / 60 Hz) magnetic field according to IEC 61000-4-8	3 A/m	3 A/m	Power-frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Note: U_T is the a. c. mains voltage prior to application of the test level.

**Guidance and manufacturer's declaration according to IEC 60601-1-2,
 clause 6.8.3.201 b) Table 204: Electromagnetic immunity**

The device is intended for use in an electromagnetic environment as specified below. The user should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – Guidance
Conducted RF according to IEC 61000-4-6 Radiated RF according to IEC 61000-4-3	3 V _{eff} 150 kHz to 80 MHz 3 V/m 80 MHz to 2.5 GHz	3 V _{eff} 150 kHz to 80 MHz 3 V/m 80 MHz to 2.5 GHz	Portable and mobile RF communications equipment should be used no closer to any part of the device (including cables) than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = 1.17\sqrt{P}$ $d = 1.17\sqrt{P} \quad \text{for 80 MHz to 800 MHz}$ $d = 2.33\sqrt{P} \quad \text{for 800 MHz to 2.5 GHz}$ where <i>P</i> is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and <i>d</i> is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^a, should be less than the compliance level in each frequency range.^b Interference may occur in the vicinity of equipment marked with the following symbol. 
Note 1:	At 80 MHz and 800 MHz, the higher frequency range applies.		
Note 2:	These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.		

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio stations and AM and FM radio and TV stations, cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level specified above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.

^b Over the frequency range of 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

**Recommended separation distances between portable and mobile
RF telecommunications equipment and the KLS Martin NervScout
according to IEC 60601-1-2 clause 6.8.3.201 b), Table 206**

The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The user can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF telecommunications equipment (transmitters) and the device as recommended below, depending on the maximum output power of the communications equipment.

Rated maximum output P of transmitter in watts	Separation distance according to frequency of transmitter (in m)		
	150 kHz to 80 MHz $d = 1.17\sqrt{P}$	80 MHz to 800 MHz $d = 1.17\sqrt{P}$	800 MHz to 2.5 GHz $d = 2.33\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.37	0.37	0.74
1	1.17	1.17	2.33
10	3.7	3.7	7.4
100	11.7	11.7	23.3

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1:	At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.
Note 2:	These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

12 Ecological Information

12.1 Disposal of Packing

Gebrüder Martin will, as a matter of course, take back the full packaging if so desired, to recycle as many parts of the packaging as possible.

If you do not wish to make use of this offer, you can dispose of the packaging with the normal paper or domestic garbage.

12.2 Ecological Aspects of Operation


If the treatment is interrupted and the unit is not being used for an extended period, we recommend you to switch it off for safety as well as economic reasons (energy saving).

If disposable products are used in a treatment, please note that such articles must first be carefully cleaned, disinfected and – where applicable – sterilized before they are disposed of as domestic or hazardous waste. Infected sharp parts of single-use products must be handled like any other „sharps“ (cannulas, needles and scalpels) in accordance with valid regulations (disposal via germ-proof and puncture-proof containers).

12.3 Disposal of the Unit

In designing the unit, we tried to avoid using composite materials wherever possible. This allows a high degree of recycling after the lifetime of the unit. We therefore offer to take the unit back for proper disposal and recycling.

Be sure to observe your national/local rules and regulations governing disposal!

	Marking of electric and electronic equipment in accordance with Directive 2002/96/EC (WEEE Directive) and the German Electrical and Electronic Equipment Act (ElektroG) This symbol on the product or its packaging indicates that the product may not be disposed of as normal household garbage.
---	---

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Rudolf Buck GmbH
KLS Martin UK Ltd.
KLS Martin Australia Pty Limited