



English

Electrosurgical Unit ME MB2



Instructions for Use

REF 90-131-52-40 Revision No. 5


Date of Release: 2015-12





Quick Instructions for Use ME MB2

For the assignment of position numbers, see section 4.2, page 12 and section 4.3, page 14.





Symbol explanation

\wedge	Safety alert sym	bol
	CAUTION	Indicates a situation which, if not avoided, could result in minor or moderate injury.
	WARNING	Indicates a situation which, if not avoided, could result
	DANGER	in death or serious injury. Indicates a situation which, if not avoided, will result in death or serious injury.
	Observe instructio	ns for use
REF	Item number (Iter	m no.)
	Manufacturer	
NON WOOD		non-wood", i.e. wood packaging material that is not and therefore suitable for international shipment.
A	Warning: Dangero	us high voltage!
	Class I equipment	acc. to IEC 60601-1
C € 0297	CE mark of confor	mity
	This product may	not be disposed of as normal household garbage



Contents

1 1.1	Product Liability and Warranty	
1.2	Scope of Delivery	
1.3	Intended Use	
1.4	Warranty	
1.5	User's Inspection	
1.6	Hotline	
1.0		
2	Notices concerning this Document	9
2.1	Symbols Used in this Document	9
2.2	Service Manual	9
3	Functional Principle	10
3.1	Monopolar Functional Principle	
3.2	Bipolar Functional Principle	
3.2	Bipolai Functional Finicipie	10
4	Commissioning	11
4.1	Technical Specification	11
4.2	Functions of the Control Elements, Connections and Displays of the ME MB2m (standard KLS Martin version)	12
4.3	Functions of the Control Elements, Connections and Displays of the ME MB2i (international version)	
4.4	Mains Connection	
4.5	Switching On / Off	16
4.6	Self-Test	16
4.7	Monopolar Operation	17
4.7.1	Connection of the Neutral Electrode	17
4.7.2	Connection of Handles	17
4.7.3	Selecting the Current Type	18
4.7.4	Power Setting	
4.7.5	Operation	
4.7.6	Functional Check	
4.8	Bipolar Operation	
4.8.1	Function of the Neutral Electrode	
4.8.2	Automatic Operating Mode	
4.8.3	Foot Switch Operation	
4.8.4	Power Setting	
4.8.5 4.8.6	Operation Functional Check	
5	Safaty Maacuros	26
	Safety Measures	
5.1	General Information	26

Instructions for Use Electrosurgical Unit ME MB2



5.2	Patient Placement	26
5.3	Application of the Neutral Electrode	27
5.4	Working with the Active Electrode	28
5.5	Work on Patients with Pacemakers or Other Implants	28
5.6	Cable Routing on the Patient	28
5.7	Operation with KLS Martin Argon Beamer System	29
5.8	Putting Down HF Instruments	29
5.9	Accessories	29
5.10	Use of Two Electrosurgery Devices on One Patient	
5.11	Inadvertent emission of HF energy	
5.12	Additional Safety Notes	30
6	Cleaning and Disinfection	31
7	Safety Checks	32
8	Accessories	33
9	Technical Data	34
10	Diagrams	36
10.1	Performance Diagrams	
10.2	Voltage Diagrams	
11	Guidelines and Manufacturer's Declaration Regarding	
	Electromagnetic Compatibility (EMC)	38
12	Ecological Information	
12.1	Disposal of Packing	
12.2	Ecological Aspects of Operation	
12.3	Disposal of the Unit	42



1 Product Liability and Warranty

1.1 General Information

We thank you for having decided to buy a KLS Martin product. This product carries the CE mark, which means that it satisfies the essential requirements laid down in the EC Directive concerning medical devices.

We are the manufacturer of this product:

Gebrüder Martin GmbH & Co. KG

A company of the KLS Martin Group KLS Martin Platz 1 · D-78532 Tuttlingen · Germany Postfach 60 · D-78501 Tuttlingen · Germany Tel. +49 7461 706-0 · Fax +49 7461 706-193 info@klsmartin.com · www.klsmartin.com

1.2 Scope of Delivery

HF generator

- Electrosurgical unit ME MB2
- Mains cable
- Instructions for Use

1.3 Intended Use

The unit is used for electrosurgical cutting or coagulation of live human tissue.

The operator may operate the device only if an on-site functional test has previously been performed by Gebrüder Martin or a person authorized by Gebrüder Martin. In addition, a responsible person designated by the operator must have been instructed in the proper handling, application and operation of the unit, as well as in its permissible combination with other medical devices, objects and accessories. This duly instructed officer shall subsequently be responsible for familiarizing the operator's staff with the unit as the need arises.

We recommend documenting all user instructions in a medical device logbook. A copy of the logbook is available from Gebrüder Martin.

The operational safety of the unit must be verified at regular intervals, see section 7 "Safety Checks", page 7.

If the device is not functionally reliable and / or safe to operate, it must be marked as such and withdrawn from service. A technical check is mandatory in any such case.



1.4 Warranty

Our Standard Terms and Conditions of Sale effective at the time shall apply. Agreements diverging from these Standard Terms and Conditions do not restrict the legal rights of the buyer.

Any warranty exceeding the above provisions shall require a contractual form and shall exclude component-related vandalism, software updates and consumables.

Important Notices

The product may only be repaired by Gebrüder Martin or a qualified person or firm expressly authorized by Gebrüder Martin to perform such work.

If the repair is carried out by a person or firm specially authorized by Gebrüder Martin, the operator of the product is required to obtain from the repairer a certificate with details about the nature and scope of the repair work done. This certificate must show the date of the repair and the details of the person or firm carrying out the work and must be signed.

In all cases where a party other than the product manufacturer performed the work, repaired products must be additionally marked with the repairer's ID label.

Improper interventions or alterations performed by third parties during the period of limitation shall void any and all warranty claims. Unauthorized actions performed on the product shall invalidate any liability claims against Gebrüder Martin.

1.5 User's Inspection

Immediately upon receipt, the goods must be checked for completeness and potential damage in transit. Notice of any such damage must be given immediately.



1.6 Hotline

• Should you have any questions on how to handle the product or use it for clinical applications, please do not hesitate to contact the Product Management:

Tel: +49 7461 706-243 Fax: +49 7461 706-190

 Should you have any technical questions, please do not hesitate to turn to our Martin Service Center:

Tel: +49 7461 706-343
Fax: +49 7461 706-408
E-mail: service@klsmartin.com

 Should you have any questions concerning maintenance contracts or training courses, please contact our Technical Service Manager:

Tel: +49 7461 706-332 E-mail: service@klsmartin.com

NOTICE

To answer your technical questions as efficiently as possible, our service technicians require the serial number of the product. Therefore, please have this number at hand when contacting our hotline. It is part of the information provided on the rating plate; see section 4.2 "Functions of the Control Elements, Connections and Displays of the ME MB2m (standard KLS Martin version)", page 12.



2 Notices concerning this Document



Non-observance of this document can lead to serious or even lethal patient injury!

Improper handling and care as well as non-intended use can lead to premature wear and / or pose a risk to patients and users!

Be sure to read, understand and follow the instructions given below!

- Every user is required to read this document completely and follow them carefully.
- In particular, be sure to heed all cautions, warnings and danger notices.
- Keep this document accessible to users at all times.

2.1 Symbols Used in this Document

Throughout this document, important information (such as general or safety-related notices) is marked with the following symbols and signal words:



Danger of death or serious injury!

Indicates a situation which, if not avoided, **could** result in death or serious injury!



Danger of minor injury!

Indicates a situation which, if not avoided, could result in minor or moderate injury!

NOTICE

Risk of material damage!

Indicates a situation which, if not avoided, could lead to material damage (loss of time, data loss, device / machine failure, etc.)!

2.2 Service Manual

Upon request, we will supply a service manual. It helps persons or companies authorized by us to maintain / repair the unit. The service manual also contains order information for all components that are available as spare parts.



3 Functional Principle

3.1 Monopolar Functional Principle

The electrosurgical unit ME MB2 is a generator that converts electrical energy from the grid into a high-frequency current. This high-frequency current is fed to a point-shaped active electrode via a supply line or a handle. This leads to high field line concentration in a small area at the application site of the active electrode in the tissue surroundings of the point of contact. The desired electrosurgical effect is achieved by this energy concentration in a small area in the area of the active electrode. During further energy transport through the patient up to the neutral electrode that is applied to a larger area, the concentration of the current is further reduced. Therefore, there is, in accordance with its intended purpose, no thermal effect in the area of the applied neutral electrode. The circuit is closed via the supply cable of the neutral electrode.

The HF generator can be activated via either the foot switch or the finger switch on the surgical handle.

3.2 Bipolar Functional Principle

Special constructional measures (insulation) permit design of bipolar instruments where the active electrode and neutral electrode are positioned directly to each other. Here the path of the high-frequency current only leads to one instrument tip to the other. Thus, there result very short current paths and surrounding coagulation areas with low power requirements.



4 Commissioning

4.1 Technical Specification


The electrosurgical unit is controlled by a micro-processor, can be used universally and has excellent performance data. The unit meets the latest safety standards.


The unit is characterized by the following constructional characteristics:

- The characteristics of the HF generator result in automatic power adjustment for different operation electrodes, tissue types and cutting speeds over wide ranges.
- Power adjustments for particularly precise settings in the lower power range and high power reserves in the upper adjustment range.
- Advanced safety concept with two micro-processors.
- Easy handling thanks to the clear arrangement of all control elements and their labeling with easy-to-remember function symbols. The power is adjusted via rotary knobs. All connections for active electrodes, neutral electrodes and foot switches are located on the front side.
- Two current types for cutting: Smooth cutting and cutting with eschar formation.
- A pulsed current type for cutting, used e. g. in polypectomy.
- Two current types for coagulation: One coagulation current with high power for contact coagulation and a coagulation current with high crest factor for spray coagulation.
- The unit additionally includes a generator for bipolar coagulation.
- Bipolar application can also be conveniently and safely be realized without a foot switch,
 via the bipolar automatic button.
- Self-test after each switching on of the unit.
- Optical and acoustic signals to indicate HF activation with activation LEDs in different colors and signal tones for cutting and coagulation.
- Combined connection sockets for KLS Martin coax cables and disposable products.
- The integrated Patient Control System (PCS) by KLS Martin performs continuous application monitoring using a two-surface neutral electrode.
- Safety system according to IEC 60601-2-2 for cases of a unit-related output error.
- Acoustic warning in case of high-frequency activation for more than 15 seconds (switch-on time).
- No opening or cooling vents, no cooling fan. Low power loss thanks to high efficiency of the generators.
- Hygiene-friendly foil front plate.
- Practice-oriented, high-quality accessory components.



4.2 Functions of the Control Elements, Connections and Displays of the ME MB2m (standard KLS Martin version)







- 1 Mains switch
- 2 Activation LED monopolar cutting
- **3** Selection button monopolar cutting 1 (pure)
- 4 Selection button monopolar cutting 2 (blend)
- 5 Display monopolar cutting
- 6 Selection button endo mode
- **7** Power controller monopolar cutting
- 8 Activation LED monopolar coagulation
- **9** Selection button contact coagulation
- **10** Selection button spray coagulation
- 11 Display monopolar coagulation
- 12 Power controller monopolar coagulation
- 13 Activation LED bipolar coagulation
- 14 Display bipolar coagulation
- 15 Power controller bipolar coagulation
- 16 Selection button AUTOCOAG bipolar
- 17 Connection socket for bipolar instruments
- 18 Switch button for blue foot switch pedal (bipolar)
- **19** Switch button for blue foot switch pedal (monopolar)
- **20** Connection socket for double-pedal foot switch for monopolar cutting and monopolar coagulation, switchable to bipolar coagulation
- 21 Monopolar connection socket
- 22 Indicator LED for neutral electrode
- 23 Connection socket for neutral electrode
- 24 Speakers
- 25 Connection socket for beamer
- **26** Mains fuses
- 27 Connection socket for mains cable
- 28 Connector for equipotential bonding
- 29 Type plate



Connection for the neutral electrode.

Neutral electrode insulated against ground (floating).



Symbol for classification of the unit (CF). The unit is defibrillation-proof.




WARNING! PLEASE READ INSTRUCTIONS FOR USE!




WARNING! HIGH-FREQUENCY CURRENTS! CAUTION! HIGH VOLTAGE!



4.3 Functions of the Control Elements, Connections and Displays of the ME MB2i (international version)







- 1 Mains switch
- 2 Activation LED monopolar cutting
- **3** Selection button monopolar cutting 1 (pure)
- 4 Selection button monopolar cutting 2 (blend)
- 5 Display monopolar cutting
- 6 Selection button endo mode
- **7** Power controller monopolar cutting
- 8 Activation LED monopolar coagulation
- **9** Selection button contact coagulation
- **10** Selection button spray coagulation
- 11 Display monopolar coagulation
- 12 Power controller monopolar coagulation
- 13 Activation LED bipolar coagulation
- 14 Display bipolar coagulation
- **15** Power controller bipolar coagulation
- 16 Selection button AUTOCOAG bipolar
- 17 Connection socket for bipolar instruments
- 18 Switch button for blue foot switch pedal (bipolar)
- **19** Switch button for blue foot switch pedal (monopolar)
- **20** Connection socket for double-pedal foot switch for monopolar cutting and monopolar coagulation, switchable to bipolar coagulation
- 21 Monopolar connection socket
- 22 Indicator LED for neutral electrode
- 23 Connection socket for neutral electrode
- 24 Speakers
- 25 Connection socket for beamer
- **26** Mains fuses
- 27 Connection socket for mains cable
- 28 Connector for equipotential bonding
- 29 Type plate



Connection for the neutral electrode.

Neutral electrode insulated against ground (floating).



Symbol for classification of the unit (CF). The unit is defibrillation-proof.



WARNING! PLEASE READ INSTRUCTIONS FOR USE!



WARNING! HIGH-FREQUENCY CURRENTS! CAUTION! HIGH VOLTAGE!



4.4 Mains Connection

For the assignment of the numbering, see section 4.2 "Functions of the Control Elements, Connections and Displays of the ME MB2m (standard KLS Martin version)", page 12.

Before initial switching on of the unit, ensure that your grid complies with the voltage setting of the unit indicated on the label (next to the power input). If the voltages should not match, please contact our Martin Service Center.

Any adjustment to another supply voltage on site may only be carried out by a person authorized to do so by Gebrüder Martin. Unauthorized persons may not open the unit!

The unit is available in two versions:

100 - 127 VAC, T 8 A 220 - 240 VAC, T 4 A

Connect the unit to the power outlet using the mains cable. Set mains switch (1) on the front side of the unit to position "ON".

4.5 Switching On / Off

After switching on at the mains switch (1), the unit is set to operating mode and performs a self-test.

4.6 Self-Test

After switching on, the unit performs a self-test. All five current type buttons are briefly lit up consecutively. The following simultaneous lighting up of the five current type buttons, the figures 888, 888 and 188 shown on the power display fields and a short acoustic signal indicate successful completion of the self-test.

The device is operational.

If there are no confirmation signals, the unit is defective and must not be operated. Please contact the Martin Service Center (see section 1.6 "Hotline", page 8).



4.7 Monopolar Operation

For the assignment of the numbering, see section 4.2 "Functions of the Control Elements, Connections and Displays of the ME MB2m (standard KLS Martin version)", page 12.

4.7.1 Connection of the Neutral Electrode

If the neutral electrode is not connected, the red indicator LED of the neutral electrode (22) flashes. If you try to activate the unit in this state via the finger or foot switch, there will be an additional acoustic warning signal. The HF current cannot be activated.

The indicator LED of the neutral electrode (22) goes dark immediately after connecting a onesurface neutral electrode. Therefore ensure that the neutral electrode is safely and fully attached to the patient.

The indicator LED of the neutral electrode (22) goes dark after reaching a safe application state when a multi-surface neutral electrode is connected. As in this case individual run-in times are to be expected, a lead time for the application must be considered.

4.7.2 Connection of Handles

For monopolar cutting and coagulation, a handle with finger switch or a foot switch in combination with a handle without finger switch can be connected.

Handles with finger switch are connected to the monopolar connection socket (21). Handles without finger switch or surgical instruments (e. g. resectoscopes) are also connected to the monopolar connection socket (21), the corresponding foot switch to the connection socket (20).

Handles with coax plug as well as handles with 3-pin plug can be connected to the monopolar connection socket (21). For further information, please see the accessories catalog of KLS Martin.

Surgical HF instruments without switch function, such as resectoscopes, must be connected to the right connection (coax) of the monopolar connection socket (21), as the foot switch acts only on this port.

The desired active electrode must be inserted into the hexagonal guide of the handle until the hexagon of the electrode snaps into place. This prevents rotation of the electrode during the application.



Risk of electric shock!

When inserting the electrode and while replacing an electrode, the high-frequency current must not be activated.



4.7.3 Selecting the Current Type

Three different current types are available in the control field **Monopolar Cut** that is marked **yellow** on the front plate.



Cutting 1 (pure)

Clean cut without eschar formation



Cutting 2 (blend)

Cutting current with low eschar formation



Endo mode

Time-controlled cutting mode that offers fractionated and hence controlled cutting for special applications, e. g. endoscopic papillotomy and endoscopic polypectomy.

This briefly intensive pulse at the beginning of the cutting process with a sufficient amount of coagulation for hemostasis is repeated in constant intervals. Thus, pedunculated growths also known as polyps are removed endoscopically by diathermic loop excision (DLE).

Two different current types are available in the control field **Monopolar Coag.** that is marked **blue** on the front plate.



Contact coagulation

Coagulation with high depth and direct contact between electrode and tissue



Spray coagulation

Coagulation current with high voltage for surface coagulation (fulguration).


Particularly suitable for hemostasis in TUR with small-area instruments such as loops




4.7.4 Power Setting


Use the power controller (7) for setting the cutting power. Use the power controller (12) for setting the coagulation power. The power setting is adjusted from a preset minimum value to a maximum value depending on the selected current type (see section 9 "Technical Data", page 34). The control characteristic line allows a very precise power increase in the lower range, and an evenly increasing power setting in the high power range. This enables precise adjustment of the lower power range for demanding operations without need for special functions, and to provide sufficient power reserves for applications with high power demands.

Please note that due to this advantage the controller setting in the middle of the adjusting range does not correspond to half of the maximum output power.





^{2 =} Cutting 2 (blend)



1 = Contact coagulation

^{3 =} Endo mode

^{2 =} Spray coagulation



4.7.5 Operation



Failure to comply with the patient and user rules for the application of electrosurgery explained in section 5 "Safety Measures", page 26, may result in severe and even lethal injury!

- Ensure safe application of the neutral electrode and correct positioning of the patient.
- After prolonged application of HF current at high power, the surface of the unit can heat up significantly.
- Only accessories in flawless working condition may be used.
- The accessory equipment must be suitable for at least the voltages indicated in section 9 "Technical Data", page 34.
- During simultaneous activation of two control elements (e. g. blue and yellow button on the handle), the generator is blocked, and there will be no HF power.

Depending on the selected application purpose, the HF current is activated by pressing the yellow activation button for the cutting current and the blue activation button for the coagulation current on the handle or foot switch, respectively. The HF current is provided in accordance with the preset current types and the set power. By actuating an activation button on the handle or foot switch, the selected working channel is activated. At the same time, an acoustic activation signal is sounded that is different for cutting and coagulation, and the corresponding activation LED (2) or (8) lights up. In case of HF activation for more than 15 s, the volume of the signal tone is increased.

NOTICE

During monopolar applications with foot switch, the switch button (19) (monopolar) must be activated for the blue foot switch pedal.

If the switch button (18) (bipolar) is activated, the bipolar HF output port is activated by pressing the blue foot switch pedal!

The yellow foot switch pedal always activates the monopolar cutting current. This is independent of the switch buttons (18) and (19).

These activation buttons of a connected handle always activate the monopolar current types! This is independent of the switch buttons (18) and (19).



4.7.6 Functional Check

All functions of the unit must be checked before using the unit.

The following functional checks must be performed:

- Remove plug of the connection cable for the neutral electrode from the connection socket (23).
 - The red indicator LED (22) flashes.
 - Upon attempts to activate the HF current, the acoustic warning signal is sounded.
 - Activation of the HF current is blocked.
- Re-insert the plug of the connection cable for the one-surface neutral electrode into the connection socket (23).
 - The red indicator LED (22) must not flash any longer.
 - In case of a two-part neutral electrode, this must be correctly applied to the patient for the alarm to disappear.
- Connect connection cable with electrode handle to the monopolar connection socket (21).
- Select the individual current types by actuating the finger switch on the handle or the foot switch.
 - The activation LEDs (2) or (8) assigned to the individual current types must light up, and the HF activation signal must sound.

ACAUTION

Risk of malfunction, risk of injury!

If the HF activation signal can be heard when no double foot switch or electrode handle is connected, the unit is defective.

- The unit must not be operated.
- A technical inspection must be performed.

A CAUTION

Risk of malfunction, risk of injury!

If the HF activation signal can be heard while the foot switch or electrode handle is connected but none of the control elements is triggered, one of these accessory parts is defective.

- The defective accessory must not be used.
- It must be replaced.



4.8 Bipolar Operation

For the assignment of the numbering, see section 4.2 "Functions of the Control Elements, Connections and Displays of the ME MB2m (standard KLS Martin version)", page 12.

4.8.1 Function of the Neutral Electrode

For bipolar application no neutral electrode is required. Even though the indicator LED (22) for the neutral electrode flashes, the HF generator can be activated in this mode.

As a changeover from bipolar to monopolar application during the course of an operation is to be expected, we recommend application of a neutral electrode under all circumstances.

4.8.2 Automatic Operating Mode

Function "Bipolar automatic operation":

One non-adjustable coagulation current type is available in the operating field **Bipolar** coagulation that is marked **blue** on the front plate.



When actuated, the selection button **AUTOCOAG** bipolar (16) has the function to automatically activate the HF current upon tissue contact and to deactivate it after the coagulation.

The ME MB2 has an additional advantageous feature:

After the beginning of the automatic bipolar coagulation, the tissue dries up more and more in the **main coagulation phase**. Further drying results in an increase in tissue resistance and therefore in the automatic shutdown of the generator.

Shortly after that, however, the tissue will be ready again for further energy input. The subsequent **pulse-type post-coagulation phase** produces optimum coagulation results.

The post-coagulation phase is ended by opening the forceps.

The post-coagulation phase is also indicated visually by the flashing of the blue activation LED (13). This signal indicates that the main coagulation phase is completed.



Unintentional HF activation!

In the function **Bipolar automatic operation** of the unit, there is the possibility of unintended HF activation.

The bipolar instrument used must never come into unintentional contact with tissue or low-resistance parts.

The switching-on process can also be triggered by contact with metal parts (e. g. trocars) or similarly low-resistance materials (e. g. liquids). Unintentional HF activation may also occur if the bipolar instrument is placed on the patient and thus comes into contact with tissue.



4.8.3 Foot Switch Operation

The function **bipolar coagulation** is available on the unit. The following accessory components can be connected:

All bipolar connection cables with small coax plugs (e. g. for bipolar forceps, bipolar HF pliers) can be connected to the connection socket (17). Bipolar coagulation is activated via the foot switch.

To activate the bipolar coagulation current, the blue foot switch pedal must be switched to the bipolar generator. To do this, press the switch button (18) (bipolar). The switch button (18) flashes. The setting remains stored even when the unit is switched off.

When the switch button (19) (monopolar) is activated, the monopolar HF output port is activated by pressing the blue foot switch pedal!


The yellow foot switch pedal always activates the monopolar cutting current. This is independent of the switch buttons (18) and (19).

• These activation buttons of a connected handle always activate the monopolar current types. This is independent of the switch buttons (18) and (19).

4.8.4 Power Setting

Use the power controller (15) for setting the bipolar coagulation power. The power setting is adjusted from a preset minimum value to a maximum value (see Technical Data). The control characteristic line allows a very precise power increase in the lower range, and an evenly increasing power setting in the high power range. This enables precise adjustment of the lower power range for demanding operations without need for special functions, and to provide sufficient power reserves for applications with high power demands.

Please note that due to this advantage the controller setting in the middle of the adjusting range does not correspond to half the output power.



1 = Bipolar coagulation



4.8.5 Operation

The switch button (18) must be switched to the current type **bipolar coagulation**!

Upon actuation of the blue foot switch pedal, the current type **bipolar coagulation** is activated.

The HF current is provided in accordance with the selected power. At the same time, an acoustic signal can be heard, and the corresponding activation LED (13) lights up. In case of HF activation for more than 15 seconds, the volume of the signal tone is increased.



Risk of malfunction, risk of injury!

After prolonged application of HF current at high power, the surface of the unit can heat up significantly.



4.8.6 Functional Check

The following functional checks must be performed:

- The neutral electrode does not have to be connected for bipolar operation. As described in section 4.8.1 "Function of the Neutral Electrode", page 22, we recommend connecting the neutral electrodes even for bipolar mode.
- Connect connection cable with instrument for bipolar coagulation to connection socket (17).
- Assign the blue foot switch pedal to the current type bipolar coagulation using the switch button (18).
- Activate the current type by pressing the foot switch. The activation LED (13) must light up, and the HF activation signal must sound.
- Press the selection button AUTOCOAG bipolar (16).
- Touch a saline-soaked swab with the clean electrodes of the bipolar instrument. The activation LED (13) must light up, and the HF activation signal must sound.



Risk of malfunction, risk of injury!

If the HF activation signal can be heard when no foot switch is connected, the unit is defective.

- The unit must not be operated.
- A technical inspection must be performed.



Risk of malfunction, risk of injury!

If the HF activation signal can be heard while the foot switch is connected but the blue pedal is not pressed, the foot switch is defective.

- The foot switch must not be operated.
- It must be replaced.



5 Safety Measures

5.1 General Information

Electrosurgical devices are high-frequency generators designed to create high voltages and currents. In order to avoid hazards to the patient, the operators and third parties, the procedure must be applied carefully, and the operating and safety notes must be strictly observed.

5.2 Patient Placement



Risk of burns from stray currents!

- The patient must be placed isolated from grounded metal parts. Particular care is required to ensure that the patient's limbs do not touch any metal structures either.
- Ensure the required high-frequency insulation against the operating table by a sufficient number of layers (insulating blankets). Since during the operation moisture, perspiration, etc. are to be expected, a waterproof foil must be used to prevent wetting of these layers which serve as high-frequency insulation.
- Fluid accumulation under the patient must be avoided under all circumstances. Use further dry cloth layers where appropriate.
- Keep areas with stronger perspiration, extremities touching the trunk or skin-on-skin contact apart from each other and dry (arm-trunk, leg-leg, breasts) using insulating blankets.
- The above requirements for insulation must also be fulfilled if the patient is repositioned during the operation.



5.3 Application of the Neutral Electrode

AWARNING

Danger of burns from missing warning signal in case of one-piece neutral electrodes!

A one-part neutral electrode cannot be monitored. In case of insufficient contact, there will be no warning signal! The control of the HF output power can be limited for some monopolar cutting currents by one-piece, large-area neutral electrodes that do not come into direct contact with the patient's skin (so-called capacitive neutral electrodes). This can lead to the full power being emitted to the patient.

- We therefore generally recommend using a two-piece disposable neutral electrode that ensures the continuous monitoring of the patient.
- If a monopolar HF current is to be used, then a neutral electrode must be placed on the patient. Also in cases where bipolar application is intended but additional application of monopolar electricity is likely in case of complications, a neutral electrode should at least be kept ready.

Electrodes and cables must be applied carefully. Pay particular attention to the following:

- The neutral electrode must butt against the patient's body as close to the operation field as possible, reliably, and with its whole surface.
- Safe contact of the neutral electrode must be ensured for the total duration of the high frequency application.
- Application of the neutral electrode to a limb must not impair perfusion.
- Electrodes and cables must be applied carefully. The supply cables to the high-frequency electrodes must be laid out without loops and so that they cannot touch the patient or other cables. Use only cables that are intended for use with the device according to the manufacturer's instructions.
- The current paths in the body should be as short as possible and run in the longitudinal or diagonal direction of the body, not across it, the latter particularly not on the chest. Any metal parts in or on the body should be removed if possible, insulated, or paid special attention to.
- After repositioning of the patient, the electrodes and cables must be controlled for proper attachment.
- Do not apply the neutral electrode above implants or other metal parts, nor above bone protrusions or scarred tissue. If necessary, prepare the application site by cleaning and degreasing it. Strong hair growth is to be removed. For removal do not use substances that desiccate the skin (e. g. alcohol).
- For removal of the neutral electrode, do not pull at the cable or the connector strap.



5.4 Working with the Active Electrode

Observe the following rules when using high-frequency surgery:

- When simultaneously using electrosurgery and monitors on the same patient, use only
 monitoring electrodes whose leads comprise protective resistors or HF reactors. Needle
 electrodes must not be used for monitoring. The active surgical electrode must not be used
 near ECG electrodes (minimum distance 15 cm).
- The high-frequency power should be set to the lowest possible setting for the respective application.

NOTICE

An insufficient effect while using the standard setting can, for example, be caused by poorly applied neutral electrodes, insufficient contact in plug connection, cables that are damaged under the insulation or encrusted electrodes. This must be checked, and defective parts must be replaced if necessary.

5.5 Work on Patients with Pacemakers or Other Implants

If patients have metallic implants, the HF current paths must, as a general rule, not lead through these implants. This must be considered when applying the active and neutral electrode, i. e. it is not allowed to apply the neutral electrode over endoprostheses.

In patients with active implants such as pacemakers or implanted electrodes, the application of the electrosurgical unit poses a risk. Application could result in irreparable damage to the active implant or impairment of its function. The following guidelines are to be observed.

- Such patients are to be monitored using suitable technology.
- A defibrillator as well as an external pacemaker should be kept ready for use.
- The output power selected for the electrosurgical unit should be as low as possible.
- Do not use the active electrode of the electrosurgical unit closer than 15 cm to the implant or its electrodes. Wherever possible, use bipolar technology.
- The rules of application such as application of the neutral electrode must be carefully followed.

In case of metal parts present in the body, such as hip stems, ensure that they are not exposed to HF currents under any circumstances. HF currents can lead to surface damage (fusing) resulting in a notching effect that can cause fatigue bending fractures of the metal part, even several years later.

5.6 Cable Routing on the Patient

During simultaneous use of the monopolar and bipolar application ensure that the respective lines have a minimum distance of 10 cm to each other. The lines must be laid out loosely and without loops.



5.7 Operation with KLS Martin Argon Beamer System

If the unit is used together with an Argon Beamer System, the safety notes in the Instructions for Use of the MABS must be observed. The correct connection cable can be ordered as item no. 80-181-51-04.

5.8 Putting Down HF Instruments

In the intervals between use, never place HF instruments on the patient.

5.9 Accessories

Consider the following notes when using accessory equipment:

- The accessory equipment must have sufficient dielectric strength. The accessory equipment must be suitable for at least the voltages indicated in section 9 "Technical Data", page 34.
- The accessory equipment must be in flawless working condition.

5.10 Use of Two Electrosurgery Devices on One Patient

Basically, due to the increased risk of accidental burns by high-frequency currents two electrosurgical units should be used on one patient only if medical requirements necessitate this.

The following rules are to be observed in the concomitant operation of two electrosurgical devices:

- Only electrosurgical units with type CF application parts by Gebrüder Martin may be used.
- Each unit to be operated in monopolar application requires a neutral electrode separately attached to the patient. It is mandatory to comply with the instructions for proper application of the neutral electrode.
- Exclusively bipolar application of a unit obviates the need for the application of the neutral electrode for this unit.
- The expected current paths from each active electrode to the corresponding neutral electrode must not overlap or intersect. To this end, always place each neutral electrode in the immediate vicinity of the operation field.
- The total power of the simultaneously applied currents must not exceed 400 W.

5.11 Inadvertent emission of HF energy

Certain foot switches have magnetically activated reed contacts and must therefore be kept out of the immediate vicinity of magnetic fields (especially of MRTs).

Otherwise this can lead to unintentional activation of the foot switches. For the same reason, magnets must be kept out of the immediate vicinity of the foot switches.

Further information can be found in the respective instructions for use of the foot switch.



5.12 Additional Safety Notes

- During operations on low-diameter body parts, bipolar technology can be used to avoid unintentional coagulation in other areas.
- The use of flammable anesthetics, nitrous oxide (N²O) and oxygen should be avoided. Sparking on the active electrode will always occur when using electrosurgical devices. Combustible substances that are used for cleaning or disinfection or as solvents must be evaporated before application of electrosurgery. There is the risk of accumulations of combustible liquid under the patient or in depressions such as the navel or in body cavities such as the vagina. Liquid that has accumulated in such areas must be removed before the electrosurgical unit may be used. Always be aware of the possible presence of combustible endogenous gases. Materials saturated with oxygen, such as cotton or gauze, may catch fire from the sparks that occur during the intended use of the electrosurgical device.
- Transurethral resection of the prostate (TUR):
 It is generally known that during transurethral resection of the prostate (TUR), in particular during continuous or long HF application, flammable gases can easily form. These gases may accumulate in the upper area of the bladder. The addition of atmospheric air further increases the risk of an explosion. The contact of HF current with such gas accumulations may result in the potential risk of an unintentional explosion. The gas-air pocket in the bladder should therefore always be kept to a minimum.
- The unit may be combined with other units only by the manufacturer or with the manufacturer's consent.
- The electrosurgical unit may interfere with other electromedical devices.



6 Cleaning and Disinfection

The operator / user is responsible for cleaning and disinfecting the device. National regulations regarding cleaning and hygiene are to be complied with.



Risk of serious injury from defective accessories!

A non-functional unit or defective accessories can endanger the patient or user, as well as affect the intended functions of the unit.

- The unit and its accessories must always be kept in flawless working condition. Intact insulation, cleanliness and potential damage must be visually ascertained.
- Units and accessories unfit for further use are to be scrapped.
- The unit is not sterilizable.
- For cleaning and disinfection, the unit must be disconnected from the grid.
- When handling cleaning agents and disinfectants, no liquid may enter the unit. Use of spray is not allowed!
- Never clean the unit with scouring agents, disinfectants, or solvents which may scratch the casing or damage the unit.
- For surface cleaning and disinfection, follow the protocols recommended by the hospital or apply another nationally recognized and approved procedure.
- If accessories are to be disinfected (surface disinfection or immersion), the information of the manufacturer of the disinfectants concerning material compatibility, dosage, and exposure time must be respected.
- Non-sterilizable accessories, such as foot switches, must be subjected to disinfecting cleaning.
- Before start-up of the unit, all residual disinfectant must have been reliably removed.
- If accessory parts are to be disinfected or sterilized, the instructions for use of the accessory parts must be observed.



7 Safety Checks

At least every 24 months the device must be subjected to the following controls to be performed by people who, by virtue of their training, knowledge and experience gained from practical activities are able to properly perform such safety checks and who are not obliged to adhere to instructions with regard to such control activities.

NOTICE

The unit may be inspected only by qualified service staff of Gebrüder Martin or by companies expressly authorized to do so.

- Perform a visual check of the unit and its accessories for mechanical damage that may impair the function of the unit.
- Verify readability of safety-relevant labels.
- Check fuse inserts of the device fuse for nominal current and fuse characteristics.
- Test function according to Instructions for Use.
- Check the linear increase or decrease of the output power upon corresponding position of the rotary knob.
- Check the actual maximum output power with the maximum values (±20 %) specified in section 9 "Technical Data", page 34 on all HF outputs for all available operating modes under consideration of the indicated nominal load resistors.
- Check acoustic and visual signals during power output.
- Perform electrical check according to test report for periodic safety checks.
- The maximum leakage current value may not be more than 1.5× the initially measured values and may not exceed the limit value either.
- The initially measured values are indicated in the attached test reports from the initial commissioning.
- We recommend entering the safety checks in a medical product log and documenting the test results.
- If the unit is not intact and / or not safe to operate, it must be repaired, or the operator must informed about the danger the unit constitutes.

NOTICE

A service manual that contains technical details and a description of the safety checks is available upon request. The service manual helps persons or companies authorized by us to maintain / repair the unit. The service manual also lists all components available as spare parts.



8 Accessories

The unit may be used only with accessories, wear parts and single-use articles whose safety of use has been demonstrated by a declaration of conformity. Use of uncertified accessories by other manufacturers may constitute a source of danger. In case of doubt, consult the manufacturer.

The accessories certified for electrosurgical units by KLS Martin can be taken from the accessories catalog of KLS Martin which can be requested or directly downloaded from www.klsmartin.com.



9 Technical Data

Description	Data electrosurgical unit ME MB2		
Mains Connection	100 – 127 V / 220 – 240 V; 50 – 60 Hz adjustable via blind plug inside the unit		
Power consumption	without HF output approx. 26 VA at max. output power approx. 800 VA		
Protection class	I		
Classification acc. to German Medical Devices Act	II b		
Leakage currents LF / HF	acc. to IEC 60601-2-2		
Туре	CF; defibrillation-proof		
Nominal frequency	450 kHz		
Modulation frequency	30 kHz		
Connection socket	Beamer system		





Description	Data electrosurgical unit ME MB2			
HF output variables				
Current type	Output	Crest f	actor	Voltage
Cutting 1 (pure)	max. 400 W at 300 Ω	1.6 at 3	00 Ω	max. 2,300 VSS
Cutting 2 (blend)	max. 300 W at 300 Ω	1.9 at 3	00 Ω	max. 2,500 VSS
Endo mode	max. 100 W at 200 Ω	1.6 at 2	200 Ω	max. 2,800 VSS
Contact coagulation	max. 250 W at 200 Ω	3.4 at 2	200 Ω	max. 3,200 VSS
Spray coagulation	max. 120 W at 300 Ω	5.6 at 3	Ω 00	max. 9,000 VSS
Bipolar coagulation	max. 100 W at 100 Ω	2.1 at 2	200 Ω	max. 600 VSS
Operating mode	Intermittent INT 10 s / 30 s corr. to 25 % duty factor			
Mains fuses	100 – 127 V: T 8 A 220 – 240 V: T 4 A G 5 x 20 mm			
Signal level	HF display: 55 dB(A) (adjustable 50 – 60 dB; contact the Martin Service Center) Alarm: 65 dB(A)			
Weight	8.6 kg			
Radio interference suppression	Limit values acc. to EN 550 Radio interference immunit			
Dimensions W x H x D	405 x 135 x 380 mm			
Environmental conditions for transport and storage	Ambient temperature $-25 ^{\circ}\text{C}$ to $+70 ^{\circ}\text{C}$ (-13 $^{\circ}\text{F}$ to $+158 ^{\circ}\text{F}$) Relative humidity $10 - 100 ^{\circ}\text{M}$ Atmospheric pressure $500 - 1,060 \text{hPa}$			
Environmental conditions for operation	Ambient temperature $+10 ^{\circ}\text{C}$ to $+40 ^{\circ}\text{C}$ (-13 $^{\circ}\text{F}$ to $+1 ^{\circ}\text{C}$ to $+40 ^{\circ}\text{C}$ to $+40 ^{\circ}\text{C}$ to $+40 ^{\circ}\text{C}$ (-13 $^{\circ}\text{F}$ to $+1 ^{\circ}\text{C}$ to $+40 ^{\circ}\text{C}$ to $+40 $,	
C € 0287	Compliant with 93/42/EEC			




10 Diagrams

10.1 Performance Diagrams







Cutting 2 (blend)



- 1 = Power controller in maximum position
- 2 = Power controller in middle position




Contact coagulation




- 1 = Power controller in maximum position
- 2 = Power controller in middle position

Spray coagulation



Bipolar coagulation





- 1 = Power controller in maximum position
- 2 = Power controller in middle position





10.2 Voltage Diagrams


HF output voltage depending on the setting
















Guidelines and Manufacturer's Declaration Regarding Electromagnetic Compatibility (EMC)

Guidance and manufacturer's declaration according to IEC 60601-1-2, clause 6.8.3.201 a) 3) Table 201: Electromagnetic emission

The device ME MB2 is intended for use in an electromagnetic environment as specified below. The user should assure that it is used in such an environment.

Emission test	Compliance level	Electromagnetic environment – Guidance
RF emissions according to EN 55011	Group 2	The device uses electromagnetic energy exclusively for the performance of its internal functions. This may cause interference in nearby electronic equipment.
RF emissions according to EN 55011	Class B	The limits of this class are kept only in standby mode (without RF current activation)!
Harmonic emissions according to IEC 61000-3-2	Class A	The device is suitable for use in all establishments, including domestic establishments and those directly connected to
Voltage fluctuations / flicker emissions according to IEC 61000-3-3	Complies	the public low-voltage power supply network that supplies buildings used for domestic purposes.



Guidance and manufacturer's declaration according to IEC 60601-1-2, clause 6.8.3.201 a) 3) Table 201: Electromagnetic emission

The device is intended for use in an electromagnetic environment as specified below. The user should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – Guidance	
Electrostatic discharge (ESD) according to IEC 61000-4-2	±6 kV contact discharge ±8 kV air discharge	±6 kV contact discharge ±8 kV air discharge	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.	
Fast electrical transients / bursts according to IEC 61000-4-4	±2 kV for mains power lines ±1 kV for input and output lines	±2 kV for mains power lines ±1 kV for input and output lines	Mains power quality should be that of a typical commercial or hospital environment.	
Surge according to IEC 61000-4-5	±1 kV differential- mode voltage ±2 kV common- mode voltage	±1 kV differential- mode voltage ±2 kV common-mode voltage	Mains power quality should be that of a typical commercial or hospital environment.	
Voltage dips, short power interruptions and supply voltage fluctuations according to IEC 61000-4-11	<5% U_T (>95% dip in U_T) for V_2 period 40% U_T (60% dip in U_T) for 5 periods 70% U_T (30% dip in U_T) for 25 periods <5% U_T (>95% dip in U_T) for 5 seconds	$0\% \ U_T$ (100% dip in U_T) for V_2 period V_3 period V_4 (60% dip in V_4) for 5 periods V_4 (30% dip in V_4) for 25 periods V_4 (100% dip in V_4) for 5 seconds	Mains power quality should be that of a typical commercial or hospital environment. If the user requires continued operation during mains power interruptions, it is recommended to power the device from an uninterruptible power supply.	
Power frequency (50 / 60 Hz) magnetic field according to IEC 61000-4-8	3 A/m	3 A/m	Power-frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.	
Note: U_T is the a. c. mains voltage prior to application of the test level.				



Guidance and manufacturer's declaration according to IEC 60601-1-2, clause 6.8.3.201 b) Table 204: Electromagnetic immunity

The device is intended for use in an electromagnetic environment as specified below. The user should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – Guidance
Conducted RF according to IEC 61000-4-6 Radiated RF according to IEC 61000-4-3	3 V _{eff} 150 kHz to 80 MHz 3 V/m 80 MHz to 2.5 GHz	3 V _{eff} 150 kHz to 80 MHz 3 V/m 80 MHz to 2.5 GHz	Portable and mobile RF communications equipment should be used no closer to any part of the device (including cables) than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d=1.17\sqrt{P}$ $d=1.17\sqrt{P}$ for 80 MHz to 800 MHz $d=2.33\sqrt{P}$ for 800 MHz to 2.5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site surveya, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol.
Note 1:	At 80 MHz and 800 MHz, the higher frequency range applies.		
Note 2:	These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.		

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio stations and AM and FM radio and TV stations, cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level specified above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.

^b Over the frequency range of 150 kHz to 80 MHz, field strengths should be less than 3 V/m.



Recommended separation distances between portable and mobile RF telecommunications equipment and the KLS Martin ME MB2 according to IEC 60601-1-2 clause 6.8.3.201 b), Table 206

The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The user can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF telecommunications equipment (transmitters) and the device as recommended below, depending on the maximum output power of the communications equipment.

	Separation distance according to frequency of transmitter (in m)				
Rated maximum output <i>P</i> of transmitter in watts	150 kHz to 80 MHz $d = 1.17\sqrt{P}$	80 MHz to 800 MHz $d = 1.17\sqrt{P}$	800 MHz to 2.5 GHz $d = 2.33\sqrt{P}$		
0.01	0.12	0.12	0.23		
0.1	0.37	0.37	0.74		
1	1.17	1.17	2.33		
10	3.7	3.7	7.4		
100	11.7	11.7	23.3		

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1:	At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.
Note 2:	These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.



12 **Ecological Information**

12.1 Disposal of Packing

Gebrüder Martin will, as a matter of course, take back the full packaging if so desired, to recycle as many parts of the packaging as possible.

If you do not wish to make use of this offer, you can dispose of the packaging with the normal paper or domestic garbage.

12.2 Ecological Aspects of Operation

If the treatment is interrupted and the unit is not being used for an extended period, we recommend you to switch it off for safety as well as economic reasons (energy saving).

If disposable products are used in a treatment, please note that such articles must first be carefully cleaned, disinfected and – where applicable – sterilized before they are disposed of as domestic or hazardous waste. Infected sharp parts of single-use products must be handled like any other "sharps" (cannulas, needles and scalpels) in accordance with valid regulations (disposal via germ-proof and puncture-proof containers).

When vaporizing tissue, ensure that you do not inhale the combustion fumes that form by design over a longer period in concentrated form. During regular use of the unit, no pollutants other than these combustion products will be formed.

For removal of such combustion products, a fume suction system can be used.

12.3 Disposal of the Unit

In designing the unit, we tried to avoid using composite materials wherever possible. This allows a high degree of recycling after the lifetime of the unit. We therefore offer to take the unit back for proper disposal and recycling.

Be sure to observe your national/local rules and regulations governing disposal!



Marking of electric and electronic equipment in accordance with Directive 2002/96/EC (WEEE Directive) and the German Electrical and Electronic Equipment Act (ElektroG)

This symbol on the product or its packaging indicates that the product may not be disposed of as normal household garbage.

KLS Martin Group

KLS Martin France SARL

68200 Mulhouse France Tel. +33 3 89 21 66 01 france@klsmartin.com KLS Martin UK Ltd.

Reading RG1 3EU United Kingdom Tel. +44 1189 000 570 info.uk@klsmartin.com KLS Martin do Brasil Ltda. CEP 04.531-011 São Paulo

Brazil Tel. +55 11 3554 2299 brazil@klsmartin.com Gebrüder Martin GmbH & Co. KG

Representative Office, 121471 Moscow Russia Tel. +7 499 792 76 19 russia@klsmartin.com

Martin Italia S.r.l.

20864 Agrate Brianza (MB) Italy Tel. +39 039 605 67 31 italia@klsmartin.com Nippon Martin K.K.

Osaka 541-0046 Japan Tel. +81 6 62 28 90 75 nippon@klsmartin.com **KLS Martin Australia Pty Limited**

Artarmon NSW 2064 Australia Tel. +61 2 9439 5316 australia@klsmartin.com Gebrüder Martin GmbH & Co. KG

Representative Office, 201203 Shanghai China Tel. +86 21 5820 6251 china@klsmartin.com

Martin Nederland/Marned B.V.

1270 AG Huizen The Netherlands Tel. +31 35 523 45 38 nederland@klsmartin.com KLS Martin L.P.

Jacksonville, FI 32246 USA Tel. +1 904 641 77 46 usa@klsmartin.com KLS Martin Asia Snd. Bhd.

14100 Simpang Ambat, Penang Malaysia Tel. +604 505 7838 malaysia@klsmartin.com Gebrüder Martin GmbH & Co. KG

Representative Office, Dubai United Arab Emirates Tel. +971 4 454 1655 middleeast@klsmartin.com



Gebrüder Martin GmbH & Co. KG A company of the KLS Martin Group

KLS Martin Platz 1 · 78532 Tuttlingen · Germany Postfach 60 · 78501 Tuttlingen · Germany Tel. +49 7461 706-0 · Fax 706-193 info@klsmartin.com · www.klsmartin.com