
Gebrauchsanweisung
Operating Instructions
Instrucciones de uso
Mode d'emploi
Istruzioni d'uso

Operationsleuchten	ML 1001/702/701/501/300/200
Operating Lights	ML 1001/702/701/501/300/200
Lámparas de quirófano	ML 1001/702/701/501/300/200
Lampes opératoires	ML 1001/702/701/501/300/200
Lampade scialitiche	ML 1001/702/701/501/300/200

© by Martin Medizin-Technik

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© by Martin Medizin-Technik

This document is under the copyright protection of **Martin Medizin-Technik**.

This document must not be modified, extended or reproduced or transmitted to a third party, without prior written permission by **Martin Medizin-Technik**.

Martin Medizin-Technik

Ludwigstaler Str. 132
D-78532 Tuttlingen

Telephone: + 49-74 61/70 60
Telefax: + 49-74 61/70 61 93

In case of technical questions or replacement parts are required, please contact our hotline service.

Martin Medizin-Technik

Hotline:
Telephone: + 49-74 61/70 63 43
Telefax: + 49-74 61/70 61 90

Replacement parts:
Telephone: + 49-74 61/70 61 99
Telefax: + 49-74 61/70 61 90

Hotline service available on workdays from 8:00 a.m. to 5:00 p.m.

We reserve the right to make changes in designs and devices.

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1 Introduction

1.1 Foreword

In their construction and technical overall conception in the field of lighting engineering, the **Martin** Operating Lights* outstandingly fulfill the high requirements which are set by today's operating technology for the lighting of the operating area.

The lights bear the EC-sign, i.e. they meet the basic requirements stipulated in the EC-directive on medical products.

In addition, the lights ML 1001 N, ML 701 N and ML 501 are UL approved.

These instructions for use will provide detailed information as to how to put the operating lights in service and use them for the intended purpose.

Being a user or a service technician authorized by the manufacturer, please read these instructions carefully prior to using the unit in the operating room, in order to make yourself acquainted with the handling and the operation of the operating light.

The manufacturer warrants that this operating lights have been properly and carefully tested before leaving the factory premises. In addition, we are prepared to assume continued responsibility for the safety, reliability and proper functioning of the operating lights as long as all checks and repairs required are carried out by technical staff that has been duly authorized by Martin to perform such work, and in as much as the operating lights are used only for the intended purposes.

Martin will assume no liability whatsoever for any damage that has been caused by making improper use of the product!

* The terms "operating light" or "light" respectively used in these operating instructions not only refer to the luminous element, but also to the mounting frame, the tracking and spring arms, the stands, the control boxes, etc.!

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1.2 Symbols within these Operating Instructions

Texts marked by this symbol and set in bold letters or italics contain very important information, like for example, the prevention of hazards to health!

Please observe these texts under any condition!

Texts marked by this symbol and set in bold letters or italics contain very important information, like for example, the prevention of material damage!

Please observe these texts under any condition!

Texts marked by this symbol and set in bold letters or italics contain comments/ hints and other useful information.

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1.3 General Information

1.3.1 Validity

The descriptions within these operating instructions refer to the **Martin** operating lights with the designation

- ML 1001;
- ML 702; *
- ML 701;
- ML 501;
- ML 300;
- ML 200 and

their respective versions and options.

Note the designation on the light, which is imprinted on the light head, and look for its specific description on the respective pages of these instructions.

The other chapters are valid for all lights.

*** Lights with the designation ML 702 are exclusively mounted in pairs! However, the correct designation ML 702/ML 702 is reduced to ML 702 in these instructions for use!**

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1.3.2 Proper Use for the Intended Purposes

The lights may only be used for purposes being in accordance with their design:
i.e., for lighting an operating or examination area by surgeons, doctors, or appropriately trained paramedics.

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1.3.3 Entry Inspection

Please check the light(s) and the accessories for shipping damage and completeness immediately after having received them.

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1.3.4 Complaints

Damage may only be claimed for if the vendor or forwarding company are immediately being informed.

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Immediately establish a certificate of damage and send it to the nearest agent of the manufacturer or directly to the manufacturer.

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If possible, use the original packages when returning goods.

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Enclose the following accompanying documents with the returned goods:

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- name and address of the owner
- type and device number
- description of the defect

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1.3.5 Warranty

For the operating lights, we grant a warranty of 12 months, starting from the day they were delivered to the user.

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Within this warranty period, our competent service agencies or the manufacturer remedy all defects which were demonstrably caused by manufacturing faults or faulty material for free.

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Parts subject to wear like lamps and replaceable hand grips are excluded from this guarantee.

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The general guarantee and warranty terms of **Martin Medizin-Technik** are valid.

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enthalten.

2 Safety Instructions

2.1 Safety Instructions for the Lights ML 1001/702/701/501/300/200

All Martin operating lights are high-quality products that have been designed and manufactured in accordance with technical standards. They fully meet, at the time when they leave the factory, all safety-related requirements.

In order to preserve this condition, you as the user are requested to refrain or avoid all actions that may have a negative influence on the safety of the lights, and to have repaired damaged parts immediately or notify the manufacturer of the damage!

Damaged lights may be a life- or health-threatening hazard for the user and/or the patient!

For this reason never use damaged lights!

Only handle, switch on or off the lights in accordance with the information provided in these instructions for use!

The mounting of the permanently installed lights and the repair of all lights may be carried out either by the manufacturer themselves or by a person or company explicitly authorized by the manufacturer!

In case the repair is executed by a person or company authorized by the manufacturer, the operator of the light is asked to demand a certificate containing the type and the extent of the repair by the person or company carrying out the repair.

This must contain

- date of the service,
- type and extent of the service,
- name of the company, and
- the signature of the service technician.

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
If repair work has been carried out by a third party, the equipment or parts repaired must be labeled with the repairer's identification mark.

After having carried out the mounting and service tasks, the service technician has to check all functions of the lights, spring and tracking arms, and the power supply units for their perfect functioning and operational safety!

The mobile stand lamps, however, may also be assembled and erected by non-experts, proceeding in accordance with the assembly instructions given in section 5.10 - 5.10.6.

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Laser pilot (Option)

The built-in laser pilot is a Class 2 laser that meets all requirements set by EN 60825/IEC 825/VDE 837. These requirements as well as pertinent regulations for the prevention of accidents (VBG 93, as far as Germany is concerned) must be observed in using the laser pilot.

The illustrated warning labels are located on the control panel where the laser pilot is activated.

The laser beam exit opening is located on the central handle. Whenever the operating light is repositioned with the help of this handle (or, in the case of the ML 1001, by using the hand rail), the laser beam will be switched on automatically. After repositioning has been completed, the laser pilot light will persist for another 2 s.

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Notes concerning prevention of accidents:

- **Due to the eyelid closure reflex, the laser radiation poses no danger to non-narcotized persons.**
- **In the case of narcotized patients, however, take care that the laser radiation will not hit the opened eye for more than 0.25 s. Cover the patient's eyes by using suitable means.**

The user has to assume full responsibility for the observation of the regulations mentioned above.

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2.2 Supplementary Safety Instructions for the Light ML 702

In contrast to the other lights, the ML 702 light is not operated with low voltage, but with 230 V AC! This voltage is led into the light head.

Even if the light is switched off, 230 V AC are available on the lines between the control box and the light head due to the accumulator-driven emergency power supply.

When triggering the lamp within the light, 4,500 V AC are generated at the inside of the light head and fed via the focus of the lamp!

This high-voltage may even be generated within the light head if the focus has been removed and the light been switched off!

In order to avoid a life-threatening electrical shock

- *never dismantle, repair or modify any parts of the ML 702 light!*
- *never reach into the light when the focus has been removed!*

In case the light or any of its parts have been damaged,

- *secure the lights or the parts against persons who might touch the light, if necessary, by putting up appropriate danger signs!*
- *cut off the power supply of both lights (note that there may nevertheless be some voltage in the lines leading from the control box to the lamp in the light head)!*
- *immediately call the service technician!*

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3 Versions and Combinations of the Lights

3.1 Versions of the Lights

The lights are being delivered in different versions. These versions differ

- in terms of the **diameter of the lamp housing**, e.g.:
 - ML 1001 corresponds to lamp housing diameter of 1000 mm
 - ML 701 corresponds to a housing diameter of 700 mm, etc.
- in terms of the **handling**, e.g.:
 - The basic version ML 701 contains an electronical control system with control panel at the hand grip of the light.
 - The version ML 701 N is a manual version with a control panel at the light head.
- in terms of the **type of mounting**, e.g.:
 - ceiling mounting
 - wall mounting
 - stand mounting
- in terms of the **options** by means of which the lights may be equipped, like for example:
 - laser pilot
 - MedTV-Pro-System
- in terms of the **lamp technology**, e.g.:
 - halogen lamp
 - metal halide lamp

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- The operating lights with the designation

- ML 1001
- ML 701
- ML 501
- ML 300
- ML 200

are equipped with one or two halogen lamp(s) respectively.

- The operating lights with the designation

- ML 702

are equipped with a metal halide lamp each.

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The lights with the designation

- ML 1001
- ML 701
- ML 501
- ML 300
- ML 200

may be mounted as individual light but also as two- or three-unit aggregate (with identical or different equipment); the combination has no impact on the handling of the lights. The handling of the individual light is light-specific.

The lights with the designation ML 702 are exclusively mounted in two-unit-aggregates with two identical lamps.

The handling of these lamps has to be
- as a matter of principle and due to the requirements of the system
- carried out in pairs!

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4 Description of the Lights

By using the new **Martin** FK 157 filter technology, these lights produce the **Martin** Super Cold Light.

4.1 Description of the Light ML 1001

The **Martin** ML 1001 operating light is available in two versions:

- The basic model ML 1001, where all functions can be electronically controlled via a central control panel in the tube connection piece.
- The ML 1001 N is the manual version with a control panel in the tube connection piece.

The **Martin** ML 1001 operating light provides the following accessory equipment:

Martin Anti-Drift System (ADS)

In the case of the **Martin** ML 1001, the operating light is arrested and kept in position by means of the built-in magnetic brakes. As soon as the user touches either the hand rail or the handle of the operating light, the magnetic brakes located in the rotary joints are released and the operating light can be repositioned with ease and precision.

Upon releasing the rail or the central handle, the operating light will be locked securely in its new position by the **Martin** Anti-Drift System.

Martin Laser pilot

The laser pilot built into the operating light exactly marks the center of the light field. This makes it particularly easy for the user to focus the light field with utmost precision even if the operating area is very small. Whenever the operating light is moved/repositioned, the laser pilot will be activated. When the operating light has come to a standstill in the course of repositioning, the laser pilot continues to be in operation for another two seconds. The laser pilot can be deactivated via the control panel.

Martin MedTV-Pro

The **Martin** MedTV-Pro is a medical television system to be used for documentation, consultation or further training/education purposes. The high-performance camera may be mounted either on the operating light itself or on a separate camera support.

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4.2 Description of the Light ML 702

There is only one version of the **Martin** ML 702 operating light available: a two-unit aggregate, with another, identical ML 702.

For reasons of construction, the lights are not equipped with a replacement lamp, and may thus only be operated in combination with a second operating light.

By applying a new kind of lamp technology, i.e. a metal halide lamp, a high light efficiency is achieved while there is only a minor heat build-up in the wound area. Compared to the ordinary operating lights with halogen bulb technology, the costs for operation and service are thus relatively low.

The handling of the light is carried out by means of a wall-mounted control box.

This control box contains a switching unit for switching on and off the two lights. Two operating hours counters for detecting the service life of the lamp belong to this switching unit, a logic control for a 15-minute-afterglowing of the most recently switched-off lamp, and an energy storage unit for bridging short-time power failures.

In contrast to other lights, this light is equipped with a metal halide lamp.

The characteristic technical features of this lamp are:

- A temperature similar to that of daylight, amounting to 4,500 Kelvin.
- A light efficiency of about 45 % of the absorbed energy.
- An intensity of light of 150,000 lx approx. for a power input of 70 W.
- A very low emission of infrared rays.
- With a service life of more than 5,000 hours, a five to ten times longer service life than that of halogen bulbs.

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4.3 Description of the Light ML 701

The **Martin** ML 701 operating light is available in two versions:

- The basic model ML 701, where all functions can be electronically controlled via a central control panel located on the handle of the operating light.
- The ML 701 N, a manual version with a control panel located on the body of the operating light.

Optionally, the ML 701 can be combined with the following accessory equipment:

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4.4 Description of the Light ML 501

The **Martin** ML 501 operating light is available in two versions:

- The basic model ML 501, which is a manual version with the control panel located on the operating light.
- The ML 501 E, where all functions can be electronically controlled via a central control panel located on the handle of the operating light.

Optionally, the ML 501 can be combined with the following accessory equipment:

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4.5 Description of the Lights ML 300/ML 200

The examination lamps ML 300 and ML 200 are characterized by a relatively high luminous intensity. Thus, they are particularly useful for operations and examinations where a concentrated and uniform light field is required.

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4.6 Options for the Lights ML 1001/701/501

These operating lights have optionally, at the ceiling tube, a card for the commutation mains/battery.

A double LED on the soffit indicates optically the power source actually in use.

Green LED: Power supply (mains)

Red LED: Power supply (battery)

Note:

If neither of the LED's is lit, the battery does not supply any power.

The LEDs are connected to the emergency relais with a three-pin-plug. For combined lights, the LED does always concern the main light.

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5 Control and Display Elements & Operating the Operating Lights

5.1 Operating the Light ML 1001

All functions of the **Martin** ML 1001 operating light can be easily controlled via a separate control panel that has been arranged in a straightforward manner. The various functions available have been marked with a symbol. To activate a function via the respective pushbutton, no more than a touch of the finger is required.

The symbols on the panel represent the following functions:

- 1 Operating light OFF/ON
(or use wall switch, as the case may be)
- 2 Laser pilot (option)
Laser stand-by button
(the green LED will come on)
- 3 Malfunction indicator lamps:
 - Main lamp I defective
 - Main lamp II defective
- 4 Decreasing light intensity
- 5 Light intensity indicator
- 6 Increasing light intensity
- 7 Reducing light field diameter
- 8 Light field indicator
- 9 Increasing light field diameter

Please note:


If one of the malfunction indicator lamps (3) has come alight, the respective halogen lamp should be replaced as soon as possible. While the halogen lamps have a rated life of between 600 and 1000 burning hours, this average service life may be significantly reduced in consequence of mechanical vibrations.

The modern **Martin** operating lights are equipped with microprocessor-controlled electronics. In case an error occurs, one of the LEDs in the luminous intensity indicator (5) will be flashing. At the same time, either an LED display moves back and forth within the light field indicator (8) or all LEDs are flashing at the same time. In case of under- or overvoltage, the two malfunction indications (3) will be flashing, thus signalling that a change of lamp is due.

In all the cases mentioned above please contact the technical service since the repair work required may only be carried out by service technicians.

Note:

In general, this does not lead to a failure of the lights. It is not necessary to switch off the lights and the operation may be terminated. In certain cases, the nominal luminous intensity will be reduced to 70 % approx.

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
Martin Anti-Drift System (ADS)

The **Martin** Anti-Drift System (ADS) allows easy and high-precision positioning of the operating light. As soon as the operating light has come to a standstill, it will be fixed in position by the built-in magnetic brakes. These brakes will be automatically released whenever the hand rail surrounding the operating light (or the central handle) is touched. Thus, repositioning is possible with utmost ease.

Accordingly, the operating light will again be arrested in its new position as soon as the user withdraws his hand from the rail or handle.

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Martin Laser pilot

To put the laser pilot into operation, press the pushbutton 2 on the control panel. As soon as the central handle or the hand rail is touched, the laser will be activated.

Upon releasing the rail or handle, the laser light will persist for another 2 seconds.

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5.3 Operating and Display Elements of the Light ML 702

The operating and display elements of the **Martin** ML 702 operating lights are placed in a wall-mounted control box.

- 1/2 On/Off switch for the lamps of the lights 1 and 2
- 3/4 Green control lamps for the lights 1 and 2
- 5/6 Red warning lights for the lights 1 and 2
- 7 Electronical beeper (integrated in the control box)
- 8 Brief instructions for use

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5.3.1 Description of Functions of the Light ML 702

In contrast to the other lights described in these operating instructions, the lights of the ML 702 version contain a metal halide lamp each (with the outstanding advantages described above).

The metal halide lamp needs a high-voltage of 4,500 V for triggering. The operating voltage after triggering ranges from 70 to 90 V approximately.

The high-voltage and the operating voltage are generated in a device integrated in the light head. The power supply of the device or the lamp respectively is provided for via the control box.

Besides the operating and display elements, the control box contains an UPS (uninterruptable power system), two operating hours counters for the two lamps and a logic control.

Switching on the lights

The lights are switched on or off by means of the Lamp1 and Lamp2 switches at the control box.

After having switched on the lights the corresponding **green** pilot lamp on the control box is lighted.

Additionally, the red warning lamp flashes three times (operating hours counter test).

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The lamps in the lights may possibly be lit up several seconds after having switched on the lights and may reach the maximum light and color temperature of 4,500 Kelvin after four minutes approximately!

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To make sure that the operating area is still running to capacity, both lights have to be switched on and the operating area has to be adjusted!

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Switching off the lights

For reasons of safety, the most recently switched-off lamp will be afterglowing for 15 minutes.

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The most recently switched-off lamp needs this time to cool down unless the lamp may be triggered again.

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Display elements

The display elements of the lights (at or in the control box respectively) are optical (control/ warning lamps) as well as acoustic devices (electronical beeper).

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Optical signals

- **Green control lamps are flashing**
If the green control lamps are flashing after having switched on the light, this means that the end of the service life of the accumulator in the UPS mode has been reached. In this case, the UPS has to be replaced as soon as possible by an authorized service technician. The service life of the accumulator depends on the constant ambient temperature. A low temperature increases and a high temperature decreases the service life. Previous experience has shown that the UPS has to be replaced every three to six years!

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- **Red control lamps are flashing three times**
When switching on the lights, the red warning lamps flash three times signaling the proper operation of the operating hours counters of the lamps.

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- Red warning lamps are flashing permanently

After 4,500 operating hours, the red warning lamp of the respective lamp starts to flash.

The flashing frequency increases automatically, until, after 6,000 operating hours, the light is lit up permanently. Within this interval, the respective lamp has to be replaced.

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Emergency power supply

In order to bridge a mains failure within the building, there is the UPS with an accumulator in the control box, whose maximum charge is sufficient to supply both lights with energy for 5 minutes approximately. Within this interval, the emergency supply of the building (stand-by power supply) has to be used!

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Acoustic warning signals

- One beep tone every ten seconds

To signal the operation of the accumulator in case of a mains failure, the UPS in the control box emits **a** beep tone every ten seconds.

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- Two beep tones every ten seconds

If the accumulator has been discharged for the most part and thus there is the danger of the interruption of the power supply of the lights, the UPS in control box, emits **two** beep tones every ten seconds approximately.

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To secure the operability of the emergency power supply of the lights by means of the accumulator, the accumulator has to be charged at least every 24 hours!

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5.3.2 Putting the Lamps in Service or Switching off of the Light ML 702

We presume that

- *the lights ML 702 and the control box were mounted by a service technician authorized by the manufacturer (according to the specifications in the Martin service manual ML 702) and the initial putting in service and the function test have been carried out!*
- *the accumulator of the UPS has been charged permanently for at least 24 hours, thus attaining its maximum charge capacity!*
- *emergency power supply of the building will be triggered after one minute of mains power supply failure within the building at the latest.*
- *you have acquainted yourself with the control elements of the lights and especially with the optical and acoustic signals of the display elements and know their meaning (see chapter 5.3, 5.3.1, and 8.2 to this end).*

- Switch on **both** lights by means of the On/Off switches (switch in position 1 corresponds to **ON**).
- Mind the control lamps and if necessary the warning lamps and/or the acoustic signals!
- Align both lights to the operating area.

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- Focus the light fields of the lights to the desired diameter by means of the sterilizable handles (1)

The light field diameter and thus the light intensity may be adjusted by turning the sterilizable handle to the right or to the left respectively.

The adjusted intensity (in kilolx) may be read out at the scale factor (2), which is located above the notch joint bow of the sterilizable handle.

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- Switch off the lights by means of the On/Off switches at the control box (Position 0 corresponds to **OFF**).

The most recently switched-off lamp will persist for another 15 minutes! Within this interval, the switched-off lamp is not ready for use at first! It is only after this interval that the lamp has sufficiently cooled down as to be triggered again!

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**5.4 Operating the Lights
ML 701/501 E**

All functions of the **Martin** ML 701/501 E operating light can be easily controlled via a separate control panel that has been arranged in a straightforward manner. The various functions available have been marked with a symbol. To activate a function via the respective pushbutton, no more than a touch of the finger is required.

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- 1** Operating light OFF/ON
(or use wall switch, as the case may be)
- 2** Laser pilot (option)
Laser stand-by button
(the green LED will come on)
- 3** Malfunction indicator lamps:
- Main lamp I defective
- Main lamp II defective
- 4** Decreasing light intensity
- 5** Light intensity indicator
- 6** Increasing light intensity
- 7** Reducing light field diameter
- 8** Light field indicator
- 9** Increasing light field diameter

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5.5 Operating the Lights ML 701 N/501

All functions of the **Martin** ML 701 N/501 Operating Lights can be easily controlled via a separate control panel that has been arranged in a straightforward manner. The various functions available have been marked with a symbol. To activate a function via the respective pushbutton, no more than a touch of the finger is required.

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- 1 Operating light OFF**
(or use wall switch, as the case may be)
- 2 Operating light ON**
(first switch on the wall switch if provided)
- 3 Indicator lamp for monitoring the two light sources.**
If one of the main light sources has become defective, this indicator lamp will light up (red).

To increase or decrease the light field diameter, rotate the sterilizable handle clockwise or anticlockwise as appropriate.

Please note:

If the malfunction indicator lamp (3) has come alight, the respective halogen lamp should be replaced as soon as possible.

While the halogen lamps have a rated life of between 600 and 1000 burning hours, this average service life may be significantly reduced in consequence of mechanical vibrations.

5.6 Operating the Lights ML 300/ML 200

In the case of the ML 300 and ML 200 examination lamps, the ON/OFF switch is located on the handle.

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
5.7 Operating the Wall-Mounted Operating Lights ML 501/501 E/300/200

All wall-mounted operating lights are provided with a wall-mounting control box containing a power switch. This switch must be set to ON prior to operating the operating light via the control panel located on the body of the operating light.

The wall-mounting control box also contains the line fuses.

	220 - 240 V	100 - 125 V
ML 501	2,0 AT slow-blow.	4,0 AT slow-blow.
ML 300	1,25 AT slow-blow.	2,5 AT slow-blow.
ML 200	1,25 AT slow-blow.	2,5 AT slow-blow.

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5.8 Replacing the Sterilizable Handle

To remove the sterilizable handle, just press the locking bar located on the upper end of the handle to release the lock. The handle can then be withdrawn with ease.

When inserting the handle, be sure that it clicks in place audibly. To achieve this, it might be necessary to rotate the handle just a little clockwise or anticlockwise.

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5.9 Replacing the Halogen Lamp

5.9.1 Replacing the Metal Halide Lamp of the Light ML 702

In contrast to the halogen lamps of the remaining operating lights, the replacing of the metal halide lamp of the light and the consequently required reset of the operating hours counter may only be carried out by an authorized service technician!

The instructions for carrying out these service tasks are to be found in the **Martin service manual ML 702**.

Please note the chapters

- 2 Safety Instructions
- 10 Maintenance and Repair of the Light Combination
- 10.1 Removing the Focus
- 10.2 Replacing the lights
- 10.4 Installing the Focus
- 11 Maintenance and Repair of the Control Box
- 11.3 Resetting the Operating Hours Counter

In order to avoid a life-threatening electrical shock, please observe the safety instructions mentioned in these descriptions under any circumstances!

5.9.2 Replacing the Halogen Lamps of the Lights ML 1001/701/501

The procedure for replacing a defective halogen lamp is identical for the ML 1001, ML 701 and ML 501 operating lights.

WARNING!

Never perform work on the operating light while it is switched on.

First remove the three countersunk-head screws located on the black disk on the underside of the operating light; for this, use (for example) a coin and rotate the screws counterclockwise.

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
The entire focusing unit may now be withdrawn from the operating light.

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The electronic lamp monitoring system ensures that the defective lamp is automatically moved aside with the help of a servo-motor provided for this purpose. At the same time, the second light source is moved into the focus of the reflector.

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Accordingly, always replace the defective lamp located on the edge of the lamp holder.

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Be sure to use only original **Martin** halogen lamps with the following ratings:

Leuchten- typ	Halogen- lampe	Kat.-Nr. (6 Stck.)
ML 1001	24 V/ 250 Watt	86-910-02
ML 701	24 V/ 150 Watt	86-907-01
ML 501	24 V/ 120 Watt	86-905-01

Type of operating light	Halogen bulb	Cat/Order No. (6 pcs.)
ML 1001	24 V/250 W	86-910-02
ML 701	24 V/150 W	86-907-01
ML 501	24 V/120 W	86-905-01

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WARNING!

When replacing a lamp, always check the cables for scorches in the lamp holder area. If such defects are found, please call your serviceman.

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Never touch the halogen lamps with your bare hands! Use the protective covering provided to insert the lamp.

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
When mounting the focusing unit, be sure to insert the plug (1) into the rectangular cutout (2) of the assembly support.

To do this, manually insert the focusing unit in such a manner that it is flush with the assembly support. Then lock the unit in place by rotating the three countersunk-head screws clockwise.

Please switch on the operating light to check it for proper functioning. The indicator lamp (see chapter 5.4, position 3 on the control panel) must now be off.

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5.9.3 Replacing the Halogen Lamps of the Lights ML 300/200

Prior to opening the operating light, switch it off and let it cool down.

Remove the cover located on the upper side of the operating light.

Replace the defective halogen lamp (rating: 24 V/50 W). Be careful not to touch the lamp with your bare fingers.

Reinstall the cover on the operating light.

To check it for proper functioning, switch on the operating light.

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
5.10 Assembling the Stand Lamps & Putting the Lamps Into Service

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5.10.1 Assembly/Erection Instructions Relating to ML 701 SO/ML 701 N SO/ML 501 SO

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
Insert the stand tube into the foot/base of the stand, then fix it in place by tightening the M 12 hexagon socket screw provided (use 6-mm hexagon-socket wrench).

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Connect the plug contacts of the stand tube with the plug contacts of the adjustable swivel arm.

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
WARNING!

First connect the (green-and-yellow) PE conductors with each other, then the non-interchangeable power supply lines. Insert the swivel arm into the stand tube.

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
Lock the swivel arm in place by tightening the four M 5 x 8 recessed-head screws provided.

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Disengage the spring element of the locking sleeve by using a screw driver, then slide the sleeve backward and remove the nose element.

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Insert the pin of the lamp bracket into the swivel arm, then insert the nose element and slide the sleeve forward until it fully engages.

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

5.10.3 ML 501 SG/ML 300 SG Assembly Instructions

Insert the stand tube into the foot of the stand. Use the retaining bar, the snap ring and the M 8 x 35 hexagon head cap screw (wrench size 13) provided to lock/secure the tube on the underside of the assembly.

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Insert the completely assembled operating light arm into the stand tube.

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Netzstecker einstecken.

5.10.4 Putting the ML 501 SG/ ML 300 SG Into Service

Connect the power plug to a suitable mains socket.

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
ML 501:

This operating light is switched on/off via the power switch on the control panel located on the operating light.

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	220 - 240 V	100 - 125 V
ML 300 SG	1,25 AT	2,5 AT
ML 501 SG	2,0 AT	4,0 AT

Line fuses:

	220 - 240 V	100 - 125 V
ML 300 SG	1.25 A, slow-bow.	2.5 A, slow-bow.
ML 501 SG	2.0 A, slow-bow.	4.0 A, slow-bow.

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

Slip on the profile leg.

Insert the screw and spring washer into the bore, then tighten screw.

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
Use the key provided to tighten the screws until the aluminum profile underneath gets plate-shaped.

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Put the plastic sleeve in place on the stand foot assembly, then insert the stand tube into the foot. Lock/secure the tube on the underside of the assembly.

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To fix the cable on the stand tube, tighten the cable a little and attach it to the cable clips.

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5.10.6 Putting the ML 200 SG Into Service

Insert the power plug into a suitable mains socket.

Switch on the operating light via the mains switch.

Line fuses:

	220 - 240 V	100 - 125 V
ML 200 SG	1.25 A, slow-bow.	2.5 A, slow-bow.

The operating light can be switched on/off by using the switch provided on the handle.

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6 Accessories

	ML 1001
Order No.	Description
86-910-02	Halogen bulb, 24 V/250 W (6 pcs.)
86-900-00	Sterilizable handle for operating lights not equip- ped with a laser pilot
86-990-04	Sterilizable handle for operating lights equip ped with a laser pilot

	ML 702
Art. Nr.	Bezeichnung
86-970-00-04	Metalldampf- Halogenlampe (1 Stck.)
86-900-00	Sterilisierbarer Handgriff für Leuchten

	ML 702
Order No.	Description
86-970-00-04	Metal halide bulb (1 pcs.)
86-900-00	Sterilizable handle for operating lights

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	ML 701
Order No.	Description
86-907-01	Halogen bulb, 24 V/150 W (6 pcs.)
86-900-00	Sterilizable handle for operating lights not equipped with a laser pilot
86-990-04	Sterilizable handle for operating lights equipped with a laser pilot

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7 Maintenance & Servicing of the Light

The maintenance instructions for the lights ML 702 are defined in the Martin service manual and may only be executed by an authorized service technician!

To keep the operating light in good working order, all service and maintenance work must be carried out with due care.

7.1 Cleaning

All operating lights described above may be cleaned by using common cleaning and disinfecting agents.

In order to prevent humidity from penetrating into the light and the tracking and spring arms, do **not** spray these parts but only wipe them with cloths moistened with detergents or disinfectants.

7.1.1 Cleaning of the Sterilizable Handles

Sterilization:

The replaceable handle may be autoclaved at 134 °C at a maximum, using a common sterilizer. During sterilization, the handle may not be put under a mechanical load, since otherwise plastic deformation might occur.

Cleaning:

We recommend using the following procedure:

- Main cleaning cycle: alkaline cleaning
- Second cleaning cycle: adding a neutralizer
- Intermediate rinsing, using normal water
- Final rinsing, using deionized water

Service life:

The replaceable handles are subject to natural wear and tear if sterilized frequently. If fatigue occurs (e.g., cracks forming; discoloration), the handle must be replaced.

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7.2 Maintenance

At regular intervals - but at least once a year - all functions of the operating light, the stand and the power supply unit need to be checked for proper functioning and operating reliability by an authorized person.

Whenever a halogen lamp is replaced, the lamp holder must to be checked for potential defects. Remember that lamp holders are wearing parts!

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8 Troubleshooting

8.1 Troubleshooting ML 1001/701/501/300/200

Symptom	Potential Cause	Remedy
Operating light cannot be switched on.	Wall-mounted switch not set to ON.	Set wall-mounted switch to ON.
Malfunction indicator lamps I and II alight on the control panel.	Both halogen lamps defective.	Replace both halogen lamps.
Malfunction indicator lamp I or II alight on the control panel.	One halogen lamp defective.	Replace the halogen lamp.
Laser pilot beam not visible.	Laser pilot stand-by button (on control panel) not activated. Sterilizable handle without glass insert on the handle.	Press laser pilot stand-by button Use sterilizable handle with glass insert.
Brightness control not working.	Focusing unit not properly inserted.	Insert focusing unit properly.
Light field diameter cannot be adjusted.	Focusing unit not properly inserted.	Insert focusing unit properly.

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8.2 Optical/Acoustic Signals/Malfunctions of the Light ML 702

Signal	Meaning	Measure
- Switching on/off in position 1 - Green control lamp is not lit up - Lamp is not lit up	- Mains fuse in the control box is defective - No mains or emergency power supply	- Call authorized service technician
- Light of the operating light is too dark or too red or blue	- The lamp in the light has not yet attained the operating temperature	- Wait for 4 minutes approx., until the operating temperature has been reached
	- End of the lamp service life attained	- Mind red warning lamp - Call authorized service technician
- Switching on/off in position 1 - Green control lamp is lit up	Light is switched on	_____
- Switching on/off in position 1 - Green control lamp is not lit up	Light is switched on - Control lamp defective	- Call authorized service technician!
- Switching on/off in position 1 - Green control lamp is flashing	Light is switched on - UPS accumulator error	- Call authorized service technician Warning! Possibly no energy supply in case of mains power failure!
- When switching on the light, the red warning lamp flashes three times	- Correct functioning of the operating hours counter	_____
- During the operation of the light, the red lamp slowly starts to flash	- The lamp within the light has reached 4,500 operating hours	- The lamp still has another 500 operating hours and is to be replaced soon. - Call authorized service technician
- During the operation of the light, the red warning lamp flashes faster and faster or is lit up permanently	- The lamp has reached 4,500 or 5,000 operating hours respectively	- The lamp has to be replaced as soon as possible! - Call authorized service technician!
- Every 10 seconds approximately there is a beep tone coming from the control box	- Mains power failure! The energy supply is provided for by the accumulator in the control box!	- The emergency power supply of the building has to be triggered after four minutes approx. at the latest!
- Every ten seconds approximately, there are two beep tones coming from the control box	- Mains power failure! The accumulator in the control box has nearly been discharged! Within a few seconds, the energy supply of the lights will be	- The emergency power supply of the building has to be triggered immediately! - The accumulator has to be charged at least another 24 hours!

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* depending on mechanical vibration intensity, voltage fluctuations, and duty cycles.

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9.2 Specifications for the Light ML 702

Electrical power requirements Operation with adapting equipment (UPS)	110 or 230 V AC/50 Hz
Power input of the light	80 Watt
Max. illuminance measured at a distance of 1 m and for a small light field and a color temperature of 4,500 K	150 klux
Color rendition index Ra. (1 - 8)	94
Light efficiency	278 lm/W
Related to 100 klux: Total irradiance Thermal radiation > 700 nm Infrared radiation > 780 nm	360 W/m ² 14 W/m ² 3 W/m ²
Diameter Lamp housing Small light field Large light field	700 mm 190 mm 300 mm
Constant light field diameter without refocusing ranging from 80 cm to 130 cm	500 mm
Operating area illuminance Modeling (small/large light field > 20 %) Depth modeling (small/large light field > 50 %) Depth illumination (small/large light field > 20 %) (values according to DIN 5035, Part 3)	25 / 28 % 99 / 99 % 25 / 28 %
Cold start time 50 % of the brightness 100 % of the brightness	nach 2 min. nach 4 min.
Restart time	15 min.
Nominal current of the lamp	70 Watt
Actual service life of the lamp	5000 h

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Edizione Italiana copyright © 1998 - Tutti i diritti riservati.
Copyright by Gebrüder Martin GmbH & Co. KG
Gebrüder Martin GmbH & Co. KG