Deutsch
English
Español
Français
Italiano

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Table of Contents

1	General Information	22
1.1	Manufacturer	22
1.2	Hotline	22
1.3	Event Reporting Duty	22
1.4	Notes on this Document	23
1.5	Case Planning	23
1.6	Terms & Acronyms	24
1.7	Validity of this Document	24
2	Scope of Delivery	25
2.1	Inspecting the Delivery for Completeness and Correctness	25
3	Intended Use	26
3.1	Intended Purpose	26
3.2	Indications	26
3.3	Contraindications	26
3.4	Possible Side Effects	27
3.5	Target Patient Groups	27
3.6	Users	28
3.7	Ambient Conditions During Application	28
3.8	Limitations of Use	28
4	Use / Application	29
4.1	Description of the Components	29
4.1.1	Structure, Operating Principle and Performance Specifications	29
4.1.2	Materials and Implant Materials Used	29
4.1.3	Processing of the Implants	29
4.1.4	Compatible Products and Accessories	29
4.2	Before Use	30
4.3	Intraoperative Approach	31
4.4	After Conclusion of the Surgery	32
5	Cleaning, Disinfection and Sterilization	33
5.1	Storage and Transport	33
6	Ecological Information / Disposal	34
6.1	Packaging	34
6.2	Disposal	34
6.3	National Regulations	34

1 General Information

1.1 Manufacturer

We thank you for buying one of our products.

The product described in these Instructions for Use was specifically manufactured in accordance with a written prescription for an identified patient known by name and is intended solely to suit the patient's individual condition and needs. This product is therefore a custom-made product within the meaning of the applicable European regulations for the development, manufacture, placing on the market and monitoring of medical devices and as such does not bear any CE marking. As a custom-made product, it nevertheless fulfills the basic requirements for the safety and performance of medical devices within the meaning of the applicable European regulations.

We are the manufacturer of this product:

KLS Martin SE & Co. KG

A company of the KLS Martin Group

KLS Martin Platz 1 · D-78532 Tuttlingen · Germany

Phone +49 7461 706-0 · Fax +49 7461 706-193

info@klsmartin.com www.klsmartin.com

1.2 Hotline

Should you have any questions on how to handle the device or product or use it for clinical applications, please do not hesitate to contact the Product Management:

Phone: +49 7461 706-0

E-mail: info@klsmartin.com

NOTICE

Each specially manufactured implant and its packaging are marked with an individual case number. Please always provide the case number in case of a complaint.

1.3 Event Reporting Duty

All serious adverse events occurring in connection with the product must be reported to KLS Martin and the respective authorities without delay.

Possible danger to the life of others by shipping contaminated products!

In case of return shipments, send only clean and disinfected products in sterile packaging.

1.4 Notes on this Document

WARNING

Possible danger to the life of patients, users and others when these Instructions for Use are disregarded!

Read Instructions for Use completely and follow them carefully. In particular, be sure to note all cautionary and warning notices.

This document refers to persons of all genders alike. Reference to various genders is avoided purely for reasons of improved readability.

The electronic version of these Instructions for Use can be requested at <https://www.klsmartin.com/en/services/instructions-for-use/>.

1.5 Case Planning

The implants are created following detailed individual case planning by the manufacturer in coordination with a user. Each implant is created exclusively for a previously identified patient known by name.

For case planning, the patient data and/or other case-related information is first uploaded to the IPS Gate® by the user. The IPS Gate® is a web-based platform that guides the user through the process for specially manufactured products in a structured manner.

The requirements, information and patient data provided by the user in IPS Gate® are used by the manufacturer for the individual case planning. A case number is assigned (example C1234567), and the case is thus unmistakably linked to an individual patient.

Once the planning has been completed, the user is given access to the case-specific product documentation (approval document) via IPS Gate®. If this meets his/her requirements, the user must approve the implant for production by submitting the signed document “Written Prescription and Design Approval”.

CAUTION

Possible danger of injury to the patient due to outdated or incorrect patient data!

The specially manufactured implants are produced using imaging procedures, e.g. CT images of the patient. If the patient’s anatomical situation has changed since the creation of the images, the specially manufactured implants must not be used.

1.6 Terms & Acronyms

Acronym	Description
A ₀	Standard for the killing of microorganisms in procedures with humid heat, ISO 15883
ASTM	American Society for Testing and Materials
CT	Computer tomography
EN	European Standard
EtO	Ethylene oxide
IPS®	Individual Patient Solutions
ISO	International Organization for Standardization
MRI	Magnetic resonance imaging
OR	Operating room
PE	Polyethylene
WD	Washer-disinfector
UHMWPE	Ultra-High Molecular Weight Polyethylene
TMJ	Temporomandibular joint

1.7 Validity of this Document

This document applies to:

- IPS Implants® marPOR – Custom-made product: polyethylene implants for craniomaxillofacial surgery

2 Scope of Delivery

The scope of delivery consists of

- the specially manufactured implant or, in some cases, multiple specially manufactured implants
- the corresponding anatomical model

If desired by the user and coordinated in advance as part of individual case planning via IPS Gate®, the scope of delivery optionally also includes custom-made planning aids for a specific patient and/or one or more orthognathic splint(s). Separate Instructions for Use are supplied for planning aids, anatomical models, and orthognathic splints. They are therefore not subject of this document.

Each delivery is accompanied by a delivery note, which contains the case number that was assigned uniquely during case planning via IPS Gate®.

2.1 Inspecting the Delivery for Completeness and Correctness

- Immediately upon receipt, the goods must be checked for completeness and absence of damage.
- After delivery of the product, check the original packaging and the packaging seal for damages.
- Reconcile the delivery with the existing case-specific product documentation from case planning, and check for correct fit. This is ideally done based on the case number. If the case number does not match the case-specific documentation, there has been a mix-up, and the product must not be used.
- Sterile products with a broken seal or damaged sterile packaging must be considered non-sterile and are not to be approved for use.
- The EtO indicator for the products must be checked immediately upon receipt. Products whose EtO indicator has changed must be rejected immediately and may not be used.
- Any transport damage as well as deviations in the case number indicated in the delivery note / on the product packaging of the case-specific product documentation must be reported immediately.

3 Intended Use

3.1 Intended Purpose

The implants are single-use devices intended for augmentation and reconstruction of bone tissue in craniomaxillofacial surgery.

3.2 Indications

- Treatment of craniomaxillofacial trauma
- Correction and reconstruction of craniomaxillofacial defects
- Craniomaxillofacial augmentation
- Plastic and cosmetic surgeries in the craniomaxillofacial area

3.3 Contraindications

Possible danger to the life of the patient due to non-observance of contraindications!

Before surgery, contraindications must be considered and clarified.

- Manifest infection in the vicinity of the location to be treated
- Material sensitivity or allergies to the implant materials used
- Quantitatively and/or qualitatively insufficient hard and/or soft tissue
- Patients who are unwilling or unable to follow instructions during the entire treatment including the post-operative phase
- Use in load-bearing areas such as the temporomandibular joint (TMJ)
- Use in areas exposed to the outside environment, e. g. sinuses
- Use outside the areas defined in the case planning

3.4 Possible Side Effects

Surgical procedures always carry a risk of adverse effects and complications. In most cases, these are not caused by the product used, but rather due to clinical conditions.

The possible general adverse effects that may occur in relation to the use of an implant include:

- Loosening of the implant
- Breaking of the implant
- Material sensitivity or allergic reactions
- Qualitative and/or quantitative impairment of the hard and/or soft tissue
- Infections
- Wound-healing disorders
- Dislocation of the implant
- Contour realignment (Dislocation of the implant can potentially lead to a secondary intervention for contour realignment)
- Unsatisfying aesthetic appearance

The user is responsible to obtain informed consent from the patient with respect to potential adverse effects and complications.

3.5 Target Patient Groups

The implants that are custom-made products are used exclusively for a previously identified patient known by name. They are to be used only according to an obligatory case planning defined before the procedure.

The implants intended for use in patients with a fully developed skeleton.

CAUTION
--

Possible danger of injury to the patient when using the implants in immature skeletal structures!

When implants are used in children or adolescents, whose skeletal structure is not yet fully developed, natural bone growth can lead to misalignment of the implant. The user must evaluate the advantages of the use of implants based on the potential risks for the patient.

3.6 Users

The products must only be used by qualified surgical staff.

The user bears the responsibility for:

- Thorough study of the case-specific product documentation
- Selection and proper execution of the surgical procedure
- Compilation of the product components and their proper use as defined in the case planning
- Avoiding or minimizing the general risks posed by surgical interventions
- Obtaining informed consent from, instructing, and taking post-operative care of the patient

3.7 Ambient Conditions During Application

The application is restricted to the OR under proper surgical conditions. There are no special requirements for temperature and humidity.

3.8 Limitations of Use

⚠ CAUTION

Possible danger of injury to the patient caused by potential misinterpretation of examination results when using CT or MRI!

The implants are made of polyethylene. The material is not magnetic and should not therefore cause disturbing artifacts in computer tomography (CT) and in magnet resonance imaging (MRI).

4 Use / Application

4.1 Description of the Components

4.1.1 Structure, Operating Principle and Performance Specifications

Custom-made products are created based on a detailed individual case planning in coordination with a user and only for a previously identified patient known by name. Each implant is individually designed based on the indication and the patient. Please refer to the case-specific product documentation for detailed information on the structure, operating principle and performance specifications of the implants. In IPS Gate®, the documentation is assigned to the respective case by means of a unique case number. In IPS Gate®, the user can call up the case-specific product documentation at any time.

The user's written prescription and design approval, which must be available prior to each production and delivery of an implant, refers to this case-specific product documentation (approval document). With the written approval, the user confirms that he/she has checked, understood and acknowledged the case-specific product documentation (approval document).

4.1.2 Materials and Implant Materials Used

Designation	Material(s)	Material specification(s)	Standard(s)
Implant	Polyethylene (PE)	Ultra-High Molecular Weight Polyethylene (UHMWPE)	DIN EN ISO 5834-1, ASTM F648

4.1.3 Processing of the Implants

CAUTION

Possible danger of injury to the patient when processing implants!

The size and shape of specially manufactured implants cannot be adjusted, as these are already adapted to the patient's anatomy based on the individual case planning and the anatomical needs of the patient stored therein.

4.1.4 Compatible Products and Accessories

CAUTION

Possible danger of injury to the patient caused by incompatible fixation devices!

The implants may never be combined with incompatible osteosynthesis screws or plates or any other incompatible fixation devices.

⚠ CAUTION**Possible danger of injury to the patient due to implant rupture!**

- When screws are placed too close to the edge of the implant it can lead to intraoperative tear out of the screw hole.
- Before the wound closure make sure the implant is not damaged.

The implants are attached to the patient using osteosynthesis screws from KLS Martin's standard osteosynthesis portfolio. The osteosynthesis screws to be used are defined during the case planning process for an implant in consultation with the user. The user must have these items on hand or order them separately, and they are also supplied with separate Instructions for Use which must be observed.

The implants are delivered together with the associated anatomical model.

If desired by the user and coordinated in advance as part of individual case planning via IPS Gate®, the scope of delivery optionally also includes custom-made planning aids for a specific patient and/or orthognathic splints. Separate Instructions for Use for these products are also supplied and must be complied with.

4.2 Before Use**⚠ CAUTION****Possible danger of injury to the patient due to implant failure!**

Before their use in a surgical environment, the implants must be subjected to a final inspection by a user under all circumstances. This user must be familiar with the entire planning process and the case-specific product documentation.

Polyethylene implants are **delivered sterile**. Products delivered in sterile condition are marked with the corresponding symbol on the label. Before use, the packaging must be checked for the expiration date and absence of damage. Sterile products with opened or damaged packaging or sterile products that have been taken out of packages with passed expiration date must be considered non-sterile. Dispose of non-sterile products.

Prior to implantation, the selected implant should remain in the sterile packaging, and the implant location should be determined and prepared first.

4.3 Intraoperative Approach

WARNING

Possible danger to the life of the patient due to contamination of the implant during handling!

Contamination of the implant with microorganisms and foreign material due to the large surface of porous implants can cause infections.

- Prepare the implantation location before the implant is removed from the sterile packaging and shaped.
 - Implant the implant immediately after removing it from its sterile packaging.
 - Do not place implants on operating cloths, operating clothes or other surfaces that could lead to contamination.
-

WARNING

Possible danger to the life of the patient due to implant exposure caused by incorrect placement!

- Cover implants with periosteum or other soft tissue structures.
-

CAUTION

Possible danger of injury to the patient due to dislocation of the implant!

The implant slipping can result in complications.

- The implants must be actively fixed for cranioplastic applications and in mandibular augmentation.
-

CAUTION

Possible danger of injury to the patient caused by improper handling!

Improper handling of the implant can result in implant failure and injuries. Damage reduces the load-bearing capacity and may lead to breakage of the implant.

CAUTION

Possible danger of injury to the patient due to the fragments produced!

Should the implant break during surgery, all fragments must be removed from the situs. Otherwise, there is a risk of foreign body reactions or migration with resulting soft tissue damage.

The implants are created based on an individual case planning process in coordination with a user and exclusively for a previously identified patient known by name. The products are individually designed based on the indication and the patient. Please refer to the case-specific product documentation for detailed

information on the attachment of the implants. In IPS Gate®, the documentation is assigned to the respective case.

The case-specific product documentation includes the preoperative situation, the planned positioning of the implants, and the planned postoperative situation. Depending on the respective indication, the case-specific product documentation may include further information on the intraoperative procedure.

NOTICE

The products used in the operation must be documented in the patient file with the case-specific data, in particular the case number. Clear traceability of the products used can only be ensured through documentation.

The implant card is / are in the outer packaging.

4.4 After Conclusion of the Surgery

⚠ CAUTION**Possible danger of injury to the patient due to improper patient and user behavior!**

Improper postoperative behavior of the patient can lead to implant failure or injuries of the patient.

Note the following:

- All physical activities that may adversely affect the implant or connected bone elements should be avoided during the healing period. This applies especially to the carrying of heavy loads (danger of falling) and to sporting activities that may involve impacts or blows.
 - The surgeon must advise the patient on post-surgical procedure and carry out follow-up.
-

⚠ CAUTION**Possible danger of injury to the patient due to dislocation of the implant!**

- Perform long-term follow-up after orbital implantation.
-

The implants are intended for single use. Explanted and/or contaminated implants must not be reprocessed and reused. They must be disposed of after completion of the surgical procedure.

5 Cleaning, Disinfection and Sterilization

The implants are supplied in sterile condition. They are intended for single use. Reprocessing and/or resterilization is not possible, as the product properties would be harmed by this process.

WARNING

Possible danger to the life of the patient when reprocessing and reusing single-use products!

The implants are intended for single use in one patient. They are marked with the corresponding symbol on the label and must not be reprocessed or reused.

Even if the single-use product seems intact after use, its integrity may be compromised, and minor defects and/or internal weak points may be present, possibly causing the product to fail.

Single-use products that have not been used but have come into contact with bodily fluids, blood, tissue and/or similar are to be considered used and must be discarded.

Reprocessing and reuse can increase the risk of contamination, e. g. by carryover of germs.

WARNING

Possible danger to the life of the patient when reprocessing contaminated single-use products!

Single-use products that have not been used but have come into contact with bodily fluids, blood, tissue and/or similar are to be considered used and must be discarded.

Reprocessing and reuse can increase the risk of contamination, e. g. by carryover of germs.

5.1 Storage and Transport

- Store and transport cleanly.
- Avoid severe temperature fluctuations.
- Protect the products against mechanical damage.

6 Ecological Information / Disposal

6.1 Packaging

KLS Martin will take back the full packaging upon request. Whenever possible, parts of the packaging will be reused.

When disposing of or recycling the packaging, please take applicable national regulations into account.

6.2 Disposal

In designing the product, we tried to avoid using composite materials wherever possible. This allows a high degree of recycling. We therefore offer you to take the product back for proper disposal and recycling.

6.3 National Regulations

The national regulations and disposal provisions must be observed for all disposal measures.

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