Deutsch
English
Español
Français
Italiano

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Table of Contents

1	General Information.....	24
1.1	Manufacturer	24
1.2	Hotline	24
1.3	Event Reporting Duty	24
1.4	Notes on this Document	25
1.5	Case Planning	25
1.6	Terms & Acronyms	26
1.7	Validity of this Document.....	26
2	Scope of Delivery	26
2.1	Inspecting the Delivery for Completeness and Correctness	26
3	Intended Use	27
3.1	Intended Purpose	27
3.2	Indications	27
3.3	Contraindications	27
3.4	Possible Side Effects	28
3.5	Target Patient Groups	28
3.6	Users.....	29
3.7	Ambient Conditions During Application.....	29
3.8	Limitations of Use.....	29
4	Use / Application	29
4.1	Description of the Components	29
4.1.1	Structure, Operating Principle and Performance Specifications	29
4.1.2	Materials and Implant Materials Used	30
4.1.3	Processing of the Implants.....	30
4.1.4	Compatible Products and Accessories.....	30
4.2	Before Use	31
4.3	Intraoperative Approach	31
4.4	After Conclusion of the Surgery	32
5	Cleaning, Disinfection and Sterilization	32
5.1	Limitations and Restrictions to Cleaning, Disinfecting and Sterilizing	32
5.2	Preparation before Cleaning	33
5.3	Cleaning and Disinfection.....	33
5.3.1	Machine Cleaning and Disinfection	33
5.4	Checks and Maintenance	35
5.4.1	Maintenance	35
5.5	Packaging.....	35

5.6 Sterilization35

5.7 Storage and Transport36

6 Ecological Information / Disposal 37

6.1 Packaging37

6.2 Disposal37

6.3 National Regulations.....37

1 General Information

1.1 Manufacturer

We thank you for buying one of our products.

The product described in these Instructions for Use was specifically manufactured in accordance with a written prescription for an identified patient known by name and is intended solely to suit the patient's individual condition and needs. This product is therefore a custom-made product within the meaning of the applicable European regulations for the development, manufacture, placing on the market and monitoring of medical devices and as such does not bear any CE marking. As a custom-made product, it nevertheless fulfills the basic requirements for the safety and performance of medical devices within the meaning of the applicable European regulations.

We are the manufacturer of this product:

KLS Martin SE & Co. KG

A company of the KLS Martin Group

KLS Martin Platz 1 · D-78532 Tuttlingen · Germany

Phone +49 7461 706-0 · Fax +49 7461 706-193

info@klsmartin.com www.klsmartin.com

1.2 Hotline

Should you have any questions on how to handle the device or product or use it for clinical applications, please do not hesitate to contact the Product Management:

Phone: +49 7461 706-0

E-mail: info@klsmartin.com

NOTICE

Each custom-made implant and its packaging are marked with an individual case number. Please always provide the case number in case of a complaint.

1.3 Event Reporting Duty

WARNING

Possible danger to the life of others by shipping contaminated products!

In case of return shipments, send only clean and disinfected products in sterile packaging.

All serious adverse events occurring in connection with the product must be reported to KLS Martin and the respective authorities without delay.

1.4 Notes on this Document

WARNING

Possible danger to the life of patients, users and others when these Instructions for Use are disregarded!

Read Instructions for Use completely and follow them carefully. In particular, be sure to note all cautionary and warning notices.

This document refers to persons of all genders alike. Reference to various genders is avoided purely for reasons of improved readability.

The electronic version of these Instructions for Use can be requested at <https://www.klsmartin.com/en/services/instructions-for-use/>.

1.5 Case Planning

The implants are created following detailed individual case planning by the manufacturer in coordination with a user. Each implant is created exclusively for a previously identified patient known by name.

For case planning, the patient data and/or other case-related information is first uploaded to the IPS Gate® by the user. The IPS Gate® is a web-based platform that guides the user through the process for custom-made products in a structured manner.

The requirements, information and patient data provided by the user in IPS Gate® are used by the manufacturer for the individual case planning. A case number is assigned (example C1234567), and the case is thus unmistakably linked to an individual patient.

Once the planning has been completed, the user is given access to the case-specific product documentation (approval document) via IPS Gate®. If this meets his/her requirements, the user must approve the implant for production by submitting the signed document “Written Prescription and Design Approval”.

CAUTION

Possible danger of injury to the patient due to outdated or incorrect patient data!

The custom-made implants are produced using imaging procedures, e.g. CT images of the patient. If the patient’s anatomical situation has changed since the creation of the images, the custom-made implants must not be used.

1.6 Terms & Acronyms

Acronym	Description
A ₀	Standard for the killing of microorganisms in procedures with humid heat, ISO 15883
ASTM	American Society for Testing and Materials
CT	Computer tomography
EN	European Standard
IPS®	Individual Patient Solutions
ISO	International Organization for Standardization
MRI	Magnetic resonance imaging
OR	Operating room
PEEK	Polyether ether ketone
WD	Washer-disinfector

1.7 Validity of this Document

This document applies to:

- IPS Implants® PEEK – Custom-made product: Polyether ether ketone implants for craniomaxillofacial surgery

2 Scope of Delivery

The scope of delivery consists of

- One custom-made implant or, in some cases, multiple custom-made implants
- the corresponding anatomical model

If desired by the user and coordinated in advance as part of individual case planning via IPS Gate®, the scope of delivery optionally also includes custom-made planning aids for a specific patient and/or one or more orthognathic splint(s). Separate Instructions for Use are supplied for planning aids, anatomical models, and orthognathic splints. They are therefore not subject of this document.

Each delivery is accompanied by a delivery note, which contains the case number that was assigned uniquely during case planning via IPS Gate®.

2.1 Inspecting the Delivery for Completeness and Correctness

- Immediately upon receipt, the goods must be checked for completeness and absence of damage.
- After delivery of the product, check the original packaging and the packaging seal for damages.

- Reconcile the delivery with the existing case-specific product documentation from case planning, and check for correct fit. This is ideally done based on the case number. If the case number does not match the case-specific documentation, there has been a mix-up, and the product must not be used.
- Any transport damage as well as deviations in the case number indicated in the delivery note / on the product packaging of the case-specific product documentation must be reported immediately.

3 Intended Use

3.1 Intended Purpose

The implants are single-use devices intended for augmentation and reconstruction of bone tissue in craniomaxillofacial surgery.

3.2 Indications

- Reconstruction of bone structures after destruction/deformation by, for example, tumors, ulcers, cysts, traumas, infections or rejection reactions
- Reconstruction in cases of congenital deformations
- Orthognathic surgeries
- Cerebrocranial traumata with elevated intracranial pressure
- Bone augmentation
- Cosmetic surgery

3.3 Contraindications

Possible danger to the life of the patient due to non-observance of contraindications!

Before surgery, contraindications must be considered and clarified.

-
- Manifest infection in the vicinity of the location to be treated
 - Material sensitivity or allergies to the implant materials used
 - Quantitatively and/or qualitatively insufficient hard and/or soft tissue
 - Patients who are unwilling or unable to follow instructions during the entire treatment including the post-operative phase
 - Use outside the areas defined in the case planning

3.4 Possible Side Effects

Surgical procedures always carry a risk of adverse effects and complications. In most cases, these are not caused by the product used, but rather due to clinical conditions.

The possible general adverse effects that may occur in relation to the use of an implant include:

- Loosening of the implant
- Breaking of the implant
- Material sensitivity or allergic reactions
- Qualitative and/or quantitative impairment of the hard and/or soft tissue
- Infections
- Wound-healing disorders
- Dislocation of the implant

The user is responsible to obtain informed consent from the patient with respect to potential adverse effects and complications.

3.5 Target Patient Groups

The implants that are custom-made products are used exclusively for a previously identified patient known by name. They are to be used only according to an obligatory case planning defined before the procedure.

The implants intended for use in patients with a fully developed skeleton.

⚠ CAUTION**Possible danger of injury to the patient when using the implants in immature skeletal structures!**

When implants are used in children or adolescents, whose skeletal structure is not yet fully developed, natural bone growth can lead to misalignment of the implant. The user must evaluate the advantages of the use of implants based on the potential risks for the patient.

3.6 Users

The products must only be used by qualified surgical staff.

The cleaning, disinfecting and sterilizing is to be performed by trained and qualified staff of the medical device preparation unit.

The user bears the responsibility for:

- Thorough study of the case-specific product documentation
- Selection and proper execution of the surgical procedure
- Compilation of the product components and their proper use as defined in the case planning
- Avoiding or minimizing the general risks posed by surgical interventions
- Obtaining informed consent from, instructing, and taking post-operative care of the patient

3.7 Ambient Conditions During Application

The application is restricted to the OR under proper surgical conditions.

There are no special requirements for temperature and humidity.

3.8 Limitations of Use

CAUTION

Possible danger of injury to the patient caused by potential misinterpretation of examination results when using CT or MRI!

The implants are made of PEEK. The material is not magnetic and should not therefore cause disturbing artifacts in computer tomography (CT) and in magnet resonance imaging (MRI).

4 Use / Application

4.1 Description of the Components

4.1.1 Structure, Operating Principle and Performance Specifications

Custom-made products are created based on a detailed individual case planning in coordination with a user and only for a previously identified patient known by name. Each implant is individually designed based on the indication and the patient. Please refer to the case-specific product documentation for detailed information on the structure, operating principle and performance specifications of the implants. In IPS Gate®, the documentation is assigned to the respective case by means of a unique case number. In IPS Gate®, the user can call up the case-specific product documentation at any time.

The user's written prescription and design approval, which must be available prior to each production and delivery of an implant, refers to this case-specific product documentation (approval document). With the

written approval, the user confirms that he/she has checked, understood and acknowledged the case-specific product documentation (approval document).

4.1.2 Materials and Implant Materials Used

Designation	Material	Material specification	Standard(s)
Implant	Polyether ether ketone (PEEK)	--	ASTM F 2026

4.1.3 Processing of the Implants

CAUTION

Possible danger of injury to the patient when processing implants!

The size and shape of custom-made implants cannot be adjusted, as these are already adapted to the patient's anatomy based on the individual case planning and the anatomical needs of the patient stored therein.

4.1.4 Compatible Products and Accessories

CAUTION

Possible danger of injury to the patient caused by incompatible fixation devices!

The implants may never be combined with incompatible osteosynthesis screws or plates or any other incompatible fixation devices.

CAUTION

Possible danger of injury to the patient due to implant rupture!

When fixing the implant, ensure that it is not damaged.

The implants are attached to the patient using osteosynthesis screws from KLS Martin's standard osteosynthesis portfolio. The osteosynthesis screws to be used are defined during the case planning process for an implant in consultation with the user. The user must have these items on hand or order them separately, and they are also supplied with separate Instructions for Use which must be observed.

The implants are delivered together with the associated anatomical model.

Optionally – if desired by the user and coordinated in advance as part of the individual case planning via IPS Gate® – the scope of delivery also includes custom-made planning aid(s) for a specific patient. Separate Instructions for Use for these products are also supplied and must be complied with.

4.2 Before Use

The implants are delivered in non-sterile condition. They must be removed from their packaging before use and subjected to a suitable cleaning, disinfection and sterilization process.

CAUTION

Possible danger of injury to the patient due to implant failure!

Before their use in a surgical environment, the implants must be subjected to a final inspection by a user under all circumstances. This user must be familiar with the entire planning process and the case-specific product documentation.

4.3 Intraoperative Approach

CAUTION

Possible danger of injury to the patient caused by improper handling!

Improper handling of the implant can result in implant failure and injuries. Damage reduces the load-bearing capacity and may lead to breakage of the implant.

CAUTION

Possible danger of injury to the patient caused by produced fragments!

Should the implant break during surgery, all fragments must be removed from the situs. Otherwise, there is a risk of foreign body reactions or migration with resulting soft tissue damage.

The implants are created based on an individual case planning process in coordination with a user and exclusively for a previously identified patient known by name. The products are individually designed based on the indication and the patient. Please refer to the case-specific product documentation for detailed information on the attachment of the implants. In IPS Gate®, the documentation is assigned to the respective case.

The case-specific product documentation includes the preoperative situation, the planned positioning of the implants, and the planned postoperative situation. Depending on the respective indication, the case-specific product documentation may include further information on the intraoperative procedure.

NOTICE

The products used in the operation must be documented in the patient file with the case-specific data, in particular the case number. Clear traceability of the products used can only be ensured through documentation.

The implant card is / are in the outer packaging.

4.4 After Conclusion of the Surgery

CAUTION

Possible danger of injury to the patient due to improper patient and user behavior!

Improper postoperative behavior of the patient can lead to implant failure or injuries of the patient. Note the following:

- All physical activities that may adversely affect the implant or connected bone elements should be avoided during the healing period. This applies especially to the carrying of heavy loads (danger of falling) and to sporting activities that may involve impacts or blows.
- The surgeon must advise the patient on post-surgical procedure and carry out follow-up.

The implants are intended for single use. Explanted and/or contaminated implants must not be reprocessed and reused. They must be disposed of after completion of the surgical procedure.

5 Cleaning, Disinfection and Sterilization

WARNING

Non-sterile handling can result in danger to the life of the patient!

Improper sterilization and non-sterile handling of the implants can pose serious health hazards to patients.

- The operator/preparer is responsible for the cleaning, disinfecting and sterilizing of the used products. It is essential that national regulations, including restrictions, be observed.

The preparation process must be validated before the medical device is used. The operator / preparer is responsible for this.

5.1 Limitations and Restrictions to Cleaning, Disinfecting and Sterilizing

WARNING

Possible danger to the life of the patient due to cleaning, disinfection and sterilization of contaminated single-use products!

Single-use products that have not been used but have come into contact with bodily fluids, blood, tissue and/or similar are to be considered used and must be discarded.

Even if the single-use product seems intact after use, its integrity may be compromised, and minor defects and/or internal weak points may be present, possibly causing the product to fail.

Cleaning, disinfecting, sterilization and use can actually increase the contamination risk, e.g. due to germ transfer.

There must be no mechanical load on the products throughout the entire cleaning and sterilization cycle (e.g. pressure caused by stacking).

If a non-sterile delivered implant has already been cleaned, disinfected and sterilized but could not be used directly (e.g. because the surgical procedure had to be postponed), it can be cleaned, disinfected and sterilized once again as described below.

5.2 Preparation before Cleaning

NOTICE

Pre-cleaning is not performed for IPS implants.

5.3 Cleaning and Disinfection

5.3.1 Machine Cleaning and Disinfection

When choosing a WD, it must be ensured that,

- the WD efficiency is proven (e.g. DGHM or FDA approval/clearance/registration or CE marking according to DIN EN ISO 15883),
- if possible, a tested program for thermal disinfection (A_0 -value > 3000 or – for older devices – at least 5 min. at 93 °C/194.4 °F) is used (with chemical disinfection there is a risk of disinfectant residues on the products),
- the used program is suitable for the products and has enough rinsing cycles.

Detergent residues can effectively be prevented with any of the following measures:

- a) performing at least three diluting steps after cleaning or neutralization or
- b) conductance control of the WD
- using only sterile or low-germ (max. 10 germs/ml) as well as low endotoxin load (max. 0.25 endotoxin units /ml) water (e.g. purified water/highly purified water)¹,
- ensuring that the air used for drying is filtered (oil-free, low germ- and particle load) and
- ensuring that the WD is regularly serviced, inspected and calibrated.

When choosing the cleaning agent system it must be ensured,

- that it is generally suitable for the cleaning of invasive medical devices made from metals and plastics,
- that, when no thermal disinfection is applied, an additional suitable disinfectant with tested efficacy (e.g. VAH/DGHM or FDA/EPA approval/clearance/registration or CE marking) is used and that it is compatible with the used cleaning agent used. In general, the manufacturer recommends thermal disinfection to keep the chemical load on the implants as low as possible.

¹ If a lower water quality should be deemed sufficient in light of national recommendations (e.g. in Germany KRINKO/RKI/BfArM recommendation for preparation), the responsibility falls solely on the preparer.

- and that the used chemicals are compatible with the product.

When choosing the cleaning agents and disinfectants, it must be ensured that they do **not** contain the following ingredients:

- organic, mineral and oxidizing acids (minimum permissible pH value 5.5)
- organic solvents (e.g. alcohols, ether, ketones, benzines)
- oxidants (e.g. hydrogen peroxide)
- halogens (chloride, iodine, bromine)
- aromatic/halogenated hydrocarbons

The manufacturer instructions for the cleaning agent regarding concentrations, temperature, application time as well as rinsing specifications must be observed.

Procedure:

1. Place the products into the WD. Ensure that products do not touch each other.
2. Start program.
3. Remove products from the WD after the end of program.
4. Inspect and package products as soon as possible after removal.
5. Ensure that suitable small-mesh sieve trays are used for implants (e.g., plates, screws, distractors) so that damage due to protruding tips and bending or scratching of the medical devices can be ruled out. In addition, this prevents implants that fall out from dropping into the WD and damaging the latter.

The following cleaning process has been validated:

Phase	Step	Temperature	Time	Water quality	Additional Information
1	Pre-cleaning	cold (not temperature-controlled)	1 min.	Drinking water	---
2	Cleaning	50 °C	5 min.	Drinking water	Minimum concentration according to the specifications of the cleaning agent manufacturer
3	Neutralization	---	---	---	Neutralization not performed according to information provided by the cleaning agent manufacturer
4	Interim rinsing	cold (not temperature-controlled)	1 min.	Demineralized water	---
5	Thermal disinfection	90 °C	5 min.	Demineralized water	A ₀ -value ≥ 3000
6	Drying	according to the program	according to the program	---	Drying was not used in the validation

The proof of the general suitability of the products for an effective machine cleaning and disinfection was provided by an independent, accredited and recognized test lab using the WD Miele PG 8353 (thermal disinfection, Miele & Cie. GmbH & Co., Gütersloh) and the pre-cleaning and cleaning agent neodisher MediClean forte (Dr. Weigert GmbH & Co. KG, Hamburg). The above-mentioned procedure was observed in doing so.

5.4 Checks and Maintenance

NOTICE

Further use is confirmed by the successful inspection of the product. The inspection release and the packaging into a sterile barrier system release the product for its use.

The products must be macroscopically clean, i.e., free from visible contamination, after the cleaning process.

5.4.1 Maintenance

Implants must under no circumstances be treated with care agents.

5.5 Packaging

Approved sterilization packaging (e.g. conforming to EN 868, ISO 11607) must be used for sterilization, subsequent transportation and storage.

5.6 Sterilization

Sterilization is to be done applying the following steam sterilization method. Other sterilization methods are not permissible.

Steam sterilization:

- Fractionated vacuum procedure (with sufficient product drying²)
- Steam sterilizer according to DIN EN 13060/DIN EN 285 or ANSI AAMI ST79 (for USA: FDA-Clearance)
- Must be validated according to DIN EN ISO 17665 (valid IQ/OQ (commissioning) and product-specific performance evaluation (PQ))
- Maximum sterilization temperature 134 °C (273 °F, plus tolerance according to DIN EN ISO 17665)
- The procedure has been validated at 132 °C with a holding time of 2 min (half-cycle procedure), so the following parameters are permissible:

Steam sterilization with fractionated pre-vacuum		
	Procedure 1	Procedure 2
Pre-vacuum cycles	min. 3 ×	min. 3 ×
Medium	saturated water steam	saturated water steam
Sterilization time	4 min. (or longer) ³	4 min. (or longer) ³
Sterilization temperature	132 °C, plus tolerance	134 °C, plus tolerance
Drying time	20 min ²	20 min ²

Evidence of the basic suitability of the devices for effective steam sterilization was provided by an independent, officially accredited and recognized test laboratory using the steam sterilizer Lautenschläger ZentraCert and the fractionated vacuum procedure. In doing so, the typical conditions at hospitals and medical offices as well as the above-mentioned procedure were taken into consideration.

Do not use gravitation sterilization, no flash sterilization, no hot air sterilization, no radiation sterilization, no formaldehyde or ethylene dioxide sterilization and also no plasma sterilization.

5.7 Storage and Transport

Possible danger to the life of others by shipping contaminated products!

In case of return shipments, send only clean and disinfected products in sterile packaging.

- Store and transport cleanly.
- Avoid severe temperature fluctuations.
- Protect the products against mechanical damage.

² The actual required drying time depends directly on the parameters, which are the sole responsibility of the user (loading configuration and density, sterilizer condition, ...), and must therefore be determined by the user. Drying times should, nevertheless, not be less than 20 min.

³ or extended sterilization time (e.g. 18 min.) for the prion activation according to national standards.

6 Ecological Information / Disposal

6.1 Packaging

KLS Martin will take back the full packaging upon request. Whenever possible, parts of the packaging will be reused.

When disposing of or recycling the packaging, please take applicable national regulations into account.

6.2 Disposal

In designing the product, we tried to avoid using composite materials wherever possible. This allows a high degree of recycling. We therefore offer you to take the product back for proper disposal and recycling.

6.3 National Regulations

The national regulations and disposal provisions must be observed for all disposal measures.

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