
Deutsch
English
Español
Français
Italiano

IPS Implants® Resorb Mg – Sonderanfertigung:
Resorbierbare Magnesium-Implantate

IPS Implants® Resorb Mg – custom-made:
Absorbable magnesium implants

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Table of Contents

1	General Information.....	23
1.1	Manufacturer	23
1.2	Hotline	23
1.3	Case Planning	24
1.4	Event Reporting Duty	24
1.5	Notes on this Document	24
1.6	Terms & Acronyms	25
1.7	Validity of this Document.....	25
2	Scope of Delivery	26
2.1	Inspecting the Delivery for Completeness and Correctness	26
3	Intended Use	27
3.1	Intended Purpose	27
3.2	Indications	27
3.3	Contraindications	27
3.4	Possible Side Effects	28
3.5	Target Patient Groups	28
3.6	Users.....	28
3.7	Ambient Conditions During Application.....	29
3.8	Warnings	29
3.9	Safety Instructions.....	29
4	Use / Application	30
4.1	Description of the Components	30
4.1.1	Structure, Operating Principle and Performance Specifications	30
4.1.2	Materials and Implant Materials Used	30
4.1.3	Processing of the Implants.....	30
4.1.4	Compatible Products and Accessories	31
4.2	Before Implantation	31
4.3	Intraoperative Approach	32
4.4	After Conclusion of the Surgery	33
5	Cleaning, Disinfection, and Sterilization	34
6	Storage and Transport.....	34
7	Ecological Information / Disposal.....	35
7.1	Packaging.....	35
7.2	Disposal	35
7.3	National Regulations	35

1 General Information

1.1 Manufacturer

Thank you for choosing a KLS Martin IPS Implants® product.

The product described in these Instructions for Use was specifically manufactured in accordance with a written prescription for an identified patient known by name and is intended solely to suit the patient's individual condition and needs. This product is therefore a custom-made product within the meaning of the applicable European regulations for the development, manufacture, placing on the market and monitoring of medical devices and as such does not bear any CE marking. As a custom-made product, it nevertheless fulfills the basic requirements for the safety and performance of medical devices within the meaning of the applicable European regulations.

We are the manufacturer of this product:

KLS Martin SE & Co. KG

A company of the KLS Martin Group

KLS Martin Platz 1 · D-78532 Tuttlingen · Germany

Phone +49 7461 706-0 · Fax +49 7461 706-193

info@klsmartin.com · www.klsmartin.com

1.2 Hotline

Should you have any questions on how to handle the device or product or use it for clinical applications, please do not hesitate to contact the Product Management:

Phone: +49 7461 706-0

E-mail: info@klsmartin.com

NOTICE

All packaging, and in some cases also the product itself, is marked with case-specific data. Please always provide the case number in case of a complaint.

1.3 Case Planning

The implants that are custom-made products are created following detailed individual case planning by the manufacturer in coordination with a qualified user. Each implant is created exclusively for a previously identified patient known by name.

For case planning, the patient data and/or other case-related information is first uploaded to the IPS Gate® by the qualified user. The IPS Gate® is a web-based platform that guides the user through the process for implants in a structured manner.

The requirements, information and patient data provided by the qualified user in IPS Gate® are used by the manufacturer for the individual case planning. A case number is assigned (example C1234567), and the case is thus unmistakably linked to an individual patient.

Once the planning has been completed, the user is given access to the case-specific product documentation (approval document) via IPS Gate®. If this meets his/her requirements, the user must approve the implant for production by submitting the signed document "Written Prescription and Design Approval".

CAUTION**Possible danger of injury of the patient due to outdated or incorrect patient data!**

The implants are produced using imaging procedures, e.g. CT images of the patient. If the patient's anatomical situation has changed since the creation of the images, the implants must not be used.

1.4 Event Reporting Duty

All serious adverse events occurring in connection with the product must be reported to KLS Martin and the respective authorities without delay.

1.5 Notes on this Document

WARNING**Possible danger to the life of patients, users and others when these Instructions for Use are disregarded!**

Read Instructions for Use completely and follow them carefully. In particular, be sure to note all cautionary, warning notices.

This document refers to persons of all genders alike. Reference to various genders is avoided purely for reasons of improved readability.

The electronic version of these Instructions for Use can be requested at <https://www.klsmartin.com/en/services/instructions-for-use/>.

1.6 Terms & Acronyms

Acronym	Description
ASTM	American Society for Testing and Materials
CT	Computer tomography
IPS®	Individual Patient Solutions
Mg	Magnesium
MRI	Magnetic resonance imaging
OR	Operation
REF	Catalogue number

1.7 Validity of this Document

This document applies to:

- IPS Implants® Resorb Mg – custom-made: Absorbable magnesium implants

2 Scope of Delivery

The product consists of:

- one absorbable implant or, in some cases, multiple absorbable implants.
- the corresponding anatomical model

If desired by the user and coordinated in advance via IPS Gate® as part of the individual case planning, the scope of the delivery may also optionally include custom-made planning aids for a specific patient. Planning aids and anatomical models are supplied with separate Instructions for Use. They are therefore not subject of this document.

Each delivery is accompanied by a delivery note, which contains the case number that was assigned uniquely during case planning via IPS Gate®.

2.1 Inspecting the Delivery for Completeness and Correctness

- Immediately upon receipt, the goods must be checked for completeness and damages.
- Reconcile the delivery with the existing case-specific product documentation from case planning, and check for correct fit. This is ideally done based on the case number. If the case number does not match the case-specific documentation, there has been a mix-up, and the product must not be used.
- Any transport damage as well as deviations in the case number indicated in the delivery note / on the product packaging of the case-specific product documentation must be reported immediately.
- Check the original packaging in particular and the packaging seal if applicable for intactness. Also check the protective film in the case of sterile packed products.

Implants in sterile packaging can be applied directly from their sterile packaging without further pre-treatment.

Sterile implants with opened (even if accidental) or damaged sterile packaging or implants that have been removed from packaging with a lapsed expiration date must be considered non-sterile and therefore must be discarded.

The delivery must agree with the user's documents from case planning. Use the case number in particular for inspection. If this data does not match the available documents, there has been a mix-up and the product may not be used.

The usage of implants as well as the clinical planning process are the responsibility of the user. The product must be designed/implemented in accordance with the current state of the art, e.g. material thickness. Special cases deviating from this rule, especially when due to anatomical conditions, must be approved explicitly in advance.

3 Intended Use

The Absorbable custom-made implant must be used in accordance with its indication, see chapter 3.2.

The manufacturer assumes no responsibility for complications caused by improper use and/or use outside the indication, incorrect combination of system components, surgical technique, incorrect treatment method, or insufficient hygiene.

3.1 Intended Purpose

The implants are single-use devices intended for alignment, reconstruction and stabilization of bone tissue. They are wholly or mainly absorbed.

3.2 Indications

- Reconstruction of bone structures in the craniomaxillofacial area after destruction/deformation by, for example, tumors, ulcers, cysts, traumas, infections or rejection reactions
- Reconstruction in cases of congenital deformations in the craniomaxillofacial area
- Orthognathic surgeries
- Cerebrocranial traumata with elevated intracranial pressure
- Bone augmentation in the craniomaxillofacial area
- Cosmetic surgery in the craniomaxillofacial area

3.3 Contraindications

Possible danger to life of the patient in case of non-observance of contraindications!

Bear the contraindications in mind.

- Manifest infection in the vicinity of the location to be treated
- Material sensitivity or allergies to the implant materials used
- Qualitatively and/or quantitatively insufficient hard and/or soft tissue
- Patients who are unwilling or unable to follow instructions during the entire treatment including the post-operative phase
- Use outside the areas defined in the case planning

3.4 Possible Side Effects

Surgical procedures always carry a risk of adverse effects and complications. In most cases, these are not caused by the product used, but rather due to clinical conditions.

The possible general adverse effects that may occur in connection with the application of the osteosynthesis system include:

- Loosening of the implant
- Breaking of the implant
- Material sensitivity or allergic reactions
- Qualitative and/or quantitative impairment of the hard and/or soft tissue
- Infections
- Wound-healing disorders
- Dislocation of the implant
- Unsatisfying aesthetic appearance
- Signs of osteolysis in the x-ray during the resorption process, which is generally not clinically relevant.

The user is responsible to obtain informed consent from the patient with respect to potential adverse effects and complications.

3.5 Target Patient Groups

The implants that are custom-made products are used exclusively for a previously identified patient known by name. They are to be used only according to an obligatory case planning defined before the procedure.

CAUTION

Possible danger of injury to the patient resulting from use of the implants in skeletal structures which are not fully grown!

Use of IPS Implants® Resorb Mg is not recommended in children or adolescents whose skeletal structure is not yet fully developed, as natural bone growth can lead to misalignment of the implant. The user must evaluate the advantages of using IPS Implants® Resorb Mg based on the potential risks for the patient.

3.6 Users

The products may be used only by expert medical staff who are qualified to use the products based on their specialized training, their knowledge and experience, as well as their knowledge of the relevant regulations.

Users must have adequate knowledge of the rules of medical practice, as well as of the state of the art and of the contents of pertinent scientific publications by medical authors.

The user must know and have mastered the surgical technique and the product as a whole.

The user bears the responsibility for:

- Thorough study of the complete product documentation
- Practical exercises in the handling of the product
- Selection and proper execution of the surgical procedure
- Compilation of the product components and their proper use
- Avoiding or minimizing the general risks posed by surgical interventions
- Obtaining informed consent from, instructing, and taking post-operative care of the patient

3.7 Ambient Conditions During Application

The application is restricted to the OR under proper surgical conditions.

There are no special requirements for temperature and humidity.

3.8 Warnings

If implants made of various materials (steel, titanium and alloys of cobalt and chromium) are used together with IPS Implants® Resorb Mg at the operation site, it is necessary to ensure that they are not in direct contact for an extended period of time (physical contact between the implants).

3.9 Safety Instructions

 CAUTION

Possible danger of injury of the patient caused by burns or unintentional movement of the implant as well as possible misinterpretation of the examination upon use of MRI!

Due to the continuous development and increasing energy density of MRI devices, a negative effect on the implants and patients cannot be ruled out for the future.

In case of doubt, MRI examination is not allowed before damage to the health of the patient has been excluded.

Possible dangers are heating or migration of implants possibly reacting to the magnetic field. Metal implants can furthermore cause disturbing artifacts in computer tomography (CT) as well as in magnet resonance imaging (MRI).

Light zones around the implant that are visible on radiological images may occur in individual cases. Based on experience to date, this has no impact on the success of healing. Including this possible event in the surgical report and the physician's letter is recommended. It is generally not necessary to remove the implant as IPS Implants® Resorb Mg are resorbed in the body and replaced with the body's own tissue.

4 Use / Application

4.1 Description of the Components

4.1.1 Structure, Operating Principle and Performance Specifications

Custom-made products are created based on a detailed individual case planning in coordination with a qualified user and only for a previously identified patient known by name.

For case planning, the patient data and/or other case-related information is first uploaded to the IPS Gate® by the qualified user. The IPS Gate® is a web-based platform that guides the user through the process for implants in a structured manner.

The requirements, information and patient data provided by the qualified user in IPS Gate® are used by the manufacturer for the individual case planning. A case number is assigned (example C1234567), and the case is thus unmistakably linked to an individual patient.

Once the planning has been completed, the user is given access to the case-specific product documentation (approval document) via IPS Gate®. If this meets his/her requirements, the user must approve the implant for production by submitting the signed document "Written Prescription and Design Approval".

Each implant is individually designed based on the indication and the patient. Please refer to the case-specific product documentation for detailed information on the structure, functionality and performance characteristics of the implants. In IPS Gate®, the documentation is assigned to the respective case by means of a unique case number. In IPS Gate®, the user can call up the case-specific product documentation at any time.

The user's written prescription and design approval, which must be available prior to each production and delivery of an implant, refers to this case-specific product documentation (approval document). With the written approval, the user confirms that he/she has checked, understood and acknowledged the case-specific product documentation (approval document).

4.1.2 Materials and Implant Materials Used

Material	Material specification	Standard
Magnesium alloy	WE43 MEO	ASTM F 1088

Owing to the underlying absorbable properties, IPS Implants® Resorb Mg are biodegraded in the patient's own body and excreted. Hence, these implants are single-use only.

4.1.3 Processing of the Implants

The size and shape of custom-made implants cannot be adjusted, as these are already adapted to the patient's anatomy based on the individual case planning and the anatomical needs of the patient stored therein.

4.1.4 Compatible Products and Accessories

** CAUTION****Possible danger of injury to the patient due to incompatible fastening elements!**

Never combine the implants with incompatible osteosynthesis screws or any other incompatible attachment elements.

** CAUTION****Possible danger of injury to the patient due to the implant breaking!**

Before the wound closure make sure the implant is not damaged.

IPS Implants® Resorb Mg are fixed to the patient using specially manufactured osteosynthesis screws by KLS Martin.

The fastening element to be used is defined during the case planning process for a specially manufactured implant in consultation with the user. The user must have these items on hand or order them separately, and they are also supplied with separate Instructions for Use which must be observed.

** CAUTION****Possible danger of injury to the patient due to the implant breaking!**

When fixing the implant, ensure that it is not damaged.

The special material properties of IPS Implants® Resorb Mg must be taken into account when inserting the osteosynthesis screws. Specifically, this means that when screwing the osteosynthesis screws tight it is always necessary to ensure that the magnesium implant is fixed to the bone without damaging the implant or the screws.

4.2 Before Implantation

** CAUTION****Possible danger of injury to the patient due to implant failure!**

Before their use in a surgical environment, the absorbable magnesium implants must be subjected to a final inspection by a qualified user under all circumstances. This user must be familiar with the entire planning process and the case-specific product documentation.

IPS Implants® Resorb Mg are supplied in **sterile condition**. Products delivered in sterile condition are marked with the corresponding symbol on the label. They must be removed from the packaging only immediately before use. Before use, the packaging must be checked for the expiration data and for absence of damage. Sterile products with opened or damaged packaging or sterile products that have been taken out of packages with passed expiration date must be considered non-sterile. They must be discarded.

Prior to implantation, the selected implant should remain in the sterile packaging, and the implant location should be determined and prepared first.

4.3 Intraoperative Approach

** WARNING****Possible danger to life of the patient due to implant loosening caused by incorrect positioning!**

Cover implants with periosteum or other soft tissue structures.

** CAUTION****Possible danger of injury to the patient due to improper handling!**

Improper handling of the implant can result in implant failure and injuries.

Damage reduces the load-bearing capacity and may lead to breakage of the implant.

** CAUTION****Possible danger of injury to the patient due to fragment formation!**

Should the implant break during surgery, all fragments must be removed from the situs. Otherwise, there is a risk of foreign body reactions or migration with resulting soft tissue damage.

Implants that are custom-made products are created following individual case planning in coordination with a qualified user. Each implant is created exclusively for a previously identified patient known by name. The products are individually designed based on the indication and the patient. Please refer to the case-specific product documentation for detailed information on the attachment of the implants. In IPS Gate®, the documentation is assigned to the respective case.

The case-specific product documentation includes the preoperative situation, the planned positioning of the implants, and the planned postoperative situation. Depending on the respective indication, the case-specific product documentation may include further information on the intraoperative procedure.

The products used in the operation must be documented in the patient file with the case-specific data, in particular the case number. Clear traceability of the products used can only be ensured through documentation.

** WARNING****Possible danger to life of the patient from non-sterile handling!**

Non-sterile handling of the implants can pose serious health hazards to patients.

Absorbable implants labeled with the marking **STERILE R** are supplied in sterile condition (gamma irradiated) and can be administered directly from the packaging without further pre-treatment.

In this situation, remove the implant from the secondary packaging outside of the sterile zone of the operating room. The primary packaging ensures the sterility of the implant in the sterile zone.

4.4 After Conclusion of the Surgery

CAUTION

Possible danger of injury to the patient due to improper behavior by the patient and user!

Improper postoperative behavior of the patient can lead to implant failure or injuries of the patient. Note the following:

- All physical activities that may adversely affect the implant or connected bone elements should be avoided during the healing period. This applies especially to the carrying of heavy loads (danger of falling) and to sporting activities that may involve impacts or blows.
- The surgeon must advise the patient on post-surgical procedure and carry out follow-up.
- Different absorption behaviors and times need to be assumed depending on the patient and their individual metabolism, age, etc. The user and the patient should be informed about the effects of deviating absorption.
- It is generally not necessary to remove the implant as the absorbable implant is biologically absorbed in the body and replaced with the body's own tissue.

IPS Implants® Resorb Mg are intended for single use. Contaminated implants must not be processed and reused. They must be disposed of after completion of the surgical procedure.

Should a supplied, absorbable implant has already been removed from its packaging but not used immediately (e.g. because the surgery had to be delayed), the implant must be discarded.

Reprocessing is not possible, as the product properties would be impaired by this process.

5 Cleaning, Disinfection, and Sterilization

WARNING

Possible danger to life of the patient from reuse and processing of single-use products!

IPS Implants® Resorb Mg are intended for single use in one patient. They are marked with the corresponding symbol on the label and must not be processed or reused.

Even if the single-use product seems intact after use, its integrity may be compromised, and minor defects and/or internal weak points may be present, possibly causing the product to fail.

Processing and reuse can increase the risk of contamination, e. g. by carryover of germs.

WARNING

Possible danger to life of the patient from processing of contaminated single-use products!

Unused single-use products that have come into contact with bodily fluids, blood, tissue and/or the like should be considered used and must be discarded.

Processing and reuse can increase the risk of contamination, e. g. by carryover of germs.

IPS Implants® Resorb Mg are supplied in sterile condition (R sterilization). The implants with the marking "Sterile R" can be used without additional pretreatment directly after removal from their sterile packaging.

Reprocessing is not possible, as the product properties would be impaired by this process.

Resterilization of the implants is not allowed.

6 Storage and Transport

WARNING

Possible danger to the life of others by shipping contaminated products!

In case of return shipments, send only clean and disinfected products in sterile packaging.

- Store and transport cleanly.
- Avoid severe temperature fluctuations.
- Protect the products from mechanical damage.

7 Ecological Information / Disposal

7.1 Packaging

KLS Martin will take back the full packaging upon request. Whenever possible, parts of the packaging will be reused.

When disposing of or recycling the packaging, please take applicable national regulations into account.

7.2 Disposal

Products that have not been fully resorbed must be disposed of properly after explantation.

In designing the product, we tried to avoid using composite materials wherever possible. This allows a high degree of recycling. We therefore offer you to take the product back for proper disposal and recycling.

7.3 National Regulations

The national regulations and disposal provisions must be observed for all disposal measures.

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