Deutsch
English
Español
Français
Italiano

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Foot switches for electrosurgical units

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The following symbols are either part of these Instructions for Use and/or of the product label.

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Safety alert symbol

CAUTION	Indicates a situation which, if not avoided, could result in minor or moderate injury.
WARNING	Indicates a situation which, if not avoided, could result in death or serious injury.
DANGER	Indicates a situation which, if not avoided, will result in death or serious injury.

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Table of Contents

1	General Information.....	22
1.1	Manufacturer.....	22
1.2	Hotline.....	22
1.3	Event Reporting Duty.....	22
1.4	Notes on this Document.....	23
1.5	Terms & Acronyms.....	23
2	Scope of Delivery.....	24
2.1	Inspecting the Delivery for Completeness and Correctness.....	25
3	Intended Use.....	26
3.1	Intended Purpose.....	26
3.2	Indications.....	26
3.3	Possible Side Effects.....	26
3.4	Residual Risks.....	26
3.5	Target Patient Groups.....	26
3.6	Users.....	26
3.7	Ambient Conditions During Application.....	26
4	Operation.....	27
4.1	Before Use.....	27
4.2	Installing and Connecting the Single-Pedal Foot Switch.....	28
4.3	Installing and Connecting the Double-Pedal Foot Switch.....	28
4.4	Foot Switch Functions.....	30
5	Cleaning, Disinfection and Sterilization.....	31
5.1	Preparation before Cleaning.....	31
5.1.1	Wipe Disinfection.....	32
5.1.2	Machine Cleaning and Disinfection.....	32
5.2	Storage and Transport.....	33
6	Repair.....	33
7	Ecological Information / Disposal.....	34
7.1	Packaging.....	34
7.2	Disposal.....	34
7.3	National Regulations.....	34

1 General Information

1.1 Manufacturer

We thank you for buying one of our products.

This product bears the CE mark. This means that it meets the basic requirements on the safety and the performance of medical devices set forth by the respective applicable European regulations.

We are the manufacturer of this product:

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info@klsmartin.com · www.klsmartin.com

1.2 Hotline

Should you have any questions on how to handle the device or product or use it for clinical applications, please do not hesitate to contact the Product Management:

Phone: +49 7461 706-0

E-mail: info@klsmartin.com

In case of technical questions and questions regarding maintenance contracts and training, please contact our Martin Service Center:

Phone: +49 7461 706-343

E-mail: service@klsmartin.com

NOTICE

All packaging, and in some cases also the product itself, is marked with a batch (LOT) and a catalogue number (REF). When making a product complaint please always include the LOT and REF.

1.3 Event Reporting Duty

All serious adverse events occurring in connection with the product must be reported to KLS Martin and the respective authorities without delay.

1.4 Notes on this Document

Possible danger to the life of patients, users and others when these Instructions for Use are disregarded!

Read Instructions for Use completely and follow them carefully. In particular, be sure to note all cautionary and warning notices.

This document refers to persons of all genders alike. Reference to various genders is avoided purely for reasons of improved readability.

The electronic version of these Instructions for Use can be requested at <https://www.klsmartin.com/en/services/instructions-for-use/>.


1.5 Terms & Acronyms

Acronym	Description
AP	Anesthetics check: Product is explosion-proof
HF current	High-frequency current
IP X8	International Protection (protection class definition): Protection against permanent submersion
MRI	Magnetic Resonance Imaging
NMR	Nuclear magnetic resonance spectroscopy

2 Scope of Delivery

This document applies for the following products:

Illustration (example)	Designation
	Single-pedal foot switch
 	Double-pedal foot switch

REF	
REF 80-811-50-04	REF 80-830-00-04
REF 80-811-30-04	REF 80-830-01-04
REF 80-811-32-04	REF 80-830-02-04
REF 80-812-10-04	REF 80-830-10-04
REF 80-830-20-04	REF 80-830-21-04
REF 80-824-01-04	

2.1 Inspecting the Delivery for Completeness and Correctness

Immediately upon receipt, the goods must be checked for completeness and damages.

Notice of any such damage must be given immediately.

After delivery of the product, check the original packaging and the packaging seal for damages.

The foot switches are packed in a fully hygienic but non-sterile packaging.

3 Intended Use

3.1 Intended Purpose

The reusable, non-invasive foot switch is intended to be used temporarily with an electrosurgical unit by KLS Martin and in particular enables its intended purposes and targets and directly supports its medical function with regards to these purposes.

3.2 Indications

The indications result from the intended purpose.

3.3 Possible Side Effects

There are no known product-specific adverse effects when used as intended.

The user is responsible to obtain informed consent from the patient with respect to potential adverse effects and complications.

3.4 Residual Risks

The user is responsible for obtaining informed consent from the patient with respect to potential adverse effects and complications in connection with the application.

3.5 Target Patient Groups

There are no restrictions regarding the patient target group.

3.6 Users

The products must only be used by qualified surgical staff.

The cleaning, disinfecting and sterilizing is to be performed by trained and qualified staff of the medical device preparation unit.

3.7 Ambient Conditions During Application

Use only in ORs under conditions appropriate for surgery or in medical facilities designed for this.

4 Operation

CAUTION

Danger of injuries by unintentional actuation of the foot switches!

All the following foot switches have magnetic reed contacts. Magnetic fields in the immediate vicinity can cause accidental actuation of the foot switches and thus uncontrolled emission of HF energy.

- Keep the foot switches away from magnetic fields (especially MRI and NMR spectroscopy)!
- Keep magnetic materials away from the foot switches!

Foot switch with reed contact

REF 80-811-30-04	REF 80-811-32-04
REF 80-811-50-04	REF 80-824-01-04
REF 80-812-10-04	REF 80-830-00-04
REF 80-830-01-04	REF 80-830-10-04
REF 80-830-02-04	REF 80-830-21-04
REF 80-830-20-04	

CAUTION

A failure of the electrosurgical unit may cause undesirable increases in output power!


4.1 Before Use

Clean and disinfect product before initial use and each reuse, see section 5 "Cleaning, Disinfection and Sterilization", page 31.

- Especially check the insulation of the cable for damage before each use.
- Damaged products (on cable or foot switch) must not be used anymore.
- Certain repairs can be performed, see section 6 "Repair", page 33.


4.2 Installing and Connecting the Single-Pedal Foot Switch

See Instructions for Use of the respective electrosurgical unit!


Connection socket (1) for the single-pedal foot switch and the double-pedal foot switch REF 80-830-10-04

- Place foot switch on the floor in a stable position.
- Lay out the connection cable of the foot switch in such a manner that it is not tensioned and there is no danger of tripping.


- Insert connector (3) of the foot switch into the corresponding connection socket (1) of the electrosurgical unit.
- Ensure that the contacts (2) are properly aligned!

- Secure plug connection by slightly turning the bayonet lock (4).

4.3 Installing and Connecting the Double-Pedal Foot Switch


See Instructions for Use of the respective electrosurgical unit!



Connection socket (1) for double-pedal foot switch

- Place foot switch on the floor in a stable position.
- Lay out the connection cable of the foot switch in such a manner that it is not tensioned and there is no danger of tripping.

- Insert connector (3) of the foot switch into the corresponding connection socket (1) of the electrosurgical unit.
- Ensure that the contact pins (2) are properly aligned!

- Secure plug connection by turning the cap nut (4) several times.

4.4 Foot Switch Functions

NOTICE

- Perform a function check before each use!
- Foot switches without AP protection must not be used in the presence of flammable or explosive substances! If the product has a green AP label, it belongs to class AP.

REF	Pedal yellow	Pedal blue	Pedal black	SWAP® Mode, round, black *	Bracket, retractable	AP protection/ IP X8
Function	Activate cutting	Activate coagulation	Activate cutting or coagulation, depending on function assigned on device	Switching between 2 preset programs	Positioning of the foot switch	
80-811-50-04			X			X/X
80-812-10-04			X			X/X
80-824-01-04	X	X				X/X
80-811-30-04	X	X				X/X
80-811-32-04	X	X		X		X/X
80-830-00-04	X	X				X/X
80-830-01-04	X	X		X	X	X/X
80-830-02-04			X			X/X
80-830-20-04	X	X				X/X
80-830-21-04	X	X		X	X	X/X

REF	Pedal light blue, left	Pedal dark blue, right	Pedal black	SWAP® Mode, round, black *	Bracket, retractable	AP protection/ IP X8
Function	Activation of the bipolar function ¹⁾	Activate coagulation	Activate cutting or coagulation, depending on function assigned on device	Switching between 2 preset programs	Positioning of the foot switch	
80-830-10-04	X	X				X/X

* only for maXium® from hardware version HW 06 and software version V3.398

* only for maXium® from hardware version HW 06 and software version V3.398

5 Cleaning, Disinfection and Sterilization

NOTICE

The operator/user is responsible for cleaning and disinfecting the device. It is essential that national regulations, including restrictions, be observed.

NOTICE

Potential damage from impermissible reprocessing of foot switches!

We generally recommend spray or wiping disinfection using a suitable disinfection agent. Never use cleaning agents that can attack the plastic surfaces, i.e. instrument cleaning solutions, abrasive cleaning agents or cleaners containing solvents (e.g. benzine, stain remover).

- Observe the information by the manufacturer of the disinfection agent regarding its suitability and concentration.
- Do not sterilize foot switches.

Foot switches with AP label are suitable for machine cleaning with some restrictions, see section 5.1.2 “Machine Cleaning and Disinfection”, page 32.

5.1 Preparation before Cleaning

The responsibility for cleaning and disinfection lies with the operator/user of the product. Be sure to observe your local cleaning and hygiene regulations.

WARNING

Risk of infection in case of contaminated handling!

Contaminated handling of the foot switch can cause severe health risks for the patient.

Wipe disinfection must be carried out after each use!

Ensure that devices and equipment are cleaned and disinfected to protect against infection!

NOTICE

Risk of damage to the surface of the foot switch!

The use of untested cleaning agents can lead to discoloration and damage to the surface.

Never clean the foot switch with abrasive cleaners, disinfectants or solvents that could scratch or damage the surface.

The foot switch may only be cleaned with agents that are suitable for plastic-coated/painted surfaces! In the case of any doubt, test a discreet area in advance.

Be sure to observe the instructions of the manufacturer of the disinfectant regarding material compatibility, dosage and exposure time.

Disinfectant residues must be safely removed before start up.

5.1.1 Wipe Disinfection

- Use a mild soap solution or standard rinsing agents as cleaning agents.
- The use of untested cleaning agents can lead to discoloration and damage to the surface.
- Wipe the surfaces of the product with a light, damp cloth and add a mild soap solution (detergent) if needed.
- Then dry all outer surfaces with a soft, clean cloth.
- Clean surfaces contaminated with a potentially infectious material (e.g. blood, secretion, or excrement) immediately.
- Use the correct effective concentration of the disinfectant.
- Only use the wiped surfaces once the disinfectant has dried.
- Comply with the hygiene guidelines
- Follow the instructions of the disinfectant manufacturer

5.1.2 Machine Cleaning and Disinfection

Foot switches with AP label are suitable for machine cleaning (thermo disinfectant / cleaning and disinfection machine) with the following restrictions:

- The maximum temperature of 55°C for the cleaning and drying process must not be exceeded.
- Shake off any remaining liquid after the rinsing process. Also check areas that are not visible.
- Do not exceed the drying time of 1.5 h.
- Supply line and gasket can be flushed a maximum of 20 times with a cleaning agent suitable for plastics.

5.2 Storage and Transport

- If the foot switch is stored on the trolley, you can use the appropriate foot switch fixations or storage facilities.
- When winding up the cable (e.g. around or besides the foot switch casing), ensure that it does not lie on the bend protector and is not tensioned. Frequent alternating bending ($\pm 45^\circ$) results in early cable breakage.
- For a potential further transport, the original packaging with its internal protection elements should be preferred.

Environmental conditions for transport and storage	Ambient temperature	-25 °C up to +70 °C
	Relative humidity	10 – 90%
	Air pressure	500 – 1060 hPa
Environmental conditions for operation	Ambient temperature	+10 °C up to +40 °C
	Relative humidity	30 – 75%
	Air pressure	700 – 1060 hPa

6 Repair

Contact the Martin Service Center in case of repair requirements, see section 1.2 "Hotline", page 22.

7 Ecological Information / Disposal


7.1 Packaging

KLS Martin will take back the full packaging upon request. Whenever possible, parts of the packaging will be reused.

When disposing of or recycling the packaging, please take applicable national regulations into account.

7.2 Disposal

This allows a high degree of recycling. In designing the product, we tried to avoid using composite materials wherever possible. This allows a high degree of recycling. We therefore offer you to take the product back for proper disposal and recycling.

	Marking of electric and electronic equipment in accordance with Directive 2012/19/EC (WEEE Directive) and the German Electrical and Electronic Equipment Act (ElektroG) This symbol on the product or its packaging indicates that the product may not be disposed of as normal household garbage.
--	---

7.3 National Regulations

The national regulations and disposal provisions must be observed for all disposal measures.

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Date of Release: 2025-06 · **REF** 90-604-58-10 · Printed in Germany · Copyright by KLS Martin SE & Co. KG · Alle Rechte vorbehalten ·