
Deutsch
English
Español
Français
Italiano

Resorb x[®]/xG/xC, SonicPins Rx

Osteosynthese-Implantate

Osteosynthesis Implants

Implantes de osteosíntesis

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Impianti da osteosintesi

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Instructions for Use

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Mode d'emploi

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CE 0297

REF 90-815-58-31
Revision 06
Date of Release: 2023-11

KLS martin
GROUP

Deutsch.....	7
English	20
Español	33
Français	46
Italiano	59

Symbolerklärung

Symbol Explanation

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Légende

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Contents

1	Notices Concerning this Document	21
1.1	Safety Instructions.....	21
1.2	Symbols Used in this Document.....	21
2	Product Liability and Warranty	22
2.1	General Information	22
2.2	Hotline	22
2.3	User's Inspection	23
2.4	Warranty	23
3	Intended Use	23
3.1	Intended Purpose	23
3.1.1	Indications	23
3.1.2	Patient Group/Patient Selection	24
3.1.3	Ambient conditions.....	24
3.1.4	Users	24
3.2	Limitations of Use	25
3.2.1	Contraindications	25
3.3	Potential Side Effects and Complications	26
4	Application	27
4.1	Description of the Components	27
4.1.1	Implant materials.....	27
4.1.2	Design.....	27
4.1.3	Modeling	27
4.1.4	SonicPins	29
4.1.5	Twist Drills	30
4.2	Prior to Implantation	30
4.3	Intraoperative Procedure	31
4.4	Upon completion of the Surgical Intervention.....	31
5	Processing	32
5.1	Basic Notes on Preparation.....	32
6	Storage and transportation.....	32
7	Disposal	32

1 Notices Concerning this Document

1.1 Safety Instructions

Non-observance of this document can lead to serious or even lethal patient injury!

Improper handling and care as well as non-intended use can lead to premature wear and/or pose a risk to patients and users!

It is the operator's responsibility to ensure that all personnel handling the product have understood and do observe the notes and instructions in this document.

- Every user is required to read this document completely and follow it carefully.
- In particular, be sure to heed all cautions, warnings and danger notices.
- Keep this document accessible to users at all times.
- It is the user's responsibility to inform the patient on all the application, cautionary, warning and danger notices contained in this document and to ensure that he/she has understood them.

1.2 Symbols Used in this Document

Throughout this document, important information (such as general or safety-related notices) is marked with the following symbols and signal words:

WARNING

Life hazard or serious injury!

Indicates a situation which, if not avoided, **could** result in death or serious injury!

CAUTION

Risk of injury!

Indicates a situation which, if not avoided, could result in minor or moderate injury!

NOTICE

Risk of material damage!

Indicates a situation which, if not avoided, could lead to material damage (loss of time, data loss, device/machine failure, etc.)!

2 Product Liability and Warranty

2.1 General Information

We thank you for having decided to buy one of our products. This product bears the CE-marking, which means that it satisfies the essential requirements laid down in the EC Directive concerning medical devices.

We are the manufacturer of this product:

KLS Martin SE & Co. KG

A company of the KLS Martin Group

KLS Martin Platz 1 · D-78532 Tuttlingen · Germany

Tel. +49 7461 706-0 · Fax +49 7461 706-193

info@klsmartin.com · www.klsmartin.com

2.2 Hotline

- Should you have any questions on how to handle the product or use it for clinical applications, please do not hesitate to contact the Product Management:

Tel: +49 7461 706-0

E-mail: info@klsmartin.com

<i>NOTICE</i>

Each product and its packaging are usually marked with a specific lot number (LOT) and a reference number (REF). Please always indicate LOT and REF in case of a complaint.

2.3 User's Inspection

- Immediately upon receipt, the goods must be checked for completeness and potential damage in transit.
- Notice of any such damage must be given immediately.
- Upon receipt of the product, the original packaging and its seal must be checked for integrity. Packages with a broken seal cannot be taken back.
- The temperature indicator must be checked immediately upon receipt. Products whose temperature indicators have changed their color from silver to black must be reported as defective and may not be used in surgical procedures.

Implants that are packed in sterile packing can be taken out of their packing and applied immediately without further pre-treatment.

Sterile implants whose sterile packaging has been opened or damaged for any reason, as well as implants taken from expired packages, must be considered non-sterile.

2.4 Warranty

Our Standard Terms and Conditions of Sale effective at the time shall apply. Agreements diverging from these Standard Terms and Conditions do not restrict the legal rights of the buyer.

3 Intended Use

3.1 Intended Purpose

The implants are intended only for fixation and stabilization inside the body in case of craniomaxillofacial osteosynthesis in non-load-bearing areas.

The product must be used according to indication. KLS Martin assumes no responsibility for complications caused by incorrect indication, incorrect choice of implant, incorrect combination of system components and surgical technique, or limitations of the treatment method or insufficient hygiene.

3.1.1 Indications

The Resorb x[®] implant system by KLS Martin is intended for surgical procedures which require internal fixation by absorbable implants for aligning, reconstructing and stabilizing bone tissue. For example:

- Osteosynthesis in non-load-bearing craniomaxillofacial areas
- Craniofacial corrective osteotomies (e.g. syndrome patients, fronto-orbital advancement) in the non-load-bearing craniomaxillofacial area
- Pediatric neurosurgery
- Pre-prosthetic augmentation

See the product documents for additional details on the indications.

3.1.2 Patient Group/Patient Selection

No restrictions.

3.1.3 Ambient conditions

The Resorb x[®] implant system is used in the sterile operating room.

3.1.4 Users

The implant system may only be used by competent medical staff who, due to their specialized training, experience and knowledge of the relevant regulations, have the necessary qualification for using the system correctly.

Users must have adequate knowledge of the rules of medical practice, as well as the state of the art and the contents of pertinent scientific articles published by medical authors.

The user must be thoroughly familiar with the surgical technique as well as with the set of implants and the instruments and accessories used.

The user bears the responsibility for

- the proper performance of the surgical procedure
- selecting the right implant components and implanting them properly
- avoiding or minimizing the general risks posed by surgical interventions

The user must be familiar with both the Resorb x[®] material and the associated system components SonicWeld Rx[®], Xcelsior water bath, and BOS drill. Their instructions for use are assumed to be known.

3.2 Limitations of Use

3.2.1 Contraindications

⚠ WARNING

Non-observance of the contraindications poses a danger of serious or even lethal patient injury!

- The osteosynthesis materials Resorb x[®], Resorb xG and Resorb xC must not be used in load-bearing areas (e.g. fractures of the lower jaw). The plates and meshes could otherwise break due to excess strain and jeopardize the success of the treatment.
- Resorb x[®] plates and meshes cannot replace bone tissue permanently. For this reason, they must not be used for purposes of permanent bone replacement (e.g. for cranioplasties).
- Obvious infections
- Hypersensitivity to foreign bodies
- Suspected sensitivity or allergies to the implant materials used
- Autoimmune diseases
- Circulatory problems, systemic diseases and metabolic disorders
- Insufficient or inadequate hard and soft tissue
- Bone tumors located in the implant base region
- Secondary diseases such as degenerative processes that may adversely affect the healing process
- Patients unwilling or unable to follow instructions during the postoperative phase due to their mental, neurological or physical condition
- Patients in a poor condition of health in the sense of a general contraindication (e.g. irradiation patients)
- Osteoporosis or osteomalacia or other severe structural bone damage that prevents the stable fixation of implant components
- Obvious drug or alcohol abuse

3.3 Potential Side Effects and Complications

- Implanting foreign bodies can generally cause inflammatory or allergic reactions.
- Reduced bone density and/or bone necroses due to mechanical stress
- Vascular changes
- Dislocations and partial dislocations
- Early and late infections, surface and/or deep infections
- Nerve and vascular injuries
- Hematomas and slow wound healing
- Periarticular calcification
- Reduced joint function and mobility
- Reduced joint load-bearing capacity; joint pain
- Allergic response or sensitivity to the implant material
- Nerve damage due to surgical trauma
- Insufficient ossification or delayed bone healing, osteoporosis, osteolysis, osteomyelitis, slow revascularization, infections or bone necrosis, may cause the implant system to loosen, crack, fracture or become bent.
- Even after biological decomposition of the implant system, the patient must still receive adequate post-operative care in order to avoid fracturing of the new bone tissue. Pain, discomfort, anomalous sensations caused by the implant; implant palpability
- Increased reaction of the soft tissue around the osteotomy and/or the implant system
- Poor healing
- During the absorption phase a temporary increase of volume (swelling) of the implants can occur.

Apart from the side effects described above, it should be noted that any surgical intervention can lead to complications not necessarily caused by the implant system (such as infections, nerve injuries and pain).

4 Application

4.1 Description of the Components

4.1.1 Implant materials

Component	Base material
Resorb x [®] , SonicPins Rx	Poly-D,L-lactide (PDLLA, 50% D-lactide and 50% L-lactide)
Resorb xG	85% poly-L-lactide (PLLA) and 15% polyglycolide (PGA)
Resorb xC	80% Poly-D,L-lactide (PDLLA, 50% D-lactide and 50% L-lactide) 20% Calciumcarbonate

In accordance with their fundamental property of being a resorbable implant system, Resorb x[®] implants are biologically degraded in the patient's body, and the degradation products are excreted. Therefore, the implants cannot be reused.

For reasons of better legibility, the individual basic materials are mentioned only in cases where there are differences in regard to their use. In all other cases, the term "Resorb x[®]" or "SonicPins" is used uniformly.

4.1.2 Design

The osteosynthesis systems Resorb x[®]/xG/xC consist of the SonicPins Rx and various plates, meshes, foils and membranes.

4.1.3 Modeling

Osteosynthesis implants must be accurately adapted to the contour of the patient's bone surface by modeling.

For preparation of the modeling of Resorb x[®] implants and meshes, thermo and modeling devices must be used that are approved for sterile medical environments and have an operating temperature of 70–90°C (158–194°F). The heating time of the implant can differ depending on the device. However, it should usually range from 15–30 s. We recommend the Xcelsior water bath (REF 52-400-10-04) from KLS Martin. Other thermo devices have not yet been tested for their suitability and should therefore not be used.

⚠ WARNING**Risk of patient injury from implant breakage!**

Cold deformation of the implants can result in critical stresses in the material and to implant breakage.

- Contour, model and cut the plates, meshes, foils or membranes only in the warm state which enables them to be modeled.
- Also refer to the Instructions for Use of the Xcelsior water bath.

The duration of the modelability of the heated Resorb x[®] material depends on its temperature, its material strength, and the environment conditions.

- Please note that a plate heated in a water bath cools down relatively quickly due to the liquid residue on its surface.
- Create the desired contour with the fewest possible modeling processes to rule out an increased risk of breakage and damage to the implant.
- Avoid sharp angles and small bending radii to reduce the risk of breaking the implant.
- It is not allowed to heat and model the component more than three times.
- Check the implant for damage such as deformed or torn-out screw holes after adapting it to the required contour. Damage of this type reduces the load-bearing capacity of the implant and may cause it to break.

Implants must not be stretched or squeezed.

If necessary, resorbable implants can be shortened in the middle between the holes. Sharp edges produced during cutting of the implant must be chamfered and smoothed.

Wrong selection or use of material can lead to undesired results.

⚠ WARNING**Risk of patient injury from undesired implant deformation!**

SonicPins and implants which are already pre-shaped must not be heated. Thermal heating would compromise their defined basic shape.

Bending Templates

As an option, the use of bending templates can be beneficial for modeling.

The bending templates allow the user to manually define the contour of the bone surface which can later be applied to the implant. This way, the use of implants that do not fit and frequent remodeling can be avoided.

WARNING

Danger of injury and risk of system failure from improper use!

Bending templates must not be implanted. We recommend using bending templates only once. Reprocessing after usage is not validated by the manufacturer.

4.1.4 SonicPins

CAUTION

System failure from SonicPin deformation!

SonicPins must not be heated or bent before they are inserted correctly.

SonicPins are approved solely for usage in combination with the SonicWeld Rx[®] system (REF 52-500-00-04 and REF 52-500-20-04). See the SonicWeld Rx[®] instructions for more detailed information.

- When using the SonicWeld Rx[®] system, make sure that the core hole of the SonicPins is pre-drilled exclusively with the special drills (triple ring identification marking). Tapping is not required.
 - 3 green rings = SonicPins with 1.6 mm diameter
 - 3 red rings = SonicPins with 2.1 mm diameter

Micro SonicPins Rx[®] are intended for the fixation of Resorb x[®] membranes. In contrast to Resorb x[®] plates and meshes, these membranes are inserted only after the Micro SonicPins Rx[®] have been positioned and subsequently welded together with the latter.

4.1.5 Twist Drills

Twist drills are intended for single use only.

- When twist drills are used, ensure that proper cooling is provided to avoid thermal injury and necrosis on the bone.
- Follow the manufacturer's instructions regarding the handpiece to be used with the twist drill and the maximum speeds.

Application of too much force can lead to excess strain that can damage or break the drill.

- Take precautions to prevent any breaking of the twist drill from injuring the patient, the user, or others.
- Work with appropriate protective devices (e.g. drill sleeve, drill tube) to avoid damage to tissue in the area of the drill.
- Operate the drills at a maximum speed of 1,000 rpm. Position the drill orthogonally and start the drilling process at low pressure and moderate, later increasing, speed.
- Abort the drilling process if imbalance occurs and replace the defective drill bit.

4.2 Prior to Implantation

WARNING

Do not combine different materials or systems!

Implant systems by KLS Martin must not be combined with systems by other manufacturers.

Use of instruments that are not specifically intended to be used with Resorb x[®] and SonicWeld Rx[®] technology can lead to failure of the implant system.

- The transport/storage packing of sterile implants is non-sterile. For this reason, sterile implants must be unpacked from their transport packing outside the sterile area of an operating room and transferred to the sterile area without compromising the sterility of the implant.
- All implant components, instruments and equipment must be in a sterile, ready-to-use condition before they are used. All absorbable implants are sterilized at the factory and packed in sterile packing.
- Implants from damaged or opened sterile packing and from packing with passed expiration dates must be considered non-sterile and disposed of.
- After delivery of the product, check the original packaging and the packaging seal for integrity.
- Check the expiration date, see symbol  on the external protective packaging.

NOTICE

The implant components used must be documented in the patient file with REF, article designation, and LOT. Documentation of the LOT is the only way to ensure clear traceability of the implants used.

4.3 Intraoperative Procedure

⚠ WARNING

Risk of contamination with chips and fragments formed during the processing of osteosynthesis plates!

Impurities on the implant, instruments, the surgical periphery or the surgical team (e.g. fragments) caused during the machining process must be removed by suitable measures before any contact with the patient (e.g. by washing in sterile water, changing the surgical gloves).

- Always ensure precise contouring and patient-specific adjustment of the Resorb x[®] implants or meshes on the bones.

4.4 Upon completion of the Surgical Intervention

⚠ WARNING

Risk of injury and/or system failure in case of improper behavior of the patient!

- Instruct patient and provide postoperative care:
 - All physical activities that could affect the implant system or the bone elements connected to it must be refrained from during the entire duration of the treatment. This applies especially to the carrying of heavy loads (danger of falling) and sporting activities that may involve impacts or blows.
 - If necessary, wear protective equipment (e.g. protective collar, splint, helmet, etc.) to ensure that the implant system is not excessively strained even in situations which cannot be consciously influenced (e.g. during sleep).

5 Processing

5.1 Basic Notes on Preparation

Resorbable implants by KLS Martin are delivered in sterile condition (gamma-ray sterilization). The implants with the **STERILE R** label can be taken out of their packing and applied immediately without further pre-treatment.

Reprocessing is not possible since this would damage the product properties.

Re-sterilization of the implants is not allowed.

6 Storage and transportation

Storage at room temperature. A maximum temperature of 37°C (98.6°F) must not be exceeded.

- Stop using the implant material as soon as the temperature indicator switches from silver to black!

7 Disposal

The national disposal regulations and guidelines must be observed when disposing of packaging material and potentially infectious material (e.g. after removing of metal).

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