
Deutsch
English
Español
Français
Italiano
Português

Monopolare HF-Elektrode (wiederverwendbar)

Monopolar HF electrode (reusable)

Electrodo de AF monopolar (reutilizable)

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Instructions for Use
Monopolar HF electrode
(reusable)

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All rights reserved
We reserve the right to make alterations

English

1 Scope

Monopolar HF electrodes item no.:
Maximum rated voltage of accessories:

Item number	U _{max}
80-510-06-04; 80-510-08-04; 80-510-12-04; 80-511-06-04; 80-515-06-04; 80-515-08-04; 80-515-12-04; 80-515-14-04; 80-515-16-04; 80-516-06-04; 80-520-06-04; 80-520-08-04; 80-521-04-04; 80-521-02-04; 80-522-04-04; 80-534-06-04; 80-534-12-04; 80-540-12-04; 80-542-06-04; 80-542-08-04; 80-542-12-04; 80-550-06-04; 80-560-06-04; 80-562-06-04; 80-562-08-04; 80-674-12-04; 80-674-13-04; 80-675-13-04; 80-676-13-04; 80-677-12-04; 80-678-05-04; 80-678-10-04; 80-678-11-04	4.3 kVp
80-520-09-04	2.9 kVp
80-680-12-04; 80-681-12-04; 80-682-12-04; 80-683-12-04; 80-684-12-04; 80-685-12-04; 80-686-12-04	0.5 kVp
80-610-10-04; 80-612-10-04	1.3 kVp

Please also refer to the label or catalog.

When combined with other HF accessories, the maximum rated voltage for all accessories is the lowest rated voltage among any of the accessories (see also item 2 “Intended use”).

2 Intended use

The HF electrodes may only be used by properly trained medical professionals. The HF electrodes are intended for open surgery and endoscopic surgery, and are used to cut and coagulate biological tissue. They are not intended for direct contact with the heart or the central circulatory system. The HF electrodes are connected to the monopolar outlet of the electrosurgical generator using an appropriate monopolar HF cable or an appropriate HF handle. The monopolar HF electrodes have a 2.4 mm or 4 mm connector. The 2.4 mm electrodes can be connected to the KLS MARTIN handle, item no. 80-210-23-04, or to a comparable product. The 4 mm electrodes can be connected to the KLS MARTIN handle, item no. 80-217-13-4, or to a comparable product. The voltage of the generator used may not exceed the maximum rated accessory voltage (see item 1) of the HF electrodes or a maximum frequency of 4 MHz. The electrodes are activated using a foot switch or an HF handle. The use of a smoke evacuator is recommended. Please observe the instructions for use for the specific accessories and generator.

3 Safety precaution – WARNING!

The maximum rated voltage of the HF electrodes can be found in these instructions for use, on the label, or in the current product catalog. If there is any uncertainty please contact the manufacturer.

HF electrodes must be cleaned, disinfected, and sterilized using a validated procedure (DIN EN ISO 17665) (see item 5 “Cleaning, disinfection, and sterilization”) before each use. A visual inspection must therefore be performed prior to each use (see item 6 “Visual inspection”). Ensure that the HF electrode is securely fastened in the handle. This must be performed carefully to avoid damaging the HF electrode and/or injuring the patient or the surgical staff. The tip of the electrode may be damaged if excess force is exerted.

Before use, confirm that a dedicated neutral electrode is correctly attached to the patient and properly connected to the HF generator. The HF electrode should not be activated as long as it is in contact with metal objects and/or optical devices. Care must be taken to ensure that the device does not come into proximity with flammable substances during the entire duration of use, as these pose a risk of explosion. The tip of the electrode may remain hot even after the electrosurgical power has been turned off, and may cause burns.

The use of HF power can be harmful to pacemakers and implantable cardioverter-defibrillators. Patients with these devices must consult a cardiologist prior to surgery.

4 Safety notes on precautions, preparation, and use

a) Read and observe the specifications in the instructions for use

Please read all of the instructions for use before using the electrosurgical instrument. This also applies to the instructions for use for the accessories used, including those of the HF neutral electrode and HF generator used for monopolar applications. Always observe and follow all specifications, safety notes, and warnings in all instructions for use.

b) Surroundings

Care must be taken to ensure that the device does not come into proximity with flammable substances (anesthetic agents, oxidizing gases, endogenous gases, etc.) during the entire duration of electrosurgical use, as these pose a risk of explosion. Use non-flammable disinfectants and cleaning agents, and do not use alcohol-based tinctures or similar liquids. All compounds containing oxygen must be kept closed and airtight for the duration of surgery.

c) Patient positioning

Ensure that the patient is positioned correctly, i.e. use insulating operating table dressings that are dry, absorbable, and impermeable. Insulate the patient from all conductive surfaces and items. Dry gauze must be padded in skin folds, inframammary folds, and between the extremities. Fluid accumulations in body cavities should be removed before the start of surgery.

d) Connections

Prior to use, ensure that the handle or cable used has been properly connected to the HF generator and that the proper power settings have been selected and are displayed. The instructions in the HF generator and HF handle/HF cable instructions for use must be followed.

e) Patient reactions

All electrosurgical instruments may potentially cause muscle stimulation during use. Although this product has been designed to minimize the risk of this undesirable effect, muscle stimulation may nonetheless trigger unexpected movement by the patient within the surgical field.

f) Handling electrosurgical instruments

The tip of the instrument may not be touched during use. The tip of the instrument may remain hot even after the electrosurgical power has been turned off, and may cause burns.

If the electrosurgical instrument is not being used, it should be laid on a dry, clean, non-conductive, and highly-visible surface that is not in contact with the patient. Unintentional activation of the instrument may result in burns to the patient.

The completeness of the system must be confirmed at the end of surgery.

5 Cleaning, disinfection, and sterilization

KLS Martin SE & Co. KG recommends using the validated reprocessing procedure described below. Other equivalent procedures may also be used. The user is responsible for ensuring the suitability of the procedure actually performed using appropriate measures (e.g. validation, routine monitoring, testing of material compatibility). It is not possible to define a maximum limit to the number of feasible cleaning, disinfection, and sterilization cycles due to the design, the materials used and the intended purpose. Depending on the type and duration of use, the HF electrodes are naturally subject to wear when used as intended.

A visual inspection must therefore be performed prior to each use (see item 6 “Visual inspection”).

5.1 Preparation for cleaning

Remove the HF electrodes from their packaging. Place them into a container/equipment intended for cleaning/sterilization. It is not necessary to disassemble the HF electrodes.

5.2 Pre-cleaning

The HF electrodes must be cleaned immediately prior to each use. For pre-cleaning, use only water and, if needed, aldehyde-free disinfectants that do not fix proteins. Thoroughly remove any surface debris using a soft brush or a synthetic fabric cloth, as particles or dried secretions may otherwise adhere to the device. This may render later cleaning and sterilization difficult or even impossible. Take care to ensure that areas which are difficult to reach are cleaned well and rinsed multiple times.

Hollow spaces and the inside of the electrodes must be flushed intensively (> 30 sec.) with at least 3 x 20 ml cold tap water (< 40 °C) using a flushing adapter (e.g. from Medisafe), a syringe, or a pressurized water pistol. **This pre- and post-cleaning step must be performed before further manual cleaning or cleaning with the washer/disinfector is performed.**

5.3 Manual Cleaning and Disinfection

Prepare separate cleaning and disinfection baths: Prepare an immersion bath for cleaning that contains a suitable liquid cleaning agent. Use a disinfectant that is compatible with the cleaning agent and that is also suitable for preparing an immersion bath. Please follow the instructions provided by the cleaning agent/disinfectant manufacturer and only use cleaning agents with a pH between 5.5 and 12.3 that are suitable for medical devices made of metal and plastic. The cleaning agent gigazyme® (Schülke & Mayr), 0.5%, and the disinfectant Korsolex plus, 3% are recommended. Do not use strong alkaline cleaning agents. Strong alkaline cleaning agents with a pH between 10-13 have a detrimental effect on the life span of the HF electrodes.

- Place the HF electrodes completely into the cleaning agent (e.g. 0.5% gigazyme®) in an ultrasound bath and sonicate for 5 minutes at a frequency of 35 kHz. Follow the instructions of the cleaning agent manufacturer. Take care to ensure that the products do not touch any other parts in the ultrasound bath. Take care to ensure that acoustic shadows do not occur in the ultrasound bath.
- Then clean the HF electrodes using a soft brush under running cold tap water (< 40 °C). Hollow areas and the inside of the electrodes must be flushed intensively using a pressurized water pistol (> 30 sec.) or similar device.
- Rinse the HF electrode thoroughly under tap water (< 40 °C) for at least 1 minute to remove any cleaning agent residue.
- Inspect the HF electrodes visually for any remaining contamination. If visible contamination remains, repeat the previous steps until all contamination has been completely removed.
- Place the HF electrodes completely into the disinfectant bath with (e.g. Korsolex Plus, 3%) for 15 minutes. Follow the manufacturer’s recommended incubation time. Ensure that the disinfectant reaches all parts of the HF electrode. Hollow areas and the inside of the electrodes must be flushed with a syringe several times (at least 3 x 20 ml).
- Then rinse thoroughly for at least 1 minute with demineralized water to remove all disinfectant residue. In addition, all narrow and inaccessible regions and hollow areas as well as the inside of the electrodes must be flushed with a syringe several times (at least 3 x 20 ml) with cold demineralized water.
- Dry the HF electrodes with a lint-free cloth. Also dry the hollow areas and channels with sterile compressed air.

5.4 Machine Cleaning and Disinfection

Only washers/disinfectors whose efficiency has been tested in accordance with DIN EN ISO 15883 may be used. Please follow the instructions provided by the cleaning agent manufacturer and only use cleaning agents with a pH between 5.5 and 12.3 that are suitable for medical devices made of metal and plastic. neodisher® MediClean forte (Dr. Weigert GmbH & Co. KG) is recommended. Neutralizing agents should not be used. Use a thermal disinfection program. The program sequence must adhere to the instructions from the washer/disinfector manufacturer. During automated cleaning/disinfection, the HF electrodes must be securely positioned and protected from mechanical damage.

Do not clean together with sharp or pointed objects.

Place the HF electrodes into a suitable washing basket. Start the program sequence with the following properties:

- 1 min pre-rinse with cold water
- Emptying
- 3 min pre-rinse with cold water
- Emptying
- Wash for 5 min at 55 °C with 0.5 % alkaline cleaning agent
- Emptying
- 3 min neutralization with warm mains water (>40 °C) and neutralizer (0.1 % neodisher® Z)
- Emptying
- 2 min intermediate rinsing with warm demineralized water (> 40 °C)
- Emptying

5.5 Disinfection

Perform automated thermal disinfection with an A0 value in accordance with national requirements (see ISO 15883, A0 value > 3,000): 5 min at 92 °C +/- 2 °C.

5.6 Drying

- 30 min at 90 °C
- Remove flushing adapter

After the program is complete, remove the HF electrodes and check them for any residual contamination. If any residues remain, repeat the automated cleaning and disinfection step until contamination is no longer visible. Dry hollow spaces and inadequately dried areas with sterile

compressed air at < 2 bar. Pack the instruments as soon as possible after removal (after drying further at a clean location) into disposable paper/film sterilization packaging (double-layer packaging) or into a sterilization container (as defined in ISO 11607 and EN 868).

5.7 Sterilization

Instruments may only be used if cleaned in advance. The HF electrodes are only designed for steam sterilization in an autoclave (fractionated pre-vacuum procedure with adequate product drying). The HF electrodes must be sterilized using saturated steam at a minimum of 134 °C and maximum of 137 °C with a holding time of at least 5 minutes to no more than 20 minutes, and dried under vacuum for at least 10 minutes. The sterilizer (Class B) must comply with applicable national standards and regulations (e.g. DIN EN 13060 or DIN EN 285). The sterilization process has been validated in accordance with DIN EN ISO 17665. Please follow the recommendations of the sterilizer manufacturer pertaining to loading, handling, and drying times. The HF electrodes may not be sterilized using hot air, EtO gas, gamma irradiation, or plasma.

NOTE: The HF electrodes must be cooled to room temperature before use. The operator must ensure that sterility is maintained after the sterilization process. Should the recommended chemicals and machines described above not be available for manual or automated cleaning and disinfection, the user is required to validate their own procedures appropriately. Similarly, if a different sterilization method than that described here is chosen, this different method must also be appropriately validated by the user.

Limits on reprocessing

The product life span depends on device wear, damage, and frequency of reprocessing.

6 Visual inspection

The insulation on the HF electrodes must be checked for pinched sites or damage before each use. HF electrodes with damage or dents may not be used.

7 Repair and modification

Defective HF electrodes may not be repaired. They must instead be replaced by new HF electrodes. Unauthorized modifications or repairs are strictly prohibited and void the manufacturer's warranty.

8 Packaging, storage, transportation, handling

The HF electrodes must be stored in a clean and dry environment. They should be stored individually in a protective container with individual compartments or sealed in film. The utmost care must always be used when transporting, cleaning, disinfecting, maintaining, and sterilizing the HF electrodes, especially their fine points and other sensitive areas. The operator must ensure that sterility is maintained after the sterilization process.

9 Returns

Returns will only be accepted if labeled as “no hygiene hazard” or “not contaminated” and packaged securely for shipping.

10 Disposal

The HF electrodes, packaging material, and accessories must be discarded in accordance with national laws and regulations.

11 About these instructions for use

These instructions for use must be kept accessible to all users for the duration of usage of the HF electrodes.

A version of the instructions for use can be found on www.klsmartin.com

©2020. KLS Martin SE & Co. KG, D-78532 Tuttlingen/Germany. All rights reserved. All company names, trademarks, service marks, trade names, product names, and logos in this document are the property of their respective holders and must be handled accordingly.

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