

A man in a white shirt and shorts is jumping in the air on a sandy beach, reaching for a red frisbee. A woman in a light blue shirt and dark shorts stands on the sand to the right, watching him. The background shows waves crashing on the shore under a clear sky.

IPS CaseDesigner<sup>®</sup> and DTX Studio Go<sup>™</sup> login

## IPS CaseDesigner® login

- Login system integration in IPS CaseDesigner® from version 2.3 on



**KLS martin**  
GROUP

**IPS®**

Select your username.

 christoph.krueger

 |

Log in automatically **Log in**

[Forgot your password?](#)

[Close application](#)

## IPS CaseDesigner® login

- Initial Password creation via email
  - Search in your emails for an email from [no-reply@dtxstudio.com](mailto:no-reply@dtxstudio.com)  
→ I can't find the email
  - Your **username** is displayed
  - Click on „**Complete registration**“ to set the initial password

Scroll down to see this email in other languages. Available languages are:  
English, Chinese, French, Italian, Japanese, Portuguese, Spanish

## DTX Studio account created

Congratulations! You now have access to the DTX Studio platform.  
You can log in with the following credentials:

Username: **christoph.krueger**

Click on the button below to set up your DTX Studio account.

Complete registration

If the link does not work, you may copy the following URL in your browser.

<https://uat-nbp-auth-service.dtxstudio.com/Account/SetPassword?setUserPasswordBewit=N2FhMzMxYTktYzYxM500MzM0LTk4NGENWU2NTG2YmVjOGQ0XDE2NzI1ODI3MTRcNmhtYmITMS10ck9qdmlzc2dWUVZwMEFjM3Fjc0RjaGV0OGhLT0FPaFpWRT1c&isResetLink=False>

Need help? Contact the software support team directly at [support@dtxstudio.com](mailto:support@dtxstudio.com)  
DTX Studio™ is a digital platform for dental treatments connecting your workflow from beginning to end. Dental professionals worldwide use it to securely share patient information. Learn more at [www.dtxstudio.com](http://www.dtxstudio.com)

Best regards,  
The DTX Studio team

## IPS CaseDesigner® login

- Initial password creation in IPS CaseDesigner®
  - Install version 2.3 on your device
  - Choose your **username** in the drop down box via the arrow. Depending on your account, several user names are linked.
  - Click on „**Forgot your password?**“
  - Search in your emails for an email from [no-reply@dtxstudio.com](mailto:no-reply@dtxstudio.com)
  - → Set your initial password
  - Go back to IPS CaseDesigner® and log in

KLS martin  
GROUP

IPS®

Select your username.

christoph.krueger

Log in automatically

Log in

Forgot your password?

Close application

## IPS CaseDesigner® login

- Login skipping in IPS CaseDesigner®
  - Install version 2.3 on your device
  - Tick [log in automatically](#) if you would like to skip the login mask for the next time

KLS martin  
GROUP

IPS®

Select your username.

christoph.krueger

Log in automatically


Log in

Forgot your password?

Close application


## DTX Studio™ Go

- Web-based account management platform DTX Studio™ Go accompanying IPS CaseDesigner®
  - Check your license key
  - Check your maintenance expiry date
  - User management
  - Check the installation overview and clear unused installation slots
  - Download installers for Windows and Mac




## DTX Studio™ Go

- Web-based account management platform DTX Studio™ Go accompanying IPS CaseDesigner®
- View the product documentation
  - IFU
  - Help files
  - Installation guides
  - Computer guidelines (hardware requirements)



## DTX Studio™ Go

- Platform access
  - [go.dtxstudio.com](https://go.dtxstudio.com)
- Log in access
  - Click on **Log in** (upper right corner)





## DTX Studio™ Go

- Log in
  - Type in your **username** (like in IPS CaseDesigner®)
  - Type in your **password** (like in IPS CaseDesigner®)
  - **Stay-logged-in** function
  - „**Forgot your password?**“ fuction (same function like in IPS CaseDesigner®)

DTX Studio  
Go

Anwendername  
Christoph.Krueger

Kennwort  
●●●●●●●●

Angemeldet bleiben

[Passwort vergessen?](#) [Login](#)

Du hast noch kein Konto? [Anmelden](#)