




LevelOne Fixation

Osteosynthesis 2.0 Mini



Oral and maxillofacial surgery is our passion! We also want to continue our development along with our customers. Day in, day out, we work to develop innovative products and services that satisfy the highest quality demands and contribute to the patient's well-being.

Table of contents

	Pages
Product features	4 - 5
Fields of use	6
Overview of materials	7
Product range	8 - 67
■ Implants	
■ Screws	8 - 13
■ Cranium	14 - 15
■ Midface	16 - 17
■ Mandible	18 - 21
■ Orthognathic surgery	22 - 27
■ Universal osteosynthesis plates	28 - 41
■ Universal mesh plates	42 - 43
■ Rigid Blue plates	44 - 45
■ Midface - alternative plates	46 - 47
■ Instruments	48 - 55
■ Storage	56 - 65
■ Sterile kits	66 - 67
Index	68 - 70

Feature, function and benefit







With LevelOne, KLS Martin offers a comprehensive system for plate and screw osteosynthesis. A clearly structured variety of implants, universally applicable and designed for specific applications, covers the most diverse areas of oral and maxillofacial surgery.

The storage system consists of basic modules, storage modules, and color-coded screw storage packs that can be filled individually and combined according to requirements.

This offers surgeons the option of choosing between preconfigured modules or configuring their own system. There are no limits to the configuration options.

LevelOne – a coherent approach that has become indispensable in osteosynthesis.

LevelOne – implants

	Feature and function	Benefits
maxDrive® Screws 	■ Predefined, self-centering guidance ■ Self-retaining mechanism ■ Various designs ■ Standard Ø 2.0 mm ■ Drill-free Ø 2.0 mm ■ Locking screws Ø 2.0 mm ■ Emergency Ø 2.3 mm	■ Easy pick-up and insertion of the screw ■ Can be screwed in at an angle ■ Direct force transfer from the screwdrive blade to the screw ■ Screw variance for multiple application options
Clip magazines and instruments 	■ Color-coded screw clips and instruments	■ Quick and easy identification of corresponding implants and instruments
Plates 	■ Variety of plate designs ■ Universal plate geometries ■ Specific implant designs	■ Versatile coverage of a very wide range of applications
Storage 	■ Freely combinable LevelOne storage modules	■ Compact storage for instruments and implants ■ Flexibility and adaptation to individual customer requirements

LevelOne Fixation



Fields of use

The possible areas of application with the 2.0-mm system are as diverse as the LevelOne system itself:

- Mandibular fractures
- Fixation of bone grafts
- Mandibular reconstruction with bone transplants
(however, not for bridging continuity defects)
- Fractures in load-bearing regions of the cranium



Overview of materials

		
	Plate selection	Screws
Materials	 Pure titanium in implant quality	 Titanium alloy in implant quality Ti-6Al-4V
Standards	DIN ISO 5832-2 DIN 17850 ASTM F 67	DIN ISO 5832-3 DIN 17851 ASTM F 136a

Mini screws 2.0 mm

maxDrive®  

Mini screws Ø 2.0 mm self-retaining

	Ø x L (mm)	maxDrive® 5	STERILE R 5	STERILE R 1
	2.0 x 4 mm	25-872-04-09	25-872-04-75	25-872-04-71
	2.0 x 5 mm	25-872-05-09	25-872-05-75	25-872-05-71
	2.0 x 6 mm	25-872-06-09	25-872-06-75	25-872-06-71
	2.0 x 7 mm	25-872-07-09	25-872-07-75	25-872-07-71
	2.0 x 9 mm	25-872-09-09	25-872-09-75	25-872-09-71
	2.0 x 11 mm	25-872-11-09	25-872-11-75	25-872-11-71
	2.0 x 13 mm	25-872-13-09	25-872-13-75	25-872-13-71
	2.0 x 15 mm	25-872-15-09	25-872-15-75	25-872-15-71
	2.0 x 17 mm	25-872-17-09	25-872-17-75	25-872-17-71
	2.0 x 19 mm	25-872-19-09		25-872-19-71

Angular stable screws Ø 2.0 mm self-retaining

	Ø x L (mm)	maxDrive® 5	STERILE R 5	STERILE R 1
	2.0 x 5 mm	25-882-05-09	25-882-05-75	25-882-05-71
	2.0 x 7 mm	25-882-07-09	25-882-07-75	25-882-07-71
	2.0 x 9 mm	25-882-09-09	25-882-09-75	25-882-09-71
	2.0 x 11 mm	25-882-11-09	25-882-11-75	25-882-11-71
	2.0 x 13 mm	25-882-13-09	25-882-13-75	25-882-13-71
	2.0 x 15 mm	25-882-15-09	25-882-15-75	25-882-15-71
	2.0 x 17 mm	25-882-17-09	25-882-17-75	25-882-17-71
	2.0 x 19 mm	25-882-19-09	25-882-19-75	25-882-19-71

Emergency screws Ø 2.3 mm self-retaining

	Ø x L (mm)	maxDrive® 5	STERILE R 5	STERILE R 1
	2.3 x 4 mm	25-873-44-09	25-873-44-75	25-873-44-71
	2.3 x 5 mm	25-873-45-09	25-873-45-75	25-873-45-71
	2.3 x 7 mm	25-873-47-09	25-873-47-75	25-873-47-71
	2.3 x 9 mm	25-873-49-09		25-873-49-71

Drill-free screws Ø 2.0 mm self-retaining

	Ø x L (mm)	maxDrive® 5	STERILE R 5	STERILE R 1
	2.0 x 5 mm	25-879-05-09	25-879-05-75	25-879-05-71
	2.0 x 6 mm	25-879-06-09		25-879-06-71
	2.0 x 7 mm	25-879-07-09	25-879-07-75	25-879-07-71
	2.0 x 9 mm	25-879-09-09	25-879-09-75	25-879-09-71



Explanation of icons

- Titanium alloy
- maxDrive®
- Centre-Drive®
- Units per package
- Sterile packaged implants

Centre-Drive®

Mini screws Ø 2.0 mm self-retaining



Ø x L (mm)	Centre-Drive®	
2.0 x 4 mm	25-662-04-09	25-662-04-75
2.0 x 5 mm	25-662-05-09	25-662-05-75
2.0 x 7 mm	25-662-07-09	25-662-07-75
2.0 x 9 mm	25-662-09-09	25-662-09-75
2.0 x 11 mm	25-662-11-09	25-662-11-75
2.0 x 13 mm	25-662-13-09	25-662-13-75
2.0 x 15 mm	25-662-15-09	25-662-15-75
2.0 x 17 mm	25-662-17-09	25-662-17-75
2.0 x 19 mm	25-662-19-09	25-662-19-75

Emergency screws Ø 2.3 mm self-retaining



Ø x L (mm)	Centre-Drive®	
2.3 x 5 mm	25-663-45-09	
2.3 x 7 mm	25-663-47-09	25-663-47-75
2.3 x 9 mm	25-663-49-09	25-663-49-75

Drill-free screws Ø 2.0 mm self-retaining



Ø x L (mm)	Centre-Drive®	
2.0 x 5 mm	25-669-05-09	25-669-05-75
2.0 x 7 mm	25-669-07-09	25-669-07-75

MMF screws

maxDrive® MMF screw

	Ø	Thread length	Total length	Item No.		
	2.0 mm	8 mm	12 mm	25-092-08-09	Ti	5
	2.0 mm	8 mm	12 mm	25-092-08-72	Ti	2
	2.0 mm	8 mm	12 mm	25-092-38-05	St	5
	2.0 mm	8 mm	12 mm	25-092-38-72	St	2
	2.0 mm	12 mm	16 mm	25-092-12-09	Ti	5
	2.0 mm	12 mm	16 mm	25-092-12-72	Ti	2
	2.0 mm	12 mm	16 mm	25-092-42-05	St	5
	2.0 mm	12 mm	16 mm	25-092-42-72	St	2

Amsterdam MMF screw

	Ø	Thread length	Total length	Item No.		
	2.0 mm	8 mm	12 mm	25-099-12-09	Ti	2
	2.0 mm	8 mm	12 mm	25-099-12-72	Ti	2
	2.0 mm	8 mm	12 mm	25-099-42-05	St	2
	2.0 mm	8 mm	12 mm	25-099-42-72	St	2
	2.0 mm	10 mm	14 mm	25-099-14-09	Ti	2
	2.0 mm	10 mm	14 mm	25-099-14-72	Ti	2
	2.0 mm	10 mm	14 mm	25-099-44-05	St	2
	2.0 mm	10 mm	14 mm	25-099-44-72	St	2

Amsterdam MMF
Developed in cooperation with:
Dr. Bart van den Bergh, VU Amsterdam, Netherlands



Explanation of icons

-  Titanium
-  Stainless steel
-  Units per package
-  maxDrive®
-  Centre-Drive®

STERILE IR Sterile packaged implants

FAMI MMF screw

(Fixation and adaptation in mandibular injuries)

 	Ø	Thread length	Total length	Item No.	
	2.0 mm	8 mm	10.4 mm	25-097-08-09	 
	2.0 mm	8 mm	10.4 mm	25-097-08-75	 



FAMI MMF –

developed in cooperation with:

- Lars Erikson D.D.S, Odont D. (Ph.D.), University of Lund, Sweden
- Dr. Dr. Rainer Fangmann, MSC, Wilhelmshaven, Germany
- PD Dr. Dr. R. Mischkowski, Städt. Klinikum Ludwigshafen, Germany


with clip magazine

Fast-Fix MMF System: washers, instruments and storage container



25-095-90-09

MMF washer, model Rome, short



25-095-91-09

MMF washer, model Rome, long



Note:
Used in conjunction with 2.0-mm maxDrive®, Centre-Drive® screws



Explanation of icons

St Stainless steel

1 Units per package

Fast-Fix MMF System

55-961-29-04		2.0-mm MMF module, complete, consisting of:
55-962-48-04	1	MMF module with red module rails
55-963-47-04	1	Lid for 2.0-mm MMF module

		Set recommendation
25-407-04-04	St 1	Screwdriver handle, only
25-486-97-07	St 1	maxDrive® blade, 2.0 / 2.3 mm
22-500-11-07	St 1	TC wire twisting forceps
11-865-12-07	St 1	Wire cutting scissors
39-311-20-07	St 1	Tamponade stopper



For further information, please refer to brochure "Maxillomandibular Fixation MMF System"

Implants: cranium

Blows and injuries to the cranium can cause blood congestion in the brain and a resulting increase in intracranial pressure. In such cases, the skullcap is opened using special trephine drills to drain the fluids.

For the subsequent covering of such drill holes we offer our "burr hole cover" plates.




Explanation of icons

Ti Pure titanium

 Plate profile

1 Units per package

STERILE IR Sterile packaged implants




50-311-17-09 **Ti** **1**

50-311-17-71 **Ti** **1**

Ø 17 mm

flat

 = 0.6 mm



50-311-27-09 **Ti** **1**

50-311-27-71 **Ti** **1**


Ø 17 mm

pre-shaped

 = 0.3 mm

Implants: midface

Nasal miniplates



50-399-10-09 Ti 1
50-399-10-71 Ti 1

 = 0.6 mm



50-399-12-09 Ti 1
50-399-12-71 Ti 1

 = 0.6 mm



Explanation of icons

-  Pure titanium
-  Plate profile
-  Units per package
- STERILE IR** Sterile packaged implants

Orbital plates



50-400-06-09  
 50-400-06-91  
 50-400-06-71  
 = 0.6 mm



50-400-08-09  
 50-400-08-91  
 50-400-08-71  
 = 0.6 mm



50-405-12-09  
 50-405-12-71  
 = 0.7 mm

Champy orbital plates



25-395-08-09  
 25-395-08-91  
 = 0.6 mm



25-395-10-09  
 25-395-10-91  
 = 0.6 mm



25-395-12-09  
 25-395-12-91  
 = 0.6 mm

Implants: mandible

Trapez-3D condylar fracture Plates




25-285-10-91  

14 x 14.5 mm

 = 1.0 mm

Standard



25-283-15-91  

25-283-15-71  

14 x 14.5 mm

 = 1.0 mm

Locking



Explanation of icons

-  Pure titanium
-  Plate profile
-  Units per package

STERILE IR Sterile packaged implants

Rhombus-3D condylar fracture plates



$\frac{1}{1}$

25-285-05-09  

25-285-05-71  

20 x 13 mm

 = 1.0 mm

Standard



$\frac{1}{1}$

25-283-05-09  

25-283-05-71  

20 x 13 mm

 = 1.0 mm

Locking

Magdeburg condylar fracture plates



25-380-04-09 Ti 5

25-380-04-91 Ti 1

25-380-04-71 Ti 1

 = 1.0 mm




25-382-04-09 Ti 5

25-382-04-91 Ti 1

25-382-04-71 Ti 1

 = 1.0 mm

Angular stable condylar fracture plates



25-283-25-91 Ti 1

25-283-25-71 Ti 1

26.5 x 16.1 mm, right

 = 1.0 mm



25-283-26-91 Ti 1

25-283-26-71 Ti 1

26.5 x 16.1 mm, left

 = 1.0 mm

Angular stable condylar fracture plates



25-283-20-91 Ti 1

25-283-20-71 Ti 1

35 x 18.2 mm, right

 = 1.0 mm



25-283-21-91 Ti 1

25-283-21-71 Ti 1

35 x 18.2 mm, left

 = 1.0 mm



25-288-08-09 Ti 1

21.4 x 32.5 mm, right

 = 1.0 mm




25-289-08-09 Ti 1

21.4 x 32.5 mm, left

 = 1.0 mm

Implants: orthognathic surgery

BSSO sliding plates



25-401-39-09 Ti 1

Slider

 = 1.0 mm

for fixing with standard screws
(not angular stable, 2.0 / 2.3 mm)

1/1



25-401-36-09 Ti 1

 = 1.0 mm

with 1-mm graduation



25-401-37-09 Ti 1

 = 1.0 mm

with 1-mm graduation



25-401-38-09 Ti 1

 = 1.0 mm

with 1-mm graduation




Explanation of icons


-  Pure titanium
-  Plate profile
-  Units per package
-  Sterile packaged implants

BSSO plates, twin row

Rhombus-3D condylar fracture plates



25-394-34-09  
 = 1.0 mm
 with 1-mm graduation



25-395-34-09  
 = 1.0 mm
 with 1-mm graduation



25-401-43-09  
 right
 = 1.2 mm



25-401-44-09  
 left
 = 1.2 mm

for fixing with angular stable screws (ThreadLock TS, 2.0 /2.3 mm) or standard screws (not angular stable, 2.0/2.3 mm)

Champy-Arnhem



25-392-04-09 Ti 5
 25-392-04-91 Ti 1
 25-392-04-71 Ti 1

= 1.0 mm
 with 1-mm graduation



25-394-04-09 Ti 5
 25-394-04-91 Ti 1
 25-394-04-71 Ti 1

= 1.0 mm
 with 1-mm graduation

Mandibular angle plate with sliding holes





50-405-06-91 Ti 1
 50-405-06-71 Ti 1

= 1.5 mm

BSSO sliding plate "Cologne"

Plate including slider

1/1



25-401-30-09 Ti 1
 25-401-30-71 Ti 1

= 0.8 mm



1/2

30 mm

11 mm

25-401-35-09 Ti 1
 25-401-35-71 Ti 1

= 0.8 mm



1/2

35 mm

16 mm

25-401-40-09 Ti 1
 25-401-40-71 Ti 1

= 0.8 mm




1/2

40 mm

21 mm

Chin plates



50-300-02-09 Ti 5

50-300-02-91 Ti 1

50-300-02-71 Ti 1

2 mm

= 0.6 mm



50-300-04-09 Ti 5

50-300-04-91 Ti 1

50-300-04-71 Ti 1

4 mm

= 0.6 mm



50-300-06-09 Ti 5

50-300-06-91 Ti 1

50-300-06-71 Ti 1

6 mm

= 0.6 mm



50-300-08-09 Ti 5

50-300-08-91 Ti 1

50-300-08-71 Ti 1

8 mm

= 0.6 mm



50-300-10-09 Ti 5

50-300-10-91 Ti 1

50-300-10-71 Ti 1

10 mm

= 0.6 mm



50-300-12-91 Ti 1

50-300-12-71 Ti 1

12 mm

= 0.6 mm




50-300-14-91 Ti 1

50-300-14-71 Ti 1

14 mm

= 0.6 mm

Mini chin plates



50-303-02-09 Ti 1

2 mm

= 0.6 mm



50-303-03-09 Ti 1

3 mm

= 0.6 mm



50-303-04-09 Ti 1

4 mm

= 0.6 mm



50-303-05-09 Ti 1

5 mm

= 0.6 mm



50-303-06-09 Ti 1

6 mm

= 0.6 mm



50-303-07-09 Ti 1

7 mm

= 0.6 mm



50-303-08-09 Ti 1

8 mm

= 0.6 mm



50-303-09-09 Ti 1

9 mm

= 0.6 mm

Maxilla
with 1-mm graduation

For further information please refer to the folder
"Arnett FAB Orthognathic System"

$\frac{1}{1}$



25-399-26-09 Ti 1

25-399-26-71 Ti 1

$\varnothing = 0.8 \text{ mm}$

25-399-26-89 Al 1



25-399-32-09 Ti 1

25-399-32-71 Ti 1

$\varnothing = 0.8 \text{ mm}$

25-399-32-89 Al 1



25-399-37-09 Ti 1

25-399-37-71 Ti 1

$\varnothing = 1.0 \text{ mm}$

25-399-37-89 Al 1

$\frac{1}{1}$



25-397-27-09 Ti 1

25-397-27-71 Ti 1

$\varnothing = 0.8 \text{ mm}$

25-397-27-89 Al 1



25-397-31-09 Ti 1

25-397-31-71 Ti 1

$\varnothing = 0.8 \text{ mm}$

25-397-31-89 Al 1



25-397-35-09 Ti 1

25-397-35-71 Ti 1

$\varnothing = 0.8 \text{ mm}$

25-397-35-89 Al 1



25-397-39-09 Ti 1

25-397-39-71 Ti 1

$\varnothing = 1.0 \text{ mm}$

25-397-39-89 Al 1

$\frac{1}{1}$



25-398-27-09 Ti 1

25-398-27-71 Ti 1

$\varnothing = 0.8 \text{ mm}$

25-398-27-89 Al 1



25-398-31-09 Ti 1

25-398-31-71 Ti 1

$\varnothing = 0.8 \text{ mm}$

25-398-31-89 Al 1



25-398-35-09 Ti 1

25-398-35-71 Ti 1

$\varnothing = 0.8 \text{ mm}$

25-398-35-89 Al 1



25-398-39-09 Ti 1

25-398-39-71 Ti 1

$\varnothing = 1.0 \text{ mm}$

25-398-39-89 Al 1

Mandible with 1-mm graduation

1/1



25-396-29-09 Ti 1
25-396-29-71 Ti 1
= 1.0 mm

25-396-29-89 Al 1



25-396-34-09 Ti 1
25-396-34-71 Ti 1
= 1.0 mm

25-396-34-89 Al 1



25-396-39-09 Ti 1
25-396-39-71 Ti 1
= 1.0 mm

25-396-39-89 Al 1



25-396-44-09 Ti 1
25-396-44-71 Ti 1
= 1.0 mm

25-396-44-89 Al 1



25-396-56-09 Ti 1
25-396-56-71 Ti 1
= 1.0 mm

25-396-56-89 Al 1

Chin plates

1/1




50-304-02-09 Ti 1
50-304-02-71 Ti 1
2 mm
= 0.8 mm



50-304-03-09 Ti 1
50-304-03-71 Ti 1
3 mm
= 0.8 mm



50-304-04-09 Ti 1
50-304-04-71 Ti 1
4 mm
= 0.8 mm



50-304-06-09 Ti 1
50-304-06-71 Ti 1
6 mm
= 0.8 mm



50-304-08-09 Ti 1
50-304-08-71 Ti 1
8 mm
= 0.8 mm



50-304-10-09 Ti 1
50-304-10-71 Ti 1
10 mm
= 0.8 mm

Implants: universal osteosynthesis plates

				
Standard	25-551-02-09  5 25-551-02-91  1 25-551-02-71  1  = 1.0 mm	25-550-04-09  5 25-550-04-91  1 25-550-04-71  1  = 1.0 mm	25-551-04-09  5 25-551-04-91  1 25-551-04-71  1  = 1.0 mm	25-552-04-09  5 25-552-04-91  1 25-552-04-71  1  = 1.0 mm
Low profile	50-362-02-09  5 50-362-02-91  1 50-362-02-71  1  = 0.6 mm	50-360-04-09  5 50-360-04-91  1 50-360-04-71  1  = 0.6 mm	50-362-04-09  5 50-362-04-91  1 50-362-04-71  1  = 0.6 mm	50-364-04-09  5 50-364-04-91  1 50-364-04-71  1  = 0.6 mm
		25-550-04-89  1	25-551-04-89  1	25-552-04-89  1



Explanation of icons

-  Pure titanium
-  Plate profile
-  Units per package
-  Sterile packaged implants



Standard

25-550-06-09  5
 25-550-06-91  1
 25-550-06-71  1

 = 1.0 mm

25-552-06-09  5
 25-552-06-91  1
 25-552-06-71  1

 = 1.0 mm

25-550-08-09  5
 25-550-08-91  1
 25-550-08-71  1

 = 1.0 mm

Low profile

50-366-06-09  5
 50-366-06-91  1
 50-366-06-71  1

 = 0.6 mm

50-368-06-09  5
 50-368-06-91  1
 50-368-06-71  1

 = 0.6 mm


50-370-08-09  5
 50-370-08-91  1
 50-370-08-71  1

 = 0.6 mm

25-550-06-89  1

25-552-06-89  1

25-550-08-89  1



1/1



Standard

25-565-04-09 Ti 5
 25-565-04-91 Ti 1
 25-565-04-71 Ti 1
 ⌀ = 1.0 mm



25-567-04-09 Ti 5
 25-567-04-91 Ti 1
 25-567-04-71 Ti 1
 ⌀ = 1.0 mm



25-574-05-09 Ti 5
 25-574-05-91 Ti 1
 25-574-05-71 Ti 1
 ⌀ = 1.0 mm



25-575-05-09 Ti 5
 25-575-05-91 Ti 1
 25-575-05-71 Ti 1
 ⌀ = 1.0 mm

Low profile

50-379-04-09 Ti 5
 50-379-04-91 Ti 1
 50-379-04-71 Ti 1
 ⌀ = 0.6 mm

50-381-04-09 Ti 5
 50-381-04-91 Ti 1
 50-381-04-71 Ti 1
 ⌀ = 0.6 mm

50-380-05-09 Ti 5
 50-380-05-91 Ti 1
 50-380-05-71 Ti 1
 ⌀ = 0.6 mm

50-382-05-09 Ti 5
 50-382-05-91 Ti 1
 50-382-05-71 Ti 1
 ⌀ = 0.6 mm

25-567-04-89 Al 1

25-567-04-89 Al 1

1/1



Standard

25-576-06-09 Ti 5
 25-576-06-91 Ti 1
 25-576-06-71 Ti 1
 ⌀ = 1.0 mm



25-577-06-09 Ti 5
 25-577-06-91 Ti 1
 25-577-06-71 Ti 1
 ⌀ = 1.0 mm



25-578-06-09 Ti 5
 25-578-06-91 Ti 1
 25-578-06-71 Ti 1
 ⌀ = 1.0 mm



25-579-06-09 Ti 5
 25-579-06-91 Ti 1
 25-579-06-71 Ti 1
 ⌀ = 1.0 mm

Low profile

50-401-06-09 Ti 5
 50-401-06-91 Ti 1
 50-401-06-71 Ti 1
 ⌀ = 0.6 mm

50-402-06-09 Ti 5
 50-402-06-91 Ti 1
 50-402-06-71 Ti 1
 ⌀ = 0.6 mm

50-403-06-91 Ti 1
 50-403-06-71 Ti 1
 ⌀ = 0.6 mm

50-404-06-91 Ti 1
 50-404-06-71 Ti 1
 ⌀ = 0.6 mm

1/1



Standard

25-568-04-09 Ti 5
 25-568-04-91 Ti 1
 25-568-04-71 Ti 1
 ⌀ = 1.0 mm



25-568-06-91 Ti 1
 25-568-06-71 Ti 1
 ⌀ = 1.0 mm



Low profile 50-384-04-91 Ti 1
 ⌀ = 0.5 mm

50-389-06-91 Ti 1
 50-389-06-71 Ti 1
 ⌀ = 0.6 mm

50-392-05-91 Ti 1
 50-392-05-71 Ti 1
 ⌀ = 0.6 mm

1/1



Standard 25-560-06-09 Ti 5
 25-560-06-91 Ti 1
 25-560-06-71 Ti 1
 ⌀ = 1.0 mm



25-561-05-09 Ti 5
 25-561-05-91 Ti 1
 25-561-05-71 Ti 1
 ⌀ = 1.0 mm



25-562-05-09 Ti 5
 25-562-05-91 Ti 1
 25-562-05-71 Ti 1
 ⌀ = 1.0 mm



25-563-05-09 Ti 5
 25-563-05-91 Ti 1
 25-563-05-71 Ti 1
 ⌀ = 1.0 mm

Low profile 50-386-06-09 Ti 5
 50-386-06-91 Ti 1
 50-386-06-71 Ti 1
 ⌀ = 0.6 mm

50-386-05-09 Ti 5
 50-386-05-91 Ti 1
 50-386-05-71 Ti 1
 ⌀ = 0.6 mm

25-560-06-89 Al 1

1/1



Standard

25-553-05-09 Ti 5
25-553-05-91 Ti 1
25-553-05-71 Ti 1

= 1.0 mm



25-554-04-09 Ti 5
25-554-04-91 Ti 1
25-554-04-71 Ti 1

= 1.0 mm

Low profile

50-390-05-09 Ti 5
50-390-05-91 Ti 1
50-390-05-71 Ti 1

= 0.6 mm

25-553-05-89 Al 1

50-391-05-09 Ti 5
50-391-05-91 Ti 1
50-391-05-71 Ti 1

= 0.6 mm

25-554-04-89 Al 1

1/1



Standard

25-573-06-09 Ti 5
25-573-06-91 Ti 1
25-573-06-71 Ti 1

= 1.0 mm



25-570-06-09 Ti 5
25-570-06-91 Ti 1
25-570-06-71 Ti 1

= 1.0 mm

Low profile

50-396-06-09 Ti 5
50-396-06-91 Ti 1
50-396-06-71 Ti 1


= 0.6 mm

25-573-06-89 Al 1

50-397-06-09 Ti 5
50-397-06-91 Ti 1
50-397-06-71 Ti 1

= 0.6 mm

25-570-06-89 Al 1



50-381-05-09 Ti 1
 50-381-05-71 Ti 1
 = 0.7 mm



50-383-05-09 Ti 1
 50-383-05-71 Ti 1
 = 0.7 mm




50-405-07-09 Ti 1
 50-405-07-71 Ti 1
 = 0.7 mm



50-407-07-09 Ti 1
 50-407-07-71 Ti 1
 = 0.7 mm



50-405-09-09 Ti 1
 50-405-09-71 Ti 1
 = 0.7 mm



50-407-09-09 Ti 1
 50-407-09-71 Ti 1
 = 0.7 mm



50-405-08-09 Ti 1
 50-405-08-71 Ti 1
 = 0.7 mm



50-406-08-09 Ti 1
 50-406-08-71 Ti 1
 = 0.7 mm

Compression plates





50-410-04-09 Ti 5
 50-410-04-91 Ti 1
 50-410-04-71 Ti 1
 = 1.0 mm



50-410-05-91 Ti 1
 50-410-05-71 Ti 1
 = 1.0 mm



50-412-05-09 Ti 5
 50-412-05-91 Ti 1
 50-412-05-71 Ti 1
 = 1.0 mm



1/1



25-350-04-09 Ti 5
 25-350-04-91 Ti 1
 = 1.0 mm



25-351-04-09 Ti 5
 25-351-04-91 Ti 1
 = 1.0 mm



25-352-04-09 Ti 5
 25-352-04-91 Ti 1
 = 1.0 mm



25-350-06-09 Ti 5
 25-350-06-91 Ti 1
 = 1.0 mm

1/1



25-352-06-09 Ti 5
 25-352-06-91 Ti 1
 = 1.0 mm



25-350-08-09 Ti 5
 25-350-08-91 Ti 1
 = 1.0 mm



25-350-16-09 Ti 5
 25-350-16-91 Ti 1
 = 1.0 mm



25-351-16-09 Ti 5
 25-351-16-91 Ti 1
 = 0.8 mm



25-352-05-09 **Ti** 5
 25-352-05-91 **Ti** 1
 ⌀ = 1.0 mm



25-353-05-09 **Ti** 5
 25-353-05-91 **Ti** 1
 ⌀ = 1.0 mm



25-354-04-09 **Ti** 5
 25-354-04-91 **Ti** 1
 ⌀ = 1.0 mm



25-355-05-09 **Ti** 5
 25-355-05-91 **Ti** 1
 ⌀ = 1.0 mm



25-356-04-09 **Ti** 5
 25-356-04-91 **Ti** 1
 ⌀ = 1.0 mm



25-357-04-09 **Ti** 5
 25-357-04-91 **Ti** 1
 ⌀ = 1.0 mm



25-358-06-09 **Ti** 5
 25-358-06-91 **Ti** 1
 ⌀ = 1.0 mm



25-359-06-09 **Ti** 5
 25-359-06-91 **Ti** 1
 ⌀ = 1.0 mm



25-360-06-09 **Ti** 5
 25-360-06-91 **Ti** 1
 ⌀ = 1.0 mm



25-361-05-09 **Ti** 5
 25-361-05-91 **Ti** 1
 ⌀ = 1.0 mm



25-362-05-09 **Ti** 5
 25-362-05-91 **Ti** 1
 ⌀ = 1.0 mm



25-363-05-09 **Ti** 5
 25-363-05-91 **Ti** 1
 ⌀ = 1.0 mm

1/1



25-368-04-09 Ti 5
 25-368-04-91 Ti 1
 = 1.0 mm



25-370-06-09 Ti 5
 25-370-06-91 Ti 1
 = 1.0 mm



25-371-06-09 Ti 5
 25-371-06-91 Ti 1
 = 1.0 mm



25-372-04-09 Ti 5
 25-372-04-91 Ti 1
 = 1.0 mm

1/1



25-364-04-09 Ti 5
 25-364-04-91 Ti 1
 = 1.0 mm




25-366-04-09 Ti 5
 25-366-04-91 Ti 1
 = 1.0 mm



25-365-04-09 Ti 5
 25-365-04-91 Ti 1
 = 1.0 mm



25-367-04-09 Ti 5
 25-367-04-91 Ti 1
 = 1.0 mm



25-374-04-09 Ti 5
 25-374-04-91 Ti 1
 ⌀ = 1.0 mm



25-376-04-09 Ti 5
 25-376-04-91 Ti 1
 ⌀ = 1.0 mm



25-375-04-09 Ti 5
 25-375-04-91 Ti 1
 ⌀ = 1.0 mm



25-377-04-09 Ti 5
 25-377-04-91 Ti 1
 ⌀ = 1.0 mm



25-374-07-09 Ti 5
 25-374-07-91 Ti 1
 ⌀ = 1.0 mm



25-376-07-09 Ti 5
 25-376-07-91 Ti 1
 ⌀ = 1.0 mm

1/1



25-400-06-09 Ti 2
25-400-06-91 Ti 1
25-400-06-71 Ti 1

60 mm
 = 1.4 mm



25-400-09-09 Ti 2
25-400-09-91 Ti 1
25-400-09-71 Ti 1

90 mm
 = 1.4 mm



25-400-12-09 Ti 2
25-400-12-91 Ti 1
25-400-12-71 Ti 1

120 mm
 = 1.4 mm



25-400-15-09 Ti 2
25-400-15-91 Ti 1
25-400-15-71 Ti 1

150 mm
 = 1.4 mm



25-400-18-09 Ti 2
 25-400-18-91 Ti 1
 25-400-18-71 Ti 1

180 mm
 = 1.4 mm



25-400-21-09 Ti 2
 25-400-21-91 Ti 1
 25-400-21-71 Ti 1

210 mm
 = 1.4 mm



25-400-24-09 Ti 2
 25-400-24-91 Ti 1

240 mm
 = 1.4 mm

Implants: universal mesh plates




1/1

25-006-02-09  1

25-006-02-71  1

 = 0.6 mm




1/1

25-006-09-09  1

25-006-09-71  1

 = 0.6 mm



Explanation of icons

-  Pure titanium
-  Plate profile
-  Units per package
-  Sterile packaged implants


1/1




25-020-03-09  1
 25-020-03-71  1
 = 0.6 mm




25-020-04-09  1
 25-020-04-71  1
 = 0.6 mm



25-020-06-09  1
 25-020-06-71  1
 = 0.6 mm



25-020-08-09  1
 25-020-08-71  1
 = 0.6 mm




25-020-24-09  1
 25-020-24-71  1
 = 0.6 mm


1/1




25-024-04-09  1
 25-024-04-71  1
 = 0.6 mm



25-024-06-09  1
 25-024-06-71  1
 = 0.6 mm



25-024-16-09  1
 25-024-16-71  1
 = 0.6 mm



25-024-24-09  1
 25-024-24-71  1
 = 0.6 mm



25-022-04-09  1
 25-022-04-71  1
 = 0.6 mm



25-022-08-09  1
 25-022-08-71  1
 = 0.6 mm

Implants: Rigid Blue plates

A series of mini and mesh plates using 2.0- and 2.3-mm screws. These plates are made of high strength titanium alloy (Ti-6AL-4V ASTM F136a) to provide high strength, low weight and a flat profile.

The plates can be used with both Centre-Drive® as well as maxDrive® screws in monocortical osteosynthesis without compression. The Rigid Blue plates can be identified by their blue anodized surface.

$\frac{1}{1}$



50-321-04-09 1

50-321-04-71 1

= 1.0 mm



50-322-04-09 1

50-322-04-71 1

= 1.0 mm



50-320-08-09 1

50-320-08-71 1

= 1.0 mm



50-406-17-09 1

50-406-17-71 1

= 1.0 mm



Explanation of icons

- Ti** Titanium alloy
-  Plate profile
- 1** Units per package

STERILE IR Sterile packaged implants

$\frac{1}{1}$



50-330-04-09 **Ti** **1**

 = 1.0 mm



50-330-06-09 **Ti** **1**

50-330-06-71 **Ti** **1**

 = 1.0 mm



50-331-06-09 **Ti** **1**

50-331-06-71 **Ti** **1**

 = 1.0 mm



50-332-06-09 **Ti** **1**

50-332-06-71 **Ti** **1**

 = 1.0 mm



25-022-28-09 **Ti** **1**

 = 1.0 mm

Implants:
alternative plates for the treatment of fractures
in the midface



1/1

50-735-04-09  

 = 1.0 mm

with 1-mm graduation



50-735-14-09  

 = 1.0 mm

with 1-mm graduation



50-735-24-09  

 = 1,0 mm

with 1-mm graduation



1/1

50-735-34-09  

 = 1.0 mm

with 1-mm graduation



50-735-36-09  

 = 1.0 mm

with 1-mm graduation



Explanation of icons

-  Pure titanium
-  Plate profile
-  Units per package




$\frac{1}{1}$
50-740-05-09  
 = 1.0 mm



100°
50-741-05-09  
 = 1.0 mm



100°
50-740-09-09  
 = 1.0 mm



100°
50-741-09-09  
 = 1.0 mm



$\frac{1}{1}$
50-742-04-09  
 = 1.0 mm



50-743-04-09  
 = 1.0 mm



50-756-05-09  
 = 1.0 mm



50-756-07-09  
 = 1.0 mm



50-757-07-09  
 = 1.0 mm

Instruments: screwdrivers and blades

Screwdrivers and blades

										
			Centre-Drive® maxDrive®							
Blades to fit		Length	25-402-99-07*	25-407-03-04 flattened	25-407-04-04* flattened	25-407-05-04*	25-406-99-07	50-425-05-07* rigid, flattened	25-410-00-07 Ratchet	25-411-00-07* Ratchet, flattened
25-434-98-07	 2.0/2.3 mm	80 mm	x							
25-540-98-07	 2.0/2.3 mm	94 mm		x	x	x	x	x	x	x
25-491-97-07	 2.0/2.3 mm	80 mm	x							
25-486-97-07	 2.0/2.3 mm	80 mm		x	x	x	x	x	x	x
25-486-98-07	 2.0/2.3 mm	94 mm		x	x	x	x	x	x	x
25-480-20-07	 2.0/2.3 mm	103 mm		x	x	x	x	x	x	x

* For storage in the LevelOne modules



Explanation of icons

- St** Stainless steel
- maxDrive®
- Centre-Drive®

Measuring clips



25-650-03-04
Measuring clip for screw length, red



25-651-01-04
Measuring clip for screw diameter



$\frac{1}{4}$
25-542-17-07 2.0/2.3 mm
Screwdriver, rigid, not detachable



$\frac{1}{4}$
25-486-17-07 2.0/2.3 mm
Screwdriver, rigid, not detachable




$\frac{1}{2}$
25-422-20-07 2.0 mm
25-424-20-07 2.0 mm
Bone graft screwdriver



$\frac{1}{2}$
25-480-20-07 2.0/2.3 mm
Retention sleeve for screw explanation

Instruments: twist drills

Twist drills				
	J-notch attachment			
	Ø x L (mm)	Stop (mm)	STERILE R	
	1.5 x 50 mm	5	25-449-05-07	5
	1.5 x 50 mm	5	25-449-05-91	1
	1.5 x 50 mm	5	25-449-05-71	1
	1.5 x 50 mm	7	25-449-07-07	5
	1.5 x 50 mm	7	25-449-07-91	1
	1.5 x 50 mm	7	25-449-07-71	1
	1.5 x 50 mm	9	25-449-09-07	5
	1.5 x 50 mm	9	25-449-09-91	1
	1.5 x 50 mm	9	25-449-09-71	1
	1.5 x 50 mm	11	25-449-11-07	5
	1.5 x 50 mm	11	25-449-11-91	1
	1.5 x 50 mm	21	25-449-16-07	5
	1.5 x 50 mm	21	25-449-16-91	1
	1.5 x 50 mm	21	25-449-16-71	1
	1.5 x 70 mm	10	25-458-16-07	5
	1.5 x 70 mm	10	25-458-16-91	1
	1.5 x 70 mm	–	25-461-16-07	5
	1.5 x 70 mm	–	25-461-16-71	1
	1.5 x 70 mm	–	25-461-16-91	1
	1.5 x 105 mm	5	50-515-05-07	1
	1.5 x 115 mm	–	25-469-11-07	1
	1.5 x 115 mm	–	25-469-11-71	1



Explanation of icons


- St** Stainless steel
- J** J-notch attachment
- I** Dental attachment
- 1** Units per package

STERILE IR Sterile packaged implants

Twist drills

St I	Dental attachment			
	Ø x L (mm)	Stop (mm)		
	1.5 x 34 mm	19	50-924-00-07	1
	1.5 x 20 mm	7	50-924-07-07	1
	1.5 x 22 mm	9	50-924-09-07	1
	1.5 x 45 mm	5	25-459-05-07	1
	1.5 x 45 mm	7	25-459-07-07	1
	1.5 x 45 mm	15	25-459-15-07	1


Instruments



25-496-15-09 **Ti** 1
15 cm / 6"
Plate-holding forceps



25-441-18-07 **St** 1
18 cm / 7 1/4"
Plate-holding forceps, curved



Lindorf
25-435-20-07 **St** 1
Mini
16 cm / 6 1/4"
Plate-holding instrument




Lindorf
25-435-25-07 **St** 1
Champy
16 cm / 6 1/4"
Plate-holding instrument



Explanation of icons

- St** Stainless steel
- Ti** Titanium alloy
- 1** Units per package



50-502-10-07 **St 1**

16 cm / 6 1/4"

Multifunctional plate bender for 1.0/1.5/2.0 systems
for 3D bending, 90° bending, cutting

50-502-98-07 **St 1**

Set of spare blades for 50-502-10-07



50-502-11-07 **St 1**

20 cm / 8"

Multifunctional plate cutter for 1.0/1.5/2.0 systems

50-502-99-07 **St 1**

Set of spare blades for 50-502-11 -07

Instruments



25-516-13-07 St 1

13 cm / 5 1/4"

Modeling pliers



25-510-12-07 St 1

Mini

12 cm / 4 1/4"

Bending pliers KLS Martin Mini System




25-416-15-07 St 1

Champy

15 cm / 6"

Bending pliers Champy System



1/2

50-501-40-07 St 1

17 cm / 6 1/2"

Depth gauge




1/2

25-420-18-07 St 1 TC GOLD

18 cm / 7"

Plate cutting pliers



1/2

50-325-01-07 St 1

15 cm / 6"

Fracture reduction clamp



1/2



25-471-15-07
Trocar




25-471-20-07
Drill sleeve

25-471-00-07 St 1


20 cm / 8"

Transbuccal system for 2.0- and 2.3-mm screws,
complete (incl. trocar and drill sleeve)

1/2



1/2



25-474-00-01 St 1

Supplementary instrument for transbuccal procedures,
complete set


Storage: implants



55-961-06-04
2.0-mm Mini module complete, without implants

2.0-mm Mini module

55-961-06-04	2.0-mm Mini module complete consisting of:
55-962-03-04	Insert module with red module rails
55-963-06-04	Cover for 2.0-mm Mini module
55-964-05-04	Insert for 2.0-mm Mini Universal
55-964-03-04	Insert for 2.0-mm Mini Trauma
	Insert 55-964-05-04 is prepared for the following implants:
25-550-04-09	Mini-plate, 4-hole, straight
25-550-08-09	Mini-plate, 8-hole, straight
25-552-04-09	Mini plate, 4-hole, straight, long
50-368-16-09	Mini plate, 16-hole, 0.6 mm, straight
50-375-04-09	Mini-plate, 4-hole, L-shape, right, medium, 0.6 mm
50-377-04-09	Mini-plate, 4-hole, L-shape, left, medium, 0.6 mm
50-379-04-09	Mini-plate, 4-hole, L-shape, right, long, 0.6 mm
50-381-04-09	Mini-plate, 4-hole, L-shape, left, long, 0.6 mm
50-386-06-09	Mini-plate, 6-hole, T-shape, 0.6 mm
50-390-05-09	Mini-plate, 5-hole, Y-shape, short, 0.6 mm
50-396-06-09	Mini-plate, 6-hole, double Y-shape, short, 0.6 mm
50-408-04-09	Mini-plate, 4-hole, Z-shape, left, 100°, long, 0.6 mm
50-409-04-09	Mini-plate, 4-hole, Z-shape, right, 100°, long, 0.6 mm
	Insert 55-964-03-04 is prepared for the following implants:
50-381-05-09	Mini-plate, 5-hole, L-shape, right, long
50-383-05-09	Mini-plate, 5-hole, L-shape, left, long
50-405-07-09	Mini-plate, 7-hole, L-shape, right, 100°
50-405-08-09	Mini-plate, 8-hole, Y-shape, 0.7 mm
50-405-09-09	Mini-plate, 9-hole, L-shape, right
50-405-12-09	Mini orbita plate, 12-hole, 0.7 mm
50-405-17-09	Mini-plate, 20-hole, straight
50-406-08-09	Mini-plate, 8-hole, double Y-shape, 0.7 mm
50-407-07-09	Mini-plate, 7-hole, L-shape, left, 100°
50-407-09-09	Mini-plate, 9-hole, L-shape, left



55-962-13-04

2.0-mm storage module without lid, without contents

55-963-09-04

Lid



55-962-23-04

2.0-mm screw storage module without lid



55-962-33-04

Double insert module

55-963-08-04

Lid

Storage: implants



55-961-07-04

2.0-mm Champy module complete, without implants

2.0-mm Champy module

55-961-07-04	Champy module complete consisting of:
55-962-03-04	Insert module with red module rails
55-963-07-04	Cover for 2.0-mm Champy module
55-964-18-04	Insert for 2.0-mm Champy module
55-964-17-04	Insert, universal
	Insert 55-964-18-04 is prepared for the following implants:
25-350-04-09	Champy plate, 4-hole, straight, short
25-350-06-09	Champy plate, 6-hole, straight, short
25-350-08-09	Champy plate, 8-hole, straight
25-351-04-09	Champy plate, 4-hole, straight, short
25-352-04-09	Champy plate, 4-hole, straight, long
25-352-05-09	Champy plate, Y-shape, 5-hole, long
25-352-06-09	Champy plate, 6-hole, straight, long
25-356-04-09	Champy plate, T-shape, 2 x 2 holes, short
25-357-04-09	Champy plate, T-shape, 2 x 2 holes, long
25-360-06-09	Champy plate, T-shape, 3 x 3 holes, short
25-364-04-09	Champy plate, L-shape, 4-hole, short, right
25-365-04-09	Champy plate, L-shape, 4-hole, long, right
25-366-04-09	Champy plate, L-shape, 4-hole, short, left
25-367-04-09	Champy plate, L-shape, 4-hole, long, left
25-368-04-09	Champy plate, H-shape, 6-hole
25-395-08-09	Champy orbita plate, 8-hole



55-962-13-04

2.0-mm storage module without lid, without contents

55-963-09-04

Lid




55-962-23-04

2.0-mm screw storage module without lid

55-963-08-04

Lid



55-962-13-04

Storage module without lid

Storage: implants



55-961-32-04

2.0-mm Arnett FAB module complete, without implants

2.0-mm Arnett FAB module

55-961-32-04	Arnett module complete consisting of:
55-962-03-04	Insert module with red module rails
55-963-52-04	Lid for Arnett module
55-964-69-04	Insert for Arnett plates "Type 1"
55-964-80-04	Insert for Arnett plates "Type 2"
	Insert 55-964-69-04 is prepared for the following implants:
25-396-29-09	Arnett plate, 4-hole, straight, small
25-396-34-09	Arnett plate, 4-hole, straight, medium
25-396-39-09	Arnett plate, 4-hole, straight, large
25-396-44-09	Arnett plate, 4-hole, straight, XL
25-396-56-09	Arnett plate, 6-hole, straight, large
25-399-26-09	Arnett plate, T-shape, 0.8 mm, small
25-399-32-09	Arnett plate, T-shape, 0.8 mm, medium
25-399-37-09	Arnett plate, T-shape, 1.0 mm, large
50-304-02-09	Arnett chin plate, 2 mm
50-304-03-09	Arnett chin plate, 3 mm
50-304-04-09	Arnett chin plate, 4 mm
50-304-06-09	Arnett chin plate, 6 mm
50-304-08-09	Arnett chin plate, 8 mm
50-304-10-09	Arnett chin plate, 10 mm
	Insert 55-964-80-04 is prepared for the following implants:
25-397-27-09	Arnett plate, L-shape, 0.8 mm, small
25-397-31-09	Arnett plate, L-shape, 0.8 mm, medium
25-397-35-09	Arnett plate, L-shape, 0.8 mm, large
25-397-39-09	Arnett plate, L-shape, 1.0 mm, XL
25-398-27-09	Arnett plate, C-shape, 0.8 mm, small
25-398-31-09	Arnett plate, C-shape, 0.8 mm, medium
25-398-35-09	Arnett plate, C-shape, 0.8 mm, large
25-398-39-09	Arnett plate, C-shape, 1.0 mm, XL



55-962-23-04
2.0-mm screw storage module without lid

55-963-08-04
Lid




55-962-33-04
Double insert module

Storage: instruments




55-969-26-04

Multi storage*, configured, with lid, without instruments
24 x 13 cm



55-969-16-04

Free storage, with lid, without silicone mat
24 x 13 cm



55-969-14-04

Free storage, with lid, without silicone mat
24 x 24 cm



Figure: possible combination with system 2.3 Fracture and Reconstruction

Instrument modules		
55-969-16-04	24 x 13 cm	Free storage, with lid
55-969-14-04	24 x 24 cm	Free storage, with lid
55-969-26-04	24 x 13 cm	Multi*, configured, with lid
Individual components		
55-964-15-04	24 x 13 cm	Free storage, without lid
55-964-30-04	24 x 24 cm	Free storage, without lid
55-969-15-04	24 x 13 cm	Multi*, configured, without lid
55-969-47-04	24 x 13 cm	Basic*, configured, without lid
55-969-30-04	24 x 24 cm	Multi*, configured, without lid
55-969-49-04	24 x 24 cm	Basic*, configured, without lid
55-963-15-04	24 x 13 cm	Multi lid
55-963-16-04	24 x 13 cm	Basic lid
55-963-30-04	24 x 24 cm	Lid
55-969-94-04	24 x 13 cm	Silicone mat for 55-969-16-04
55-969-93-04	24 x 24 cm	Silicone mat for 55-969-14-04
Module clips		
55-969-98-04		to hold 3 racks, left
55-969-99-04		to hold 3 racks, right

* Multi = Storage prepared for LevelOne multifunctional instruments

* Basic = Storage prepared for traditional instruments

Storage: set suggestions

Set suggestions



55-969-26-04		Instrument tray Multi, configured, 24 x 13 cm prepared for:
25-402-99-07	②	Screwdriver handle, medium
25-406-99-07	②	Screwdriver handle, large
25-435-10-07	①	Lindorf plate-holding instrument 1.0 mm
25-435-15-07	①	or Lindorf plate-holding instrument 2.0 mm
25-435-20-07	①	or Lindorf plate-holding instrument 2.0 mm
25-435-25-07	①	or Lindorf plate-holding instrument 2.0 mm Champy
50-502-10-07	②	Multifunctional pliers for 1.0/2.0/2.0 mm
50-502-11-07	①	Cutter for 1.0/2.0/2.0 mm

Set suggestions



55-969-57-04		Instrument tray Basic, configured, 24 x 13 cm prepared for:
25-402-99-07	②	Screwdriver handle, medium
25-406-99-07	②	Screwdriver handle, large
25-420-16-07	①	Plate cutting pliers
25-496-15-09	①	Plate-holding forceps
25-510-12-07	①	Bending pliers KLS Martin Mini System
25-516-13-07	②	Modeling pliers



Figure: possible combination with system 2.3 Fracture and Reconstruction

Set suggestions



55-969-24-04		Instrument tray Multi, configured, 24 x 24 cm prepared for:
25-402-99-07	②	Screwdriver handle, medium
25-406-99-07	②	Screwdriver handle, large
25-496-15-09	①	Plate-holding forceps
50-205-17-07	②	Plate-holding forceps
50-325-01-07	①	Fracture clamp
50-501-40-07	①	Depth gauge
50-502-10-07	②	Multifunctional pliers for 1.0/1.5/2.0 mm
50-502-11-07	①	Cutter for 1.0/1.5/2.0 mm
50-504-23-07	①	Multifunctional instrument 1.5/2.0/3.0 mm type A
50-504-24-07	①	Multifunctional instrument 1.5/2.0/3.0 mm type B

Set suggestions



55-969-59-04		Instrument tray Basic, configured, 24 x 24 cm prepared for:
25-402-99-07	②	Screwdriver handle, medium
25-406-99-07	②	Screwdriver handle, large
25-420-18-07	①	Plate cutting pliers
25-496-15-09	①	Plate-holding forceps
25-510-12-07	①	Bending pliers KLS Martin Mini System
25-516-13-07	②	Modeling pliers
50-205-17-07	②	Plate-holding forceps
50-325-01-07	①	Fracture clamp
50-504-23-07	①	Multifunctional instrument 1.5/2.0/3.0 mm type A
50-504-24-07	①	Multifunctional instrument 1.5/2.0/3.0 mm type B

Sterile kits

The use of sterile packaged implants with all the associated advantages entails adapting the usual workflow. The selection of implants in particular is changed and the time required for passing these on at the operating table is increased, as the individually sterile packaged implants must first be opened before use.

There are specific indications, where the majority of cases are treated in the same manner and using the same type of implants. For example, 4- or 6-hole plates with a total of four Ø 2.0 x 5 mm-screws are often used for treating fractures in the mandible.

To optimize the workflow for passing on sterile packaged implants and to maintain efficiency in the OR compared to using an open implant system, it makes sense to offer pre-defined sterile kits containing all the plates, screws and drills required for a standard restoration. For this reason, the following kits will be available in the future.



Sterile kits

25-550-45-71	Sterile kit Mini 2 x Mini plate 2.0, 4-hole, straight (25-550-04) 8 x maxDrive® screw 2.0 x 5 mm (25-872-05) 1 x Twist drill (25-449-05)
25-550-46-71	Sterile kit Mini 2 x Mini plate 2.0, 6-hole, straight (25-550-06) 8 x maxDrive® screw 2.0 x 5 mm (25-872-05) 1 x Twist drill (25-449-05)



11

11-865-12-07.....13

22

22-500-11-07.....13

25

25-006-02-09.....42

25-006-02-71.....42

25-006-09-09.....42

25-006-09-71.....42

25-020-03-09.....43

25-020-03-71.....43

25-020-04-09.....43

25-020-04-71.....43

25-020-06-09.....43

25-020-06-71.....43

25-020-08-09.....43

25-020-08-71.....43

25-020-24-09.....43

25-020-24-71.....43

25-022-04-09.....43

25-022-04-71.....43

25-022-08-09.....43

25-022-08-71.....43

25-022-28-09.....45

25-023-04-09.....35

25-023-04-71.....35

25-024-04-09.....43

25-024-04-71.....43

25-024-06-09.....43

25-024-06-71.....43

25-024-16-09.....43

25-024-16-71.....43

25-024-24-09.....43

25-024-24-71.....43

25-092-08-09.....10

25-092-08-72.....10

25-092-12-09.....10

25-092-12-72.....10

25-092-38-05.....10

25-092-38-72.....10

25-092-42-05.....10

25-092-42-72.....10

25-095-90-09.....12

25-095-91-09.....12

25-097-08-09.....11

25-097-08-75.....11

25-099-12-09.....10

25-099-12-72.....10

25-099-14-09.....10

25-099-14-72.....10

25-099-42-05.....10

25-099-42-72.....10

25-099-44-05.....10

25-099-44-72.....10

25-240-04-09.....35

25-240-06-09.....35

25-241-04-09.....35

25-241-06-09.....35

25-283-05-09.....19

25-283-05-71.....19

25-283-15-71.....18

25-283-15-91.....18

25-283-20-71.....21

25-283-20-91.....21

25-283-21-71.....21

25-283-21-91.....21

25-283-25-71.....20

25-283-25-91.....20

25-283-26-71.....20

25-283-26-91.....20

25-285-05-09.....19

25-285-05-71.....19

25-285-10-91.....18

25-288-08-09.....21

25-289-08-09.....21

25-350-04-09.....36, 58

25-350-04-91.....36

25-350-06-09.....36, 58

25-350-06-91.....36

25-350-08-09.....36, 58

25-350-08-91.....36

25-350-16-09.....36

25-350-16-91.....36

25-351-04-09.....36, 58

25-351-04-91.....36

25-351-16-09.....36

25-351-16-91.....36

25-352-04-09.....36, 58

25-352-04-91.....36

25-352-05-09.....37, 58

25-352-05-91.....37

25-352-06-09.....36, 58

25-352-06-91.....36

25-353-05-09.....37

25-353-05-91.....37

25-354-04-09.....37

25-354-04-91.....37

25-355-05-09.....37

25-355-05-91.....37

25-356-04-09.....37, 58

25-356-04-91.....37

25-357-04-09.....37, 58

25-357-04-91.....37

25-358-06-09.....37

25-358-06-91.....37

25-359-06-09.....37

25-359-06-91.....37

25-360-06-09.....37, 58

25-360-06-91.....37

25-361-05-09.....37

25-361-05-91.....37

25-362-05-09.....37

25-362-05-91.....37

25-363-05-09.....37

25-363-05-91.....37

25-364-04-09.....38, 58

25-364-04-91.....38

25-365-04-09.....38, 58

25-365-04-91.....38

25-366-04-09.....38, 58

25-366-04-91.....38

25-367-04-09.....38, 58

25-367-04-91.....38

25-368-04-09.....38, 58

25-368-04-91.....38

25-370-06-09.....38

25-370-06-91.....38

25-371-06-09.....38

25-371-06-91.....38

25-372-04-09.....38

25-372-04-91.....38

25-374-04-09.....39

25-374-04-91.....39

25-374-07-09.....39

25-374-07-91.....39

25-375-04-09.....39

25-375-04-91.....39

25-376-04-09.....39

25-376-04-91.....39

25-376-07-09.....39

25-376-07-91.....39

25-377-04-09.....39

25-377-04-91.....39

25-380-04-09.....20

25-380-04-71.....20

25-380-04-91.....20

25-382-04-09.....20

25-382-04-71.....20

25-382-04-91.....20

25-392-04-09.....24

25-392-04-71.....24

25-392-04-91.....24

25-394-04-09.....24

25-394-04-71.....24

25-394-04-91.....24

25-394-34-09.....23

25-395-08-09.....17, 58

25-395-08-91.....17

25-395-10-09.....17

25-395-10-91.....17

25-395-12-09.....17

25-395-12-91.....17

25-395-34-09.....23

25-396-29-09.....27, 60

25-396-29-71.....27

25-396-29-89.....27

25-396-34-09.....27, 60

25-396-34-71.....27

25-396-34-89.....27

25-396-39-09.....27, 60

25-396-39-71.....27

25-396-39-89.....27

25-396-44-09.....27, 60

25-396-44-71.....27

25-396-44-89.....27

25-396-56-09.....27, 60

25-396-56-71.....27

25-396-56-89.....27

25-397-27-09.....26, 60

25-397-27-71.....26

25-397-27-89.....26

25-397-31-09.....26, 60

25-397-31-71.....26

25-397-31-89.....26

25-397-35-09.....26, 60

25-397-35-71.....26

25-397-35-89.....26

25-397-39-09.....26, 60

25-397-39-71.....26

25-397-39-89.....26

25-398-27-09.....26, 60

25-398-27-71.....26

25-398-27-89.....26

25-398-31-09.....26, 60

25-398-31-71.....26

25-398-31-89.....26

25-398-35-09.....26, 60

25-398-35-71.....26

25-398-35-89.....26

25-398-39-09.....26, 60

25-398-39-71.....26

25-398-39-89.....26

25-399-26-09.....26, 60

25-399-26-71.....26

25-399-26-89.....26

25-399-32-09.....26, 60

25-399-32-71.....26

25-399-32-89.....26

25-399-37-09.....26, 60

25-399-37-71.....26

25-399-37-89.....26

25-400-06-09.....40

25-400-06-71.....40

25-400-06-91.....40

25-400-09-09.....40

25-400-09-71.....40

25-400-09-91.....40

25-400-12-09.....40

25-400-12-71.....40

25-400-12-91.....40

25-400-15-09.....40

25-400-15-71.....40

25-400-15-91.....40

25-400-18-09.....41

25-400-18-71.....41

25-400-18-91.....41

25-400-21-09.....41

25-400-21-71.....41

25-400-21-91.....41

25-400-24-09.....41

25-400-24-91.....41

25-401-30-09.....24

25-401-30-71.....24

25-401-35-09.....24

25-401-35-71.....24

25-401-36-09.....22

25-401-37-09.....22

25-401-38-09.....22

25-401-39-09.....22

25-401-40-09.....24

25-401-40-71.....24

25-401-43-09.....23

25-401-44-09.....23

25-402-99-07.....48, 64, 65

25-406-99-07.....48, 64, 65

25-407-03-04.....48

25-407-04-04.....13, 48

25-407-05-04.....48

25-410-00-07.....48

25-411-00-07.....48

25-416-15-07.....54

25-420-16-07.....64

25-420-18-07.....54, 65

25-422-20-07.....49

25-424-20-07.....49

25-434-98-07.....48

25-435-10-07.....64

25-435-15-07.....64

25-435-20-07.....52, 64

25-435-25-07.....52, 64

25-441-18-07.....52

25-449-05-07.....50

25-449-05-71.....50

25-449-05-91.....50

25-449-07-07.....50

25-449-07-71.....50

25-449-07-91.....50

25-449-09-07.....50

25-449-09-71.....50

25-449-09-91.....50

25-449-11-07.....50

25-449-11-91.....50

25-449-16-07.....50

25-449-16-71.....50

25-449-16-91.....50

25-458-16-07.....50

25-458-16-91.....50

25-459-05-07.....51

25-459-07-07.....51</

25-550-04-91.....	28	25-566-04-91.....	30	25-662-19-75.....	9	25-882-05-09.....	8	50-304-04-09.....	27, 60
25-550-06-09.....	29	25-567-04-09.....	31	25-663-45-09.....	9	25-882-05-71.....	8	50-304-04-71.....	27
25-550-06-71.....	29	25-567-04-71.....	31	25-663-47-09.....	9	25-882-05-75.....	8	50-304-06-09.....	27, 60
25-550-06-89.....	29	25-567-04-89.....	31	25-663-47-75.....	9	25-882-07-09.....	8	50-304-06-71.....	27
25-550-06-91.....	29	25-567-04-91.....	31	25-663-49-09.....	9	25-882-07-71.....	8	50-304-08-09.....	27, 60
25-550-08-09.....	29, 56	25-568-04-09.....	32	25-663-49-75.....	9	25-882-07-75.....	8	50-304-08-71.....	27
25-550-08-71.....	29	25-568-04-71.....	32	25-669-05-09.....	9	25-882-09-09.....	8	50-304-10-09.....	27, 60
25-550-08-89.....	29	25-568-04-91.....	32	25-669-05-75.....	9	25-882-09-71.....	8	50-304-10-71.....	27
25-550-16-09.....	29	25-568-06-71.....	32	25-669-07-09.....	9	25-882-09-75.....	8	50-311-17-09.....	15
25-550-16-71.....	30	25-568-06-91.....	32	25-669-07-75.....	9	25-882-11-09.....	8	50-311-17-71.....	15
25-550-16-89.....	30	25-569-04-09.....	30	25-872-04-09.....	8	25-882-11-71.....	8	50-311-27-09.....	15
25-550-16-91.....	30	25-569-04-71.....	30	25-872-04-71.....	8	25-882-11-75.....	8	50-311-27-71.....	15
25-550-40-09.....	30	25-569-04-91.....	30	25-872-04-75.....	8	25-882-13-09.....	8	50-320-08-09.....	44
25-550-40-71.....	30	25-570-06-09.....	33	25-872-05-09.....	8	25-882-13-71.....	8	50-320-08-71.....	44
25-550-45-71.....	67	25-570-06-71.....	33	25-872-05-71.....	8	25-882-13-75.....	8	50-321-04-09.....	44
25-551-02-09.....	28	25-570-06-89.....	33	25-872-05-75.....	8	25-882-15-09.....	8	50-321-04-71.....	44
25-551-02-71.....	28	25-570-06-91.....	33	25-872-06-09.....	8	25-882-15-71.....	8	50-322-04-09.....	44
25-551-02-91.....	28	25-571-04-09.....	30	25-872-06-71.....	8	25-882-15-75.....	8	50-322-04-71.....	44
25-551-04-09.....	28	25-571-04-71.....	30	25-872-06-75.....	8	25-882-17-09.....	8	50-325-01-07.....	54, 65
25-551-04-71.....	28	25-571-04-89.....	30	25-872-07-09.....	8	25-882-17-71.....	8	50-330-04-09.....	45
25-551-04-89.....	28	25-571-04-91.....	30	25-872-07-71.....	8	25-882-17-75.....	8	50-330-06-09.....	45
25-551-04-91.....	28	25-573-06-09.....	33	25-872-07-75.....	8	25-882-19-09.....	8	50-330-06-71.....	45
25-552-04-09.....	28, 56	25-573-06-71.....	33	25-872-09-09.....	8	25-882-19-71.....	8	50-331-06-09.....	45
25-552-04-71.....	28	25-573-06-91.....	33	25-872-09-71.....	8	25-882-19-75.....	8	50-331-06-71.....	45
25-552-04-89.....	28	25-574-05-09.....	31	25-872-09-75.....	8			50-332-06-09.....	45
25-552-04-91.....	28	25-574-05-71.....	31	25-872-11-09.....	8			50-332-06-71.....	45
25-552-06-09.....	29	25-574-05-91.....	31	25-872-11-71.....	8			50-360-04-09.....	28
25-552-06-71.....	29	25-575-05-09.....	31	25-872-11-75.....	8			50-360-04-71.....	28
25-552-06-89.....	29	25-575-05-71.....	31	25-872-13-09.....	8			50-360-04-91.....	28
25-552-06-91.....	29	25-575-05-91.....	31	25-872-13-71.....	8			50-362-02-09.....	28
25-553-05-09.....	33	25-576-06-09.....	31	25-872-13-75.....	8			50-362-02-71.....	28
25-553-05-71.....	33	25-576-06-71.....	31	25-872-15-09.....	8			50-362-02-91.....	28
25-553-05-89.....	33	25-576-06-91.....	31	25-872-15-71.....	8			50-362-04-09.....	28
25-553-05-91.....	33	25-577-06-09.....	31	25-872-15-75.....	8			50-362-04-71.....	28
25-554-04-09.....	33	25-577-06-71.....	31	25-872-17-09.....	8			50-362-04-91.....	28
25-554-04-71.....	33	25-577-06-91.....	31	25-872-17-71.....	8			50-364-04-09.....	28
25-554-04-89.....	33	25-578-06-09.....	31	25-872-17-75.....	8			50-364-04-71.....	28
25-554-04-91.....	33	25-578-06-71.....	31	25-872-19-09.....	8			50-364-04-91.....	28
25-560-06-09.....	32	25-578-06-91.....	31	25-872-19-71.....	8			50-366-06-09.....	29
25-560-06-71.....	32	25-579-06-09.....	31	25-873-44-09.....	8			50-366-06-71.....	29
25-560-06-89.....	32	25-579-06-71.....	31	25-873-44-71.....	8			50-366-06-91.....	29
25-560-06-91.....	32	25-579-06-91.....	31	25-873-44-75.....	8			50-368-06-09.....	29
25-561-05-09.....	32	25-650-03-04.....	49	25-873-45-09.....	8			50-368-06-71.....	29
25-561-05-71.....	32	25-651-01-04.....	49	25-873-45-71.....	8			50-368-06-91.....	29
25-561-05-91.....	32	25-662-04-09.....	9	25-873-45-75.....	8			50-368-16-09.....	30, 56
25-562-05-09.....	32	25-662-04-75.....	9	25-873-47-09.....	8			50-368-16-71.....	30
25-562-05-71.....	32	25-662-05-09.....	9	25-873-47-71.....	8			50-368-16-91.....	30
25-562-05-91.....	32	25-662-05-75.....	9	25-873-47-75.....	8			50-370-08-09.....	29
25-563-05-09.....	32	25-662-07-09.....	9	25-873-49-09.....	8			50-370-08-71.....	29
25-563-05-71.....	32	25-662-07-75.....	9	25-873-49-71.....	8			50-370-08-91.....	29
25-563-05-91.....	32	25-662-09-09.....	9	25-879-05-09.....	8			50-371-04-09.....	30
25-564-04-09.....	30	25-662-09-75.....	9	25-879-05-71.....	8			50-371-04-71.....	30
25-564-04-71.....	30	25-662-11-09.....	9	25-879-05-75.....	8			50-371-04-91.....	30
25-564-04-91.....	30	25-662-11-75.....	9	25-879-06-09.....	8			50-373-04-09.....	30
25-565-04-09.....	31	25-662-13-09.....	9	25-879-06-71.....	8			50-373-04-71.....	30
25-565-04-71.....	31	25-662-13-75.....	9	25-879-07-09.....	8			50-373-04-91.....	30
25-565-04-91.....	31	25-662-15-09.....	9	25-879-07-71.....	8			50-375-04-09.....	30, 56
25-566-04-09.....	30	25-662-15-75.....	9	25-879-07-75.....	8			50-375-04-91.....	30
25-566-04-71.....	30	25-662-17-09.....	9	25-879-09-09.....	8			50-377-04-09.....	30, 56
25-566-04-89.....	30	25-662-17-75.....	9	25-879-09-71.....	8			50-377-04-71.....	30
		25-662-19-09.....	9	25-879-09-75.....	8			50-377-04-91.....	30

39

39-311-20-07.....13

50

50-205-17-07.....	65
50-300-02-09.....	25
50-300-02-71.....	25
50-300-02-91.....	25
50-300-04-09.....	25
50-300-04-71.....	25
50-300-04-91.....	25
50-300-06-09.....	25
50-300-06-71.....	25
50-300-06-91.....	25
50-300-08-09.....	25
50-300-08-71.....	25
50-300-08-91.....	25
50-300-10-09.....	25
50-300-10-71.....	25
50-300-10-91.....	25
50-300-12-71.....	25
50-300-12-91.....	25
50-300-14-71.....	25
50-300-14-91.....	25
50-303-02-09.....	25
50-303-03-09.....	25
50-303-04-09.....	25
50-303-05-09.....	25
50-303-06-09.....	25
50-303-07-09.....	25
50-303-08-09.....	25
50-303-09-09.....	25
50-304-02-09.....	27, 60
50-304-02-71.....	27
50-304-03-09.....	27, 60
50-304-03-71.....	27

50-379-04-09.....	31, 56	50-403-07-91.....	35	50-741-05-09.....	47
50-379-04-71.....	31	50-404-06-71.....	31	50-741-09-09.....	47
50-379-04-91.....	31	50-404-06-91.....	31	50-742-04-09.....	47
50-380-05-09.....	31	50-404-07-09.....	35	50-743-04-09.....	47
50-380-05-71.....	31	50-404-07-71.....	35	50-756-05-09.....	47
50-380-05-91.....	31	50-404-07-91.....	35	50-756-07-09.....	47
50-381-04-09.....	31, 56	50-405-06-71.....	24	50-757-07-09.....	47
50-381-04-71.....	31	50-405-06-91.....	24	50-924-00-07.....	51
50-381-04-91.....	31	50-405-07-09.....	34, 56	50-924-07-07.....	51
50-381-05-09.....	34, 56	50-405-07-71.....	34	50-924-09-07.....	51
50-381-05-71.....	34	50-405-08-09.....	34, 56		
50-382-05-09.....	31	50-405-08-71.....	34	55	
50-382-05-71.....	31	50-405-09-09.....	34, 56	55-961-06-04.....	56
50-382-05-91.....	31	50-405-09-71.....	34	55-961-07-04.....	58
50-383-05-09.....	34, 56	50-405-12-09.....	17, 56	55-961-29-04.....	13
50-383-05-71.....	34	50-405-12-71.....	17	55-961-32-04.....	60
50-384-04-91.....	32	50-405-17-09.....	35, 56	55-962-03-04.....	56, 58, 60
50-386-05-09.....	32	50-405-17-71.....	35	55-962-13-04.....	57, 59
50-386-05-71.....	32	50-406-04-09.....	35	55-962-23-04.....	57, 59, 61
50-386-05-91.....	32	50-406-04-71.....	35	55-962-33-04.....	57, 61
50-386-06-09.....	32, 56	50-406-04-91.....	35	55-962-48-04.....	13
50-386-06-71.....	32	50-406-08-09.....	34, 56	55-963-06-04.....	56
50-386-06-91.....	32	50-406-08-71.....	34	55-963-07-04.....	58
50-389-06-71.....	32	50-406-17-09.....	44	55-963-08-04.....	57, 59, 61
50-389-06-91.....	32	50-406-17-71.....	44	55-963-09-04.....	57, 59
50-390-05-09.....	33, 56	50-407-04-09.....	35	55-963-15-04.....	63
50-390-05-71.....	33	50-407-04-71.....	35	55-963-16-04.....	63
50-390-05-91.....	33	50-407-04-91.....	35	55-963-30-04.....	63
50-391-05-09.....	33	50-407-07-09.....	34, 56	55-963-47-04.....	13
50-391-05-71.....	33	50-407-07-71.....	34	55-963-52-04.....	60
50-391-05-91.....	33	50-407-09-09.....	34, 56	55-964-03-04.....	56
50-392-05-71.....	32	50-407-09-71.....	34	55-964-05-04.....	56
50-392-05-91.....	32	50-408-04-09.....	35, 56	55-964-15-04.....	63
50-396-06-09.....	33, 56	50-408-04-71.....	35	55-964-17-04.....	58
50-396-06-71.....	33	50-408-04-91.....	35	55-964-18-04.....	58
50-396-06-91.....	33	50-409-04-09.....	35, 56	55-964-30-04.....	63
50-397-06-09.....	33	50-409-04-71.....	35	55-964-69-04.....	60
50-397-06-71.....	33	50-409-04-91.....	35	55-964-80-04.....	60
50-397-06-91.....	33	50-410-04-09.....	34	55-969-14-04.....	63
50-399-10-09.....	16	50-410-04-71.....	34	55-969-15-04.....	63
50-399-10-71.....	16	50-410-04-91.....	34	55-969-16-04.....	63
50-399-12-09.....	16	50-410-05-71.....	34	55-969-24-04.....	65
50-399-12-71.....	16	50-410-05-91.....	34	55-969-26-04.....	63
50-400-06-09.....	17	50-412-05-09.....	34	55-969-30-04.....	63
50-400-06-71.....	17	50-412-05-71.....	34	55-969-47-04.....	63
50-400-06-91.....	17	50-412-05-91.....	34	55-969-49-04.....	63
50-400-08-09.....	17	50-425-05-07.....	48	55-969-57-04.....	63
50-400-08-71.....	17	50-501-40-07.....	54, 65	55-969-59-04.....	65
50-400-08-91.....	17	50-502-10-07.....	53, 64, 65	55-969-93-04.....	63
50-401-06-09.....	31	50-502-11-07.....	53, 64, 65	55-969-94-04.....	63
50-401-06-71.....	31	50-502-23-07.....	65	55-969-98-04.....	63
50-401-06-91.....	31	50-502-24-07.....	65	55-969-99-04.....	63
50-402-06-09.....	31	50-502-98-07.....	53		
50-402-06-71.....	31	50-502-99-07.....	53		
50-402-06-91.....	31	50-515-05-07.....	50		
50-402-07-09.....	35	50-735-04-09.....	46		
50-402-07-71.....	35	50-735-14-09.....	46		
50-402-07-91.....	35	50-735-24-09.....	46		
50-403-06-71.....	31	50-735-34-09.....	46		
50-403-06-91.....	31	50-735-36-09.....	46		
50-403-07-09.....	35	50-740-05-09.....	47		
50-403-07-71.....	35	50-740-09-09.....	47		

KLS Martin Group

KLS Martin Australia Pty Ltd.

Sydney · Australia
Tel. +61 2 9439 5316
australia@klsmartin.com

KLS Martin do Brasil Ltda.

São Paulo · Brazil
Tel. +55 11 3554 2299
brazil@klsmartin.com

KLS Martin Medical (Shanghai) International Trading Co., Ltd

Shanghai · China
Tel. +86 21 5820 6251
info@klsmartin.com

KLS Martin India Pvt Ltd.

Chennai · India
Tel. +91 44 66 442 300
india@klsmartin.com

KLS Martin Italia S.r.l.

Milan · Italy
Tel. +39 039 605 67 31
info@klsmartin.com

KLS Martin Japan K.K.

Tokyo · Japan
Tel. +81 3 3814 1431
info@klsmartin.com

KLS Martin SE Asia Sdn. Bhd.

Penang · Malaysia
Tel. +604 261 7060
malaysia@klsmartin.com

KLS Martin de México, S.A. de C.V.

Mexico City · Mexico
Tel. +52 55 7572 0944
mexico@klsmartin.com

KLS Martin Nederland B.V.

Huizen · Netherlands
Tel. +31 35 523 45 38
info@klsmartin.com

KLS Martin SE & Co. KG

Moscow · Russia
Tel. +7 499 792 76 19
russia@klsmartin.com

KLS Martin Taiwan Ltd.

Taipei · Taiwan
Tel. +886 2 2325 3169
taiwan@klsmartin.com

KLS Martin SE & Co. KG

Dubai · United Arab Emirates
Tel. +971 4 454 16 55
middleeast@klsmartin.com

KLS Martin UK Ltd.

Reading · United Kingdom
Tel. +44 118 467 1500
info.uk@klsmartin.com

KLS Martin LP

Jacksonville · Florida, USA
Tel. +1 904 641 77 46
usa@klsmartin.com

KLS Martin SE Asia Sdn. Bhd.

Hanoi · Vietnam
Tel. +49 7461 706-0
vietnam@klsmartin.com

KLS Martin SE & Co. KG

A company of the KLS Martin Group

KLS Martin Platz 1 · 78532 Tuttlingen · Germany
PO Box 60 · 78501 Tuttlingen · Germany
Tel. +49 7461 706-0 · Fax +49 7461 706-193
info@klsmartin.com · www.klsmartin.com