



smartOne®


Monopolar Micro Dissection Needles

The safe and easy Solution
for oromaxillofacial, plastic
and dental Surgery

smartOne® Electrodes – Micro Dissection Needles



The safe and easy Solution
for oralmaxillofacial, plastic and
dental Surgery



Precision as its finest

The smartOne® micro dissection needle electrodes are especially for fine and small incisions. Due to the fine tip of the electrodes, it is possible to precisely prepare the tissue before insertion of the implant.



nonStick tungsten material

All smartOne® micro dissection needles are made of tungsten. This material reduces the tissue from attaching, thus avoiding undesired bleeding.



Full power with low settings

Due to the fine tip and special coating, it is possible to work with low settings on the electrosurgery device.

The perfect solution of saving the tissue of the patient.


IPS Implants®

One Patient. One Solution.



Active electrodes

STERILE | EO

max. 3000 Vp Ø 2.4 mm 

NSTC

NSTC = nonStick tungsten

10 SU = Sales unit/pack

Micro dissection needles



80-700-40-04 10
smartOne® tungsten-coated micro dissection needle



80-700-42-04 10
smartOne® tungsten-coated micro dissection needle, 45° angled



80-700-41-04 10
smartOne® tungsten-coated micro dissection needle



80-700-43-04 10
smartOne® tungsten-coated micro dissection needle, 45° angled

Article overview

Article number	Description	SU
80-700-40-04	smartOne® tungsten-coated micro dissection needle, 20 mm	10
80-700-41-04	smartOne® tungsten-coated micro dissection needle, 30 mm	10
80-700-42-04	smartOne® tungsten-coated micro dissection needle, 40 mm, 45° angled, 5 mm	10
80-700-43-04	smartOne® tungsten-coated micro dissection needle, 40 mm, 45° angled, 10 mm	10

KLS Martin Group

KLS Martin Australia Pty Ltd.

Sydney · Australia
Tel. +61 2 9439 5316
australia@klsmartin.com

KLS Martin do Brasil Ltda.

São Paulo · Brazil
Tel. +55 11 3554 2299
brazil@klsmartin.com

KLS Martin Medical (Shanghai) International Trading Co., Ltd

Shanghai · China
Tel. +86 21 5820 6251
info@klsmartin.com

KLS Martin India Pvt Ltd.

Chennai · India
Tel. +91 44 66 442 300
india@klsmartin.com

KLS Martin Italia S.r.l.

Milan · Italy
Tel. +39 039 605 67 31
info@klsmartin.com

KLS Martin Japan K.K.

Tokyo · Japan
Tel. +81 3 3814 1431
info@klsmartin.com

KLS Martin SE Asia Sdn. Bhd.

Penang · Malaysia
Tel. +604 261 7060
malaysia@klsmartin.com

KLS Martin de México, S.A. de C.V.

Mexico City · Mexico
Tel. +52 55 7572 0944
mexico@klsmartin.com

KLS Martin Nederland B.V.

Huizen · Netherlands
Tel. +31 35 523 45 38
infonl@klsmartin.com

KLS Martin SE & Co. KG

Moscow · Russia
Tel. +7 499 792 76 19
russia@klsmartin.com

KLS Martin Taiwan Ltd.

Taipei · Taiwan
Tel. +886 2 2325 3169
taiwan@klsmartin.com

KLS Martin SE & Co. KG

Dubai · United Arab Emirates
Tel. +971 4 454 16 55
middleeast@klsmartin.com

KLS Martin UK Ltd.

Reading · United Kingdom
Tel. +44 118 467 1500
info.uk@klsmartin.com

KLS Martin LP

Jacksonville · Florida, USA
Tel. +1 904 641 77 46
usa@klsmartin.com

KLS Martin SE Asia Sdn. Bhd.

Hanoi · Vietnam
Tel. +49 7461 706-0
info@klsmartin.com

KLS Martin SE & Co. KG

A company of the KLS Martin Group

KLS Martin Platz 1 · 78532 Tuttlingen · Germany
PO Box 60 · 78501 Tuttlingen · Germany
Tel. +49 7461 706-0 · Fax +49 7461 706-193
info@klsmartin.com · www.klsmartin.com