




## **Twin Pad II**

The new reusable neutral electrode  
The best combination between economy and safe application!

## **Twin Pad II** —

The new reusable neutral electrode

The best combination between economy and safe application!



### **Reusability and safety meets economy!**

The Twin Pad II combines the safety advantages of a split neutral electrode with the economic advantages of reusability.


### **Cost reducing potential**

Due to the possibility of reusability and a high number of reprocessing cycles of the Twin Pad II, there is a high potential to reduce costs in the hospital.

### **Compatible KLS Martin devices:**

In combination with all our KLS Martin generators, like the maXium® or maXium® smart C, it is possible to use our Patient Control System (PCS).

### **NE indicator for highest safety standard**



The PCS system (Patient Control System – a monitoring system for split neutral electrodes) and the NE indicator, not only displays the status of the applied neutral electrode but also of the application quality – continuously throughout the entire application.

The PCS guarantees a high safety standard during all kinds of HF-surgery.



maxium® smart Cart



e-version

i-version

m-version

### Compatibility with all KLS Martin generators!

#### Major electrocautery units:



#### maxium®

80-042-00-04 m-version  
80-042-02-04 i-version  
80-042-04-04 e-version



#### maxium® smart C

80-043-00-04 m-version  
80-043-01-04 i-version  
80-043-02-04 e-version



#### ME MB 3

80-040-11-04 m-version  
80-040-12-04 i-version

#### Minor electrocautery units:



#### ME 102

80-010-02-04



#### miniCutter

80-008-03-04

#### Article overview

Article number	Description
80-342-07-04	Twin Pad II
80-294-40-04	Connection cable KLS Martin devices
80-294-43-04	Connection cable Erbe devices
80-294-44-04	Connection cable international

## KLS Martin Group

### **KLS Martin Australia Pty Ltd.**

Sydney · Australia  
Tel. +61 2 9439 5316  
australia@klsmartin.com

### **KLS Martin do Brasil Ltda.**

São Paulo · Brazil  
Tel. +55 11 3554 2299  
brazil@klsmartin.com

### **KLS Martin Medical (Shanghai)**

International Trading Co. Ltd.  
Shanghai · China  
Tel. +86 21 5820 6251  
china@klsmartin.com

### **KLS Martin India Pvt Ltd.**

Chennai · India  
Tel. +91 44 66 442 300  
india@klsmartin.com

### **Martin Italia S.r.l.**

Milan · Italy  
Tel. +39 039 605 67 31  
italia@klsmartin.com

### **Nippon Martin K.K.**

Tokyo · Japan  
Tel. +81 3 3814 1431  
nippon@klsmartin.com

### **KLS Martin SE Asia Sdn. Bhd.**

Penang · Malaysia  
Tel. +604 505 7838  
malaysia@klsmartin.com

### **KLS Martin de México S.A. de C.V.**

Mexico City · Mexico  
mexico@klsmartin.com

### **Martin Nederland/Marned B.V.**

Huizen · The Netherlands  
Tel. +31 35 523 45 38  
nederland@klsmartin.com

### **Gebrüder Martin GmbH & Co. KG**

Moscow · Russia  
Tel. +7 499 792-76-19  
russia@klsmartin.com

### **KLS Martin Taiwan Ltd.**

Taipei · Taiwan  
Tel. +886 2 2325 3169  
taiwan@klsmartin.com

### **Gebrüder Martin GmbH & Co. KG**

Dubai · United Arab Emirates  
Tel. +971 4 454 16 55  
middleeast@klsmartin.com

### **KLS Martin UK Ltd.**

London · United Kingdom  
Tel. +44 1189 000 570  
uk@klsmartin.com

### **KLS Martin LP**

Jacksonville · Florida, USA  
Tel. +1 904 641 77 46  
usa@klsmartin.com

### **Gebrüder Martin GmbH & Co. KG**

#### **A company of the KLS Martin Group**

KLS Martin Platz 1 · 78532 Tuttlingen · Germany  
P.O. Box 60 · 78501 Tuttlingen · Germany  
Tel. +49 7461 706-0 · Fax +49 7461 706-193  
info@klsmartin.com · www.klsmartin.com

